***Luminex***<sup>®</sup>

complexity simplified.

# **LX200 xPONENT**

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- 70-75 % 异丙醇或 70-75 % 乙醇
- 蒸馏水
- 0.1 M NaOH
- 10 % 家用含氯消毒液
- Luminex® 100/200 Calibration Kit (REF. No. LX200-CAL-K25) 仪器配套校准试剂盒
- Luminex® 100/200 Performance Verification Kit (REF.No.LX200-CON-K25) 仪器配套验证试剂盒
- Sonicator bath 超声水槽/超声波清洗器

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**Automated Maintenance - Calibration Verification**

Instructions Select the Automated Maintenance Option or Maintenance Command. Then, select the appropriate kit if applicable.

Automated Maintenance Options

① Calibration Verification Performance Verification Fluidics Prep System Shutdown

Command	Location	Reagent	Status	Information
Warmup	None	None	Pending	
Prime	None	None	Pending	
Alcohol Flush	1,E7	70% Alcohol	Pending	
Wash	1,B7	DI H2O	Pending	
Wash	1,B7	DI H2O	Pending	
CAL1	1,A10	CAL1	Pending	
MagCAL1	1,B10	MagCAL1	Pending	
CAL2	1,C10	CAL2	Pending	
Wash	1,B7	DI H2O	Pending	
CON1	1,A11	CON1	Pending	
MagCON1	1,B11	MagCON1	Pending	
CON2	1,C11	CON2	Pending	
Wash	1,B7	DI H2O	Pending	
Fluidics1	1,D11	Fluidics1	Pending	
Fluidics2	1,D11	Fluidics2	Pending	

Reagents

Calibration Kit Performance Verification Kit

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③ Run

System Status: Disconnected, Command: System State, Thursday 5/14/2020 9:01 PM

Delta Cal Temp: Laser: Sheath Pressure: Region Events: XY Status: Total Events: Warm Up

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**Automated Maintenance - Performance Verification**

Instructions Select the Automated Maintenance Option or Maintenance Command. Then, select the appropriate kit if applicable.

Automated Maintenance Options

- Calibration Verification
- Performance Verification** (①)
- Fluidics Prep
- System Shutdown

Command	Location	Reagent	Status	Information
Warmup		None	Pending	
Prime		None	Pending	
Alcohol Flush	1,E7	70% Alcohol	Pending	
Wash	1,B7	DI H2O	Pending	
Wash	1,B7	DI H2O	Pending	
CON1	1,A11	CON1	Pending	
MagCON1	1,B11	MagCON1	Pending	
CON2	1,C11	CON2	Pending	
Wash	1,B7	DI H2O	Pending	
Fluidics1	1,D11	Fluidics1	Pending	
Fluidics2	1,E11	Fluidics2	Pending	
Wash	1,B7	DI H2O	Pending	

Reagents

Calibration Kit: B7622 (②) 此处显示批号

Performance Verification Kit: B7623 (②) 此处显示批号

Legend

- AF Alcohol Flush
- Sanitize
- Wash
- Soak
- Drain
- CAL1
- MagCAL1
- CAL2
- CON1
- MagCON1
- CON2
- Fluidics1
- Fluidics2
- Multiple

Vortex each reagent vial for 30 seconds. Place 5 drops of each reagent into the designated well. Fill the designated reservoirs with DI H2O, 70% Isopropanol, 10 - 20% household bleach, etc. if applicable to the chosen routine.

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**Luminex** User: Administrator Log Off Help XONENT SOFTWARE SOLUTIONS

Home Samples Batches Results Protocols **Maintenance** Admin

**Auto Maint**

**Automated Maintenance - System Shutdown**

**Instructions** Select the Automated Maintenance Option or Maintenance Command. Then, select the appropriate kit if applicable.

Automated Maintenance Options

Calibration Verification Performance Verification Fluidics Prep **① System Shutdown**

Command	Location	Reagent	Status	Information
Sanitize	1,H7	10% Bleach	Pending	
Wash	1,B7	DI H2O	Pending	
Wash	1,B7	DI H2O	Pending	
Soak	1,B7	DI H2O	Pending	

Reagents

Calibration Kit: B76229 Performance Verification Kit: B76232

Legend

- Alcohol Flush
- Sanitize
- Wash
- Soak
- Drain
- CAL1
- MagCAL1
- CAL2
- CON1
- MagCON1
- CON2
- Fluidics1
- Fluidics2
- Multiple

Automated Maintenance Plate

DI H2O (B7) 水

70% EtOH (E7)

10% Bleach (H7)

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Report Import Kit **③ Run**

System Status: Disconnected Command: System State: Thursday 5/14/2020 9:08 PM Stop Pause Eject Delta Cal Temp: Laser: Sheath Pressure: Region Events: XY Status: Total Events: Warm Up

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



③ 点击 Support


The screenshot displays the Luminox software interface. At the top, the user is logged in as 'Administrator'. The navigation menu includes Home, Samples, Batches, Results, Protocols, Maintenance, and Admin. The 'Maintenance' menu item is highlighted with a red box. Below the navigation menu, the 'Support Utility' page is active. It features a sidebar with options like Auto Maint, Lot Management, Ccmds & Routines, Probe & Heater, System Info, System Status, and Schedule. The 'Support Utility' option in the sidebar is highlighted with a red box. The main content area shows instructions for using the support utility, a checked checkbox for 'Include Batch Info' (circled in red), and a table of batches. The table has columns for Name, Protocol, Protocol Version, Date, and Status. The 'Support' button at the bottom right of the page is also highlighted with a red box. At the bottom of the interface, the system status bar shows 'Disconnected' with a red circle around it, and the system state is 'Thursday 5/14/2020 9:09 PM'.

Name	Protocol	Protocol Version	Date	Status
AT 2020	SSAAT Assay	1	5/14/2020 9:08 PM	Pending
SAI_012_210420	LS1A04012	1	4/21/2020 11:08 AM	Complete
scov1 plate1 IgG 202004...	20200416 COV2 DuPlex...	2	4/19/2020 3:55 PM	Complete
200417-LSMU	LSMUTR009	1	4/17/2020 11:17 AM	Complete

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## Luminex Big File Transfer

Luminex customers, employees or otherwise affiliated may upload files that are too large to send via e-mail using the following form.

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### Notes

- Either you or the recipient must have an "@luminexcorp.com" email address.
- Multiple recipients should be separated with either semi-colons or commas.
- A download link for the file(s) will be emailed to you and the recipient(s).
- The download links will expire after 14 days.
- You may drag & drop files from your desktop onto this webpage with most browsers.

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**Luminex**<sup>®</sup>  
complexity simplified.