biotechne®

SIMPLE WESTERN SIZE 方法优化

萨日娜 Rina Sa


Field Application Scientist

Tel : 17319114959

E-mail : rina.sa@bio-techne.com

4/17/2020

SIMPLE WESTERN SIZE 方法优化

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Simple Western Size Assay Buffer Compatibility (Wes, Sally Sue, and Peggy Sue)

Lysis Buffer Compatibility Table

LYSIS BUFFER	USAGE	VENDOR & CATALOG	RANGE TESTED	CHEMI SIGNAL	CHEMI RESOLUTION	FLUORESCENT STANDARDS/ MW SIZING	RECOMMENDATIONS
--------------	-------	------------------	--------------	--------------	------------------	----------------------------------	-----------------

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<http://www.proteinsimple.com/antibody/antibodies.html>

Simple Western Antibody Database

The Simple Western Antibody Database is a user-interactive listing of antibodies that have been screened and tested in Simple Western Charge- and Size-based assays. This database is intended to provide assay general development guidance in identifying and selecting antibodies to test.

To bring you even more Simple Western-certified antibodies, we've partnered with our sister companies [R&D Systems](#) and [Novus Biologicals](#) so you don't have to do the work. You can find their antibodies listed below, or for more details visit their websites and search by Application (type "Simple Western"). We've given these antibodies the seal of approval they deserve after being certified in-house by our experienced scientists.

The performance of any antibody selected will need to be optimized for the test system being examined. We encourage users to contribute to the database by clicking the Submit New Antibody button and sharing specific information about the antibodies that have been tested. If you need help with or have questions about any of the information in this antibody database, please contact support@proteinsimple.com.

The database can be sorted by clicking on the column headings. In addition, keywords can be searched and identified from the database using the Search box.

ANTIBODIES (2311)

Go >

Reset

Show All

Submit New Antibody >

Page: << 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19
20 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39
40 41 42 43 44 45 46 47 48 49 50 51 52 53 54 55 56 57 58 59
60 61 62 63 64 65 66 67 68 69 70 71 72 73 74 75 76 77 78 79
80 81 82 83 84 85 86 87 88 89 90 91 92 93 > >>

Protein Target ▲	Antibody Type	Antibody Vendor	Product Number	Protein Isoform	Antibody Species	Cell Model(s)	Separated by Size or Charge	Antibody Dilution	Apparent MW in kDa on Simple Western	Matrix
AGXT	Primary	Novus Biologicals	NBP1-89200	AGXT	Rabbit Polyclonal	Liver (left), HepG2 (right)	Size-Wes, Sally Sue/Peggy Sue	1:20	47	12-230kD
ARG1	Primary	Novus Biologicals	NBP2-14787	ARG1	Mouse Monoclonal	Liver	Size-Wes, Sally Sue/Peggy Sue	1:20	42	12-230kD
ARG1	Primary	Novus Biologicals	NBP1-87490	ARG1	Rabbit Polyclonal	Liver (left), HepG2 (right)	Size-Wes, Sally Sue/Peggy Sue	1:50	41	12-230kD

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Video Article

Use of Capillary Electrophoresis Immunoassay to Search for Potential Biomarkers of Amyotrophic Lateral Sclerosis in Human Platelets

Jessica M. Sage¹, LaSharice Hall², April McVey³, Richard J. Barohn³, Jeffrey M. Statland³, Omar Jawdat³, Mazen M. Dimachkie³, Abdulbaki Agbas¹

The most critical part of CEI is optimizing the assay conditions for each antibody purchased from different vendors, type of antibody (monoclonal vs. polyclonal), optimum protein concentrations, sample preparation, sample denaturation temperature, and electrophoresis voltage applied on the capillaries. We have developed a single-assay format optimization method for the CEI that should be implemented before any new assays, which will save time and resources.

Target protein: **total and phosphorylated** transactivation response DNA/RNA binding protein (TARDP). Due to its size (43 kDa), the acronym **TDP-43** will be used.

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
Figure 2: Signal clarity depends on the sample quality. (A) Whole platelet lysate homogenate interferes with anti-TDP-43 antibody binding; therefore, a noisy electropherogram was observed. (B) Platelet cytosolic fraction was obtained from subjecting the whole lysate to centrifugation (16,000 $\times g$ for 30 min). Most of the membranous proteins were removed; hence, anti-TDP-43 antibody binding to TDP-43 protein was enhanced (blue line trace).

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


Figure 1: Assay layout. Both primary antibody and target protein sample optimization can be performed in one assay. Capillaries 2–7, 8–13, 14–19, and 20–24 represent various protein concentration and primary antibody range. Capillary 25 represents positive control. Anti-ERK antibody was used

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Antibody dilution : 1:300
Sample concentration: 0.5mg/mL

Figure 3: Linear dynamic range for platelet cytosol protein concentration. (A) Both protein concentrations and antibody titrations were optimized in one plate during the same run. (B) The linear working range (0.2–0.8 mg/mL) for protein was established. A 0.5 mg/mL protein load was labeled by α -ERK antibody as the positive control (capillary 25).

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


Figure 5: A representation of predictive phosphorylation value (PPV) of TDP-43. Absolute amount of phosphorylated TDP-43 and pan TDP-43 alone did not show much difference between the ALS and control groups. However, PPV indicated a slight increase in the ALS cohort, although there was no statistical difference between the two groups due to insufficient numbers of subjects (ALS = 25, control = 27). A low Cohen's d value between the means of ALS and control group showed a low effect size between the two groups due to small sample size (control = 25, ALS = 27).

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
FIGURE 5. Comparative data showing 14-3-3 gamma protein expression (blue bars) and the normalized expression (orange bars) in various concentrations of Jurkat cell lysate. Targets were detected by the NIR channel on Jess.




FIGURE 6. Comparative total protein data (blue bars) and the expression of β -actin (orange bars) in six human whole tissue lysates (0.3 mg/mL) tested using Jess.

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042-196 Buy Online	10X System Control Primary Antibody-Rabbit To be used with rabbit primary antibody to Antibody Diluent. Detects system controls. Match with Standard pack 1, 2, 3 or 4 (250 µl/vial)
042-191 Buy Online	10X System Control Primary Antibody-Mouse To be used with mouse primary antibody to Antibody Diluent. Detects system controls. Match with Standard pack 1, 2, 3 or 4 (250 µl/vial)

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
Table of Contents

Chapter 1:

General Information	1
Introduction to Simple Western Size Assays	2
Safety	3
User Attention Notifications	3
Chemical Hazards	4
Safety Data Sheets	4
Customer Service and Technical Support	5

Chapter 2:

Getting Started	7
Simple Western Size Assay Development Overview	8
Choose a Model System	8
General Guidelines	8
Prepare Stock Lysates	9
Prepare Simple Western Samples	9
Screen Antibodies	12
Choosing an Antibody	12
Suggested Protocol for Your Initial Experiment	14
Suggested Plate Templates for Your Initial Experiment	15
Simon	15
Wes	16
Sally and Peggy	17
Sally Sue and Peggy Sue	18
Data Acquisition and Interpretation	19

Chapter 3:

Sample Preparation Optimization	21
Buffer Compatibility for Stock Lysates	22
Preparing Samples for Simple Western Size Assays	22
Preparing Samples in Commercial Sample Buffers	22
Preparing Immunoprecipitation Samples	24
Preparing Samples with More SDS and DTT	24
Preparing Samples with an Additional Component	26
Preparing Samples with an Alternate Reducing Agent	26
Preparing Non-Reduced Samples	26
Alternate Denaturing Conditions	27

Chapter 4:

Simple Western Size Assay Optimization	29
Simple Western Separation Optimization	30
Sample and Stacking Matrix Load Time	30
Separation Time	32
Simple Western Immunoassay Optimization	33
Determining the Linear Dynamic Range	33
Optimizing Antibody Dilution and Incubation Time	35
Accelerating Antibody Binding Kinetics (Simon, Sally, and Peggy)	37

Chemi Exposure Time	39
Using Immunoassay Reagents from Other Vendors	40
Antibody Diluent/Blocking Buffer	40
Secondary Antibody	41
Luminol/Peroxide	42

Chapter 5:

Additional Protocols and Applications	45
1- Incubation Step Protocol	46
Multiplexing	47
Loading Controls	48
Biological Loading Control	49
System Control (Wes, Sally Sue, and Peggy Sue)	49
Normalizing Data in Compass	50

Chapter 6:

Troubleshooting	53
Troubleshooting Flowchart	54
Incorrect or Varying Molecular Weight	55
Signal Burn-Out	56
Separation Issues	58
Increased Assay Variability	58
Bubbles in the Capillary	60
No Sample Separation	61
Fluorescent Standards Don't Separate	62
No Chemiluminescent Signal	63

Appendix A:

Reagent and Lysis Buffer Compatibility	65
Simple Western Reagent Compatibility	66
Buffering Reagents	66
Salts	67
Reducing Agents	67

Denaturing Agents	68
Detergents	69
Fixative	70
Viscosity/Density Additives	71
Miscellaneous Reagents	71

Simple Western Commercial Lysis Buffer Compatibility	72
--	----

Appendix B:

Sample-Specific Lysis Protocols	75
Cell Lines — Guidelines and Preparation	76
Kit Components	76
Additional Consumables/Equipment Needed	76
Important Information	76
Lysis Buffer Preparation	77
Procedure for Lysis of Adherent Cells	77
Procedure for Lysis of Suspension Cells	78
Tissue Homogenization Protocol	79
Recommended Materials	79
Important Information	79
Protocol	79
Fine Needle Aspirate (FNA) — Guidelines and Preparation	81
Practical Considerations	81
Sample Procurement and Processing	81
Preclinical FNA Procedure	81
Clinical FNA Procedure	81
Sample Stability	81
Cell Yields	82
Hemoglobin Contamination	82
Effect of Residual Cell Media on Simple Western Size Assay Performance	82

THANK YOU