



Compass Software User Guide

Copyright © 2015 ProteinSimple. All rights reserved.

ProteinSimple
3001 Orchard Parkway
San Jose, CA 95134
Toll-free: (888) 607-9692
Tel: (408) 510-5500
Fax: (408) 510-5599
email: support@proteinsimple.com
web: www.proteinsimple.com

Compass Software User Guide for Wes, Sally Sue and Peggy Sue

P/N 031-109

Revision 3, August 2015

For research use only. Not for use in diagnostic procedures

Patents and Trademarks

ProteinSimple's Simple Western technology is covered by issued and pending patents in the U.S. and other countries. For more information please see http://www.proteinsimple.com/Intellectual_Property.html.

ProteinSimple and the ProteinSimple logo are trademarks or registered trademarks of ProteinSimple. Other marks appearing in these materials are marks of their respective owners.

Table of Contents

Chapter 1:

Let's Get Started..... 1

Launching Compass..... 2

Compass Overview 2

 Changing the Screen View 2

 Assay Screen 2

 Run Summary Screen 4

 Analysis Screen..... 5

 Screen Panes..... 5

 Title Bar 6

 Main Menu 6

 Instrument Status Bar..... 6

 Screen Tab 6

 View Bar 7

 Compass Status Bar..... 7

Software Menus 7

 File Menu 8

 Edit Menu..... 8

 View Menu..... 8

 Instrument Menu 9

 Window Menu 9

 Help Menu 10

Changing the Compass Main Window Layout 10

 Resizing the Main Compass Window 10

 Resizing the Screen Tab 10

 Resizing Screen Panes..... 11

 Changing the Location of Screen Panes..... 12

 Restoring the Main Window to the Default
 Layout 14

Software Help..... 14

Checking for and Installing New Versions of
Compass 14

Viewing Release Notes 15

Viewing the Software Log..... 15

Compass Version Information 15

Directory and File Information 16

 File Types..... 18

Chapter 2:

Size Assays..... 19

Assay Screen Overview..... 20

 Assay Screen Panes 20

 Software Menus Active in the Assay Screen..... 21

Reagent Color Coding 22

 Immunoassays..... 22

 Total Protein Assays..... 23

Opening an Assay 24

Immunoassays: Creating a New Assay..... 25

 Step 1 - Open a Template Assay..... 25

 Step 2 - Assign Assay Plate Reagents (Optional) . 26

Step 3 - Modifying the Assay Protocol (Optional for All Instruments)	31
Step 4 - Add Assay Notes (Optional)	34
Step 5 - Select a Schedule (Optional for Sally Sue and Peggy Sue)	35
Step 6 - Add Assay Plate Annotations (Optional)	36
Step 7 - Save the Assay	39
Step 8 - Modify Default Analysis Parameters (Optional)	39
Total Protein Assays: Creating a New Assay	41
Step 1 - Open a Template Assay	41
Step 2 - Assign Assay Plate Reagents (Optional) .	43
Step 3 - Modifying the Assay Protocol (Optional for All Instruments)	47
Steps 4 - 8	50
Making Changes to an Existing Assay	51
Switching Between Open Assays	51
Creating a Template Assay	52
Viewing and Changing the Detection Exposures	53
Copying Protocols and Templates	54
Copying an Assay Protocol	54
Copying an Assay Template	55
Printing Protocols and Templates	55
Printing an Assay Protocol	55
Printing an Assay Template	55
Importing and Exporting Protocols and Templates ..	56
Importing an Assay Protocol	56
Exporting an Assay Protocol	57
Template Export and Import	58
Exporting the Template to a CSV file	58
Template Cut and Paste	60
Chapter 3:	
Running a Size Assay	63
Starting a Run	64
Step 1 - Get Ready	64
Step 2 - Start the Run	64
Step 3 - Post-Run Procedures	71
Stopping a Run	72
Chapter 4:	
Charge Assays	73
Assay Screen Overview	74
Assay Screen Panes	74
Software Menus Active in the Assay Screen	75
Reagent Color Coding	77
Opening an Assay	78
Creating a New Assay	79
Step 1 - Open a Template Assay	79
Step 2 - Assign Assay Plate Reagents (Optional) .	80
Step 3 - Modifying the Assay Protocol (Optional)	83
Step 4 - Add Assay Notes (Optional)	88
Step 5 - Select a Schedule (Optional)	89
Step 6 - Add Assay Plate Annotations (Optional)	90
Step 7 - Save the Assay	94
Step 8 - Modify Default Analysis Parameters (Optional)	95
Making Changes to an Existing Assay	96
Switching Between Open Assays	96
Creating a Template Assay	97
Viewing and Changing the Detection Exposures	98
Copying Protocols and Templates	99
Copying an Assay Protocol	99
Copying an Assay Template	100

<i>Printing Protocols and Templates</i>	100
<i>Printing an Assay Protocol</i>	100
<i>Printing an Assay Template</i>	100
<i>Importing and Exporting Protocols and Templates</i>	101
<i>Importing an Assay Protocol</i>	101
<i>Exporting an Assay Protocol</i>	102
<i>Importing an Assay Template</i>	102
<i>Exporting an Assay Template</i>	103
Chapter 5:	
Running a Charge Assay on Peggy Sue	105
<i>Starting a Run</i>	106
<i>Step 1 - Get Ready</i>	106
<i>Step 2 - Start the Run</i>	106
<i>Step 3 - Post-Run Procedures</i>	113
<i>Stopping a Run</i>	114
Chapter 6:	
Run Status	115
<i>Run Summary Screen Overview</i>	116
<i>Run Summary Screen Panes</i>	116
<i>Software Menus Active in the Run Summary Screen</i>	117
<i>Opening Run Files</i>	119
<i>Opening One Run File</i>	119
<i>Opening Multiple Run Files</i>	119
<i>Viewing File and Run Status Information</i>	120
<i>Assay Steps: Size-based Assays</i>	121
<i>Assay Steps: Charge-based Assays</i>	123
<i>Watching Standards Separation Movies</i>	125
<i>Viewing Current and Voltage Plots</i>	127
<i>Switching Between Open Run Files</i>	128
<i>Closing Run Files</i>	128
Chapter 7:	
Controlling Wes, Sally Sue and Peggy Sue	129
<i>Instrument Control</i>	130
<i>Starting a New Run</i>	130
<i>Opening Trays (Sally Sue and Peggy Sue)</i>	130
<i>Cleaning</i>	131
<i>Self Test</i>	134
<i>Viewing and Changing System Properties</i>	134
<i>Viewing Log Files</i>	135
<i>Error Logs</i>	135
<i>Self Test Logs</i>	138
<i>Status Modes</i>	140
Chapter 8:	
Size Assay Data Analysis	141
<i>Analysis Screen Overview</i>	143
<i>Analysis Screen Panes</i>	143
<i>Software Menus Active in the Analysis Screen</i>	145
<i>Opening Run Files</i>	147
<i>Opening One Run File</i>	147
<i>Opening Multiple Run Files</i>	147
<i>How Run Data is Displayed in the Analysis Screen</i>	150
<i>Experiment Pane: Assay and Capillary Information</i>	150
<i>Graph Pane: Electropherogram Data</i>	151
<i>Image Pane: Capillary Separation Image Data</i>	152
<i>Lane Pane: Virtual Blot-Like Image Data</i>	153
<i>Std Curve Pane: Standard Curve Fit Data</i>	154
<i>Peaks Pane: Calculated Results</i>	155
<i>Capillaries Pane: User-Specified Peak Names</i>	157
<i>Viewing Run Data</i>	160

Switching Between Sample, Standards and Registration Data Views.....	160	Adding or Removing Sample Data.....	196
Selecting and Displaying Capillary Data.....	165	Hiding Sample Data.....	198
Switching Between Single and Multiple Views of the Capillaries.....	169	Changing Peak Names for Sample Data.....	199
Hiding Capillary Data.....	173	Displaying Sample Data for Named Peaks Only.....	201
Setting Run Data Display Filters.....	174	Changing the Virtual Blot View.....	203
Compass Run Data Notifications and Warnings....	176	Adjusting the Contrast.....	203
Checking Your Results.....	177	Inverting the Virtual Blot.....	204
Step 1 – Review the Fluorescent Sizing Standards Movie.....	177	Selecting Lane Labels.....	204
Step 2 – Checking Fluorescent Sizing Standards.....	178	Viewing the Uncorrected Sample Baseline....	206
Step 3 – Checking Capillary Registrations (Sally Sue and Peggy Sue Only).....	182	Overlaying Standards Data on Sample Lanes.	206
Step 4 – Checking the Ladder.....	183	Moving Lanes in the Virtual Blot View.....	208
Step 5 – Checking Samples.....	185	Changing the Electropherogram View.....	209
Step 6 – Assigning Peak Names (Optional).....	188	Autoscaling the Electropherogram.....	210
Group Statistics.....	188	Stacking Multiple Electropherograms.....	211
Using Groups.....	189	Overlaying Multiple Electropherograms.....	212
Viewing Statistics.....	190	Zooming.....	213
Hiding or Removing Capillaries in Group Analysis.....	192	Customizing the Data Display.....	214
Copying Data Views and Results Tables.....	194	Selecting Data Viewing Options.....	223
Copying Data Views.....	194	Adding and Removing Baseline Points.....	228
Copying Results Tables.....	194	Selecting the X-Axis Molecular Weight Range.	229
Saving the Graph View as an Image File.....	194	Closing Run Files.....	231
Exporting Run Files.....	195	Compass Analysis Settings Overview.....	232
Exporting Results Tables.....	195	Advanced Analysis Settings.....	234
Exporting Raw Sample Electropherogram Data.....	195	Standards Settings.....	235
Changing Sample Protein Identification.....	196	Sample Settings.....	235
		Image Settings.....	235
		Advanced Analysis Settings Groups.....	236
		Creating a New Analysis Group.....	236
		Changing the Default Analysis Group.....	237
		Modifying an Analysis Group.....	238

Deleting an Analysis Group	238	Changing the Default Standards Group	280
Applying Analysis Groups to Specific Run Data	238	Modifying a Standards Group.....	281
Images Analysis Settings.....	241	Deleting an Analysis Group	281
Exposure Settings	242	Applying Analysis Groups to Specific Run Data	282
Changing the Sample Data Exposure		Importing and Exporting Analysis Settings.....	284
Displayed	243	Importing Analysis Settings	284
Peak Fit Analysis Settings	244	Exporting Analysis Settings.....	284
Range Settings	245	Chapter 9:	
Baseline Settings.....	245	Charge Assay Data Analysis	287
Peak Find Settings	245	Analysis Screen Overview	288
Peak Fit Analysis Settings Groups	246	Analysis Screen Panes.....	288
Creating a New Peak Fit Group.....	247	Software Menus Active in the Analysis Screen .	290
Changing the Default Peak Fit Group	248	Opening Run Files	292
Modifying a Peak Fit Group	249	Opening One Run File.....	292
Deleting a Peak Fit Group	249	Opening Multiple Run Files.....	292
Applying Peak Fit Groups to Specific Run Data	249	How Run Data is Displayed in the Analysis Screen .	295
Peak Names Settings	252	Experiment Pane: Assay and Capillary	
Peak Names Analysis Settings Groups.....	253	Information.....	295
Creating a Peak Names Group	253	Graph Pane: Electropherogram Data	296
Adding Peak Names Groups	257	Image Pane: Capillary Separation Image Data	297
Modifying a Peak Names Group	258	Lane Pane: Virtual Blot-Like Image Data	298
Deleting a Peak Names Group	258	Peaks Pane: Calculated Results.....	299
Applying Peak Names Groups to Run Data	259	Capillaries Pane: User-Specified Peak Names..	300
Standard Curve Settings.....	263	Viewing Run Data	302
Applying Peak Names Groups to Run Standard		Switching Between Sample, Standards and	
Curve	265	Registration Data Views.....	302
System or Loading Control Settings.....	268	Selecting and Displaying Capillary Data	307
Standards Settings.....	270	Switching Between Single and Multiple Views of	
Standards Analysis Settings Groups.....	272	the Capillaries.....	311
Changing the Capillary Used for the Ladder ...	273	Hiding Capillary Data.....	315
Creating a New Standards Group.....	275	Setting Run Data Display Filters.....	316

<i>Compass Run Data Notifications and Warnings</i>	318	<i>Viewing the Uncorrected Sample Baseline</i>	345
<i>Checking Your Results</i>	319	<i>Overlaying Standards Data on Sample Lanes</i>	346
<i>Step 1 – Review the Fluorescent Sizing Standards Movie</i>	319	<i>Moving Lanes in the Virtual Blot View</i>	348
<i>Step 2 – Checking Fluorescent Sizing Standards</i>	320	<i>Changing the Electropherogram View</i>	349
<i>Step 3 – Checking Capillary Registrations</i>	324	<i>Autoscaling the Electropherogram</i>	350
<i>Step 4 – Checking Samples</i>	325	<i>Stacking Multiple Electropherograms</i>	351
<i>Step 5 – Assigning Peak Names (Optional)</i>	327	<i>Overlaying Multiple Electropherograms</i>	352
<i>Group Statistics</i>	328	<i>Zooming</i>	353
<i>Using Groups</i>	328	<i>Customizing the Data Display</i>	354
<i>Viewing Statistics</i>	329	<i>Selecting Data Viewing Options</i>	363
<i>Hiding or Removing Capillaries in Group Analysis</i>	332	<i>Adding and Removing Baseline Points</i>	369
<i>Copying Data Views and Results Tables</i>	334	<i>Selecting the X-Axis pI Range</i>	370
<i>Copying Data Views</i>	334	<i>Closing Run Files</i>	371
<i>Copying Results Tables</i>	334	<i>Compass Analysis Settings Overview</i>	372
<i>Saving the Graph View as an Image File</i>	334	<i>Advanced Analysis Settings</i>	374
<i>Exporting Run Files</i>	335	<i>Standards Settings</i>	375
<i>Exporting Results Tables</i>	335	<i>Sample Settings</i>	375
<i>Exporting Raw Sample Electropherogram Data</i>	335	<i>Image Settings</i>	375
<i>Changing Sample Protein Identification</i>	336	<i>Advanced Analysis Settings Groups</i>	376
<i>Adding or Removing Sample Data</i>	336	<i>Creating a New Analysis Group</i>	376
<i>Hiding Sample Data</i>	338	<i>Changing the Default Analysis Group</i>	377
<i>Changing Peak Names for Sample Data</i>	339	<i>Modifying an Analysis Group</i>	378
<i>Displaying Sample Data for Named Peaks Only</i>	341	<i>Deleting an Analysis Group</i>	378
<i>Changing the Virtual Blot View</i>	343	<i>Applying Analysis Groups to Specific Run Data</i>	378
<i>Adjusting the Contrast</i>	343	<i>Images Analysis Settings</i>	381
<i>Inverting the Virtual Blot</i>	344	<i>Exposure Settings</i>	382
<i>Selecting Lane Labels</i>	344	<i>Changing the Sample Data Exposure Displayed</i>	383
		<i>Peak Fit Analysis Settings</i>	384
		<i>Range Settings</i>	385
		<i>Baseline Settings</i>	385

<i>Peak Find Settings</i>	385	Chapter 11:	
<i>Peak Fit Analysis Settings Groups</i>	385	Compass Access Control and	
<i>Creating a New Peak Fit Group</i>	386	21 CFR Part 11 Compliance	423
<i>Changing the Default Peak Fit Group</i>	387	<i>Overview</i>	424
<i>Modifying a Peak Fit Group</i>	387	<i>Enabling Access Control</i>	425
<i>Deleting a Peak Fit Group</i>	388	<i>Logging In to Compass</i>	426
<i>Applying Peak Fit Groups to Specific Run Data</i>	388	<i>Resolving Log In Issues</i>	427
<i>Peak Names Settings</i>	391	<i>Saving Changes</i>	427
<i>Peak Names Analysis Settings Groups</i>	392	<i>Signing Files</i>	428
<i>Creating a Peak Names Group</i>	392	<i>Instrument Command Log</i>	429
<i>Adding Peak Names Groups</i>	396	<i>Run File History</i>	431
<i>Modifying a Peak Names Group</i>	397	<i>Troubleshooting Problems and Suggested</i>	
<i>Deleting a Peak Names Group</i>	397	<i>Solutions</i>	431
<i>Applying Peak Names Groups to Run Data</i>	398	<i>Authorization Server</i>	432
<i>Standards Settings</i>	401	<i>Server Administration</i>	433
<i>Standards Analysis Settings Groups</i>	402	<i>Adding Non-admin Users</i>	434
<i>Creating a New Standards Group</i>	403	<i>Adding Admin Users</i>	439
<i>Changing the Default Standards Group</i>	406	<i>Resetting User Passwords</i>	440
<i>Modifying a Standards Group</i>	406	<i>Encryption Details</i>	441
<i>Deleting an Analysis Group</i>	407		
<i>Applying Analysis Groups to Specific Run Data</i>	407		
<i>Importing and Exporting Analysis Settings</i>	409		
<i>Importing Analysis Settings</i>	409		
<i>Exporting Analysis Settings</i>	410		
Chapter 10:			
Setting Your Preferences	411		
<i>Custom Preference Options</i>	412		
<i>Setting Data Export Options</i>	413		
<i>Selecting Custom Plot Colors for Graph Overlay</i>	414		
<i>Setting Up Wes, Sally Sue and Peggy Sue to Send</i>			
<i>Tweets</i>	416		

Chapter 1:

Let's Get Started

Chapter Overview

- Launching Compass
- Compass Overview
- Software Menus
- Changing the Compass Main Window Layout
- Software Help
- Checking for and Installing New Versions of Compass
- Viewing Release Notes
- Viewing the Software Log
- Compass Version Information
- Directory and File Information

Launching Compass



To open Compass, double-click the icon on the computer desktop.


Compass Overview

Compass has three main screens:

- **Assay** - You'll create and review your assay.
- **Run Summary** - Check out the status of your run.
- **Analysis** - Take a look at the data from your experiment.

Changing the Screen View

To move between the Assay, Run Summary and Analysis screens, just click the button in the screen tab located in the upper right corner of the main window.



Assay Screen

The Assay screen is used to create, view, and edit assays. You can assign well locations for assay plate reagents, modify assay protocol steps, enter assay notes and add annotations for individual wells on the assay plate.

The screenshot displays the Compass software interface for an assay. At the top, the window title is "2013-10-27_14-50-57_Wes Assay - Compass". The menu bar includes "File", "Edit", "Instrument", "Window", and "Help". The status bar shows "Ready" and a green "Start" button. On the right, there are buttons for "Assay", "Run Summary", and "Analysis".

The main interface is divided into several sections:


- Layout:** A small grid showing lanes 1 through 25 and rows A through I.
- Protocol:** A table listing assay parameters and their values.

	Value
Separation Matrix	
Stacking Matrix	
Sample	
Separation Time (min)	30.0
Separation Voltage (volts)	375
Matrix Removal	
Antibody Diluent Time (min)	5.0
Primary Antibody Time (min)	30.0
Secondary Antibody Time (min)	30.0
Detection	
- Template:** A large grid showing the assay layout with 25 lanes and 5 rows (A-E).

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
A	Calyculin-treated Jurkat																								
B	Antibody Diluent																								
C	Anti...	9272	4691	4058	4060	ab81283						9272	4691	4058	4060	ab81283									
D	Secondary Antibody																								
E	Luminol/Peroxide																								


Run Summary Screen

The Run Summary screen is used to monitor status of a run in progress, watch movies of the separation in the capillaries, and view current and voltage plots for each run.



Analysis Screen

The Analysis screen is used to view data from your assay, including the graph view (electropherograms), lane view (shown below), capillary images, and a table with your results. You can also analyze your data here after the run is finished.



Screen Panes

Assay, Run Summary and Analysis screens all have multiple screen panes that let you view the individual components of a run, assay or data file. Each pane has a labeled tab and a unique icon. We'll describe panes specific to each screen later in the individual screen sections.

The active pane in a screen is blue. To view a pane, click in the pane or on its tab. The example below shows panes in the Analysis screen, and the Lane pane is active:




Title Bar

In the title bar you will see the run file name and the icons that allow the main Compass window to be minimized, maximized or closed.




Main Menu

Access to various software, instrument and screen operations is available through the main menu. More details on menu commands can be found in "Software Menus" on page 7.



Instrument Status Bar


The instrument status bar is used to start runs and cleaning protocols, relay system status and show run progress. More details on instrument control and status can be found in Chapter 7, "Controlling Wes, Sally Sue and Peggy Sue".



NOTE: You will only see the instrument status bar when Compass is connected to an instrument. There is no status bar on computer workstations that you're only using for data analysis.


Screen Tab

The screen tab lets you move between Assay, Run Summary or Analysis screens and is located in the upper right corner of the main window. Just click a button to view a screen.



View Bar

The view bar is only displayed in the Analysis screen as part of the main menu bar and allows you to switch between displaying sample chemiluminescent data, fluorescent standards or capillary registration information, data for a single capillary or all capillaries in the run, or grouped capillary data. View bar options are detailed in “Switching Between Sample, Standards and Registration Data Views” on page 160 for size assays or page 302 for charge assays, and “Using Groups” on page 189 for size assays or page 328 for charge assays.



Compass Status Bar

The status bar is located in the lower right corner of the main window. It displays active software processes and their progress.




Software Menus

A brief description of the software menus in the main menu are described in this section. Not all menus are available in every screen, and menu commands change depending on what screen is active. The menus and commands available for each screen will be detailed in the individual screen sections.


File Menu

The File menu contains basic file commands.




Edit Menu

The Edit menu contains basic editing commands, analysis and preferences options. Specific details on preferences are described in Chapter 10, "Setting Your Preferences".




View Menu

The View menu can only be seen in the Analysis screen, and allows you to change how your data is displayed.




Instrument Menu

The Instrument menu is only available when Compass is connected to directly to your instrument. Instrument control options are explained in Chapter 7, *“Controlling Wes, Sally Sue and Peggy Sue”*.



Window Menu


The Window menu allows you to switch between Assay, Run Summary or Analysis screens, and restore screens to the Compass default layout.



- **Assay** - Displays the Assay screen where you can create, view, and edit assays.
- **Run Summary** - Displays the Run Summary screen which tells you the status of a run in progress.
- **Analysis** - Displays the Analysis screen which is used to view sample electropherograms, lane data and results.
- **Default Layout** - Restores the individual panes in the current screen back to their default size and location.

Help Menu

The Help menu provides access to Help, software updates, release notes and other software information.



- **Help Contents** - Displays the Compass Help file.
- **Check for Updates** - Automatically checks to see if a new version of Compass is available.
- **Release Notes** - Displays Compass release notes for the current and prior versions.
- **Compass Log** - Displays the Compass software log file.
- **About Compass** - Displays the Compass software version and build information.

Changing the Compass Main Window Layout


You can easily resize the Compass main window and the individual panes in each screen. Screen panes can also be moved outside of the main window.

Resizing the Main Compass Window

To resize the main window, roll the mouse over a corner or border until the sizing arrow appears. Then just click and drag to resize.


Resizing the Screen Tab

The screen tab can be sized to show all or just some of the screen buttons. To resize, roll the mouse over the left edge of the tab until the sizing arrow appears, then click and drag to resize. If a screen button is hidden, a double arrow will display in the tab. Click to display and select the hidden screen.



Resizing Screen Panes

- **To resize a pane** - Roll the mouse over the pane border until the sizing arrow appears. Then just click and drag to resize.
- **To maximize a pane** - Click the maximize button in the upper right corner or double-click the tab. The other panes in the screen will automatically minimize to pane bars in the task area along the window border.




- **To restore all minimized panes** - Click **Restore** on the minimized pane bar.



- **To restore only one minimized pane** - Click the pane icon on the minimized pane bar.



- **To restore a maximized pane to its original size** - Double-click the tab or right click the tab and click **Restore**.



- **To restore all panes to their original sizes** - Select **Window** in the main menu and click **Default Layout**.

Changing the Location of Screen Panes

Panes can be moved to different locations within a screen.

- **To move a pane** - Click on its tab and drag it to the new location. As the pane is moved, area guides will display to assist you in choosing a drop location.



Area guides with a black arrow let you know that if the pane is dropped at that location, it will be resized and relocated as an individual pane in that area of the screen.




Area guides with a folder let you know that if the pane is dropped at that location, it will be added as a new tab in an area with one or more pane tabs.




Area guides with a window let you know that if the pane is dropped at that location, it will become a separate window outside the Compass main window.

The following figure shows the Analysis screen after moving the Graph pane.



- **To detach a pane from the main window** - Click on its tab and drag it outside the main Compass window or right click the tab and click **Detached**.



- **To move a detached pane back inside the main window** - Right click the tab and deselect **Detached**.
- **To restore all panes to their original locations** - Select **Window** in the main menu and click **Default Layout**.

Restoring the Main Window to the Default Layout

To restore screen pane sizes and locations to the original Compass layout, select **Window** from the main menu and click **Default Layout**.


Software Help

Select **Help** and click **Help Contents** to view the Compass Help file.

Checking for and Installing New Versions of Compass

Compass can automatically check to see if a newer version of software is available. To do this:

1. Make sure the computer being used has an active internet connection.
2. Select **Help** and click **Check for Updates**. If an update is found, the following screen will display:



3. Click **Finish** to start the download and install the update.
4. Follow the on-screen instructions to complete the software installation.
5. Reboot the computer before using the new version of Compass.

Viewing Release Notes

Select **Help** and click **Release Notes** to view feature updates and bug fixes for new and past versions of Compass. We recommend you review these notes whenever a software update is installed.




Viewing the Software Log

Select **Help** and click **Compass Log** to view the software log file.


Compass Version Information

Select **Help** and click **About Compass** to view the software version and build number information.




Directory and File Information

The main Compass directory is located in the **Program Files** folder, and also contains PDF files of the Wes and Sally Sue/Peggy Sue User Guides.



Compass assay and run files are located in the **Documents** folder in the User directory on your computer:



- **Assays Folder** - Contains all assay files that you've saved.
- **New Assays Folder** - Contains Simple Western assay template files.
- **Runs Folder** - Contains all run files. Run data is automatically written to this folder.

NOTE: When a Compass software update is performed, the template assays in the New Assays folder are overwritten. If you have customized these assays, we recommend saving them in a unique subfolder prior to updating the software, then transferring them back to the New Assays folder after the update to avoid losing your assay customizations.

File Types

The following file types are used by Compass:

- **Assay Files** - Use an *.assay file extension.
- **Run Files** - Use a *.cbz file extension. The default file format for run files is Date_Time_AssayName. An example run file name would be 2012-09-28_18-50-53_Simple Western.cbz.
- **Protocol Files** - Exported protocol files use a *.protocol file extension.
- **Template Files** - Exported template files use a *.template file extension.
- **Analysis Settings Files** - Exported analysis settings files use a *.settings file extension.

Chapter 2:

Size Assays

Chapter Overview

- Assay Screen Overview
- Reagent Color Coding
- Opening an Assay
- Immunoassays: Creating a New Assay
- Total Protein Assays: Creating a New Assay
- Making Changes to an Existing Assay
- Switching Between Open Assays
- Creating a Template Assay
- Viewing and Changing the Detection Exposures
- Copying Protocols and Templates
- Printing Protocols and Templates
- Importing and Exporting Protocols and Templates

Assay Screen Overview

The Assay screen is used to create, view, and edit assays. To access this screen, click **Assay** in the screen tab:



Assay Screen Panes

The Assay screen has four panes:

- **Layout** - Displays a map of the assay plate and shows where assay reagents will be located.
- **Protocol** - Lists individual assay protocol steps and parameters that Wes, Sally Sue or Peggy Sue will execute for each of the capillaries.
- **Notes** - Lets you enter specific assay information that is saved with the assay and can be used for future reference.
- **Template** - Enter annotations for the individual well and row reagents in the assay plate.

The screenshot shows the 'Assay: Sally Sue Size' window. It features a 'Layout' pane with a 96-well plate map, a 'Protocol' pane with a table of parameters, and a 'Template' pane with a reagent grid.

	Cycle 1	Cycle 2	Cycle 3	Cycle 4	Cycle 5	Cycle 6	Cycle 7	Cycle 8
Separation Matrix								
Stacking Matrix								
Sample								
Separation Time (min)	40.0	40.0	40.0	40.0	40.0	40.0	40.0	40.0
Separation Voltage (volts)	250	250	250	250	250	250	250	250
Matrix Removal								
Antibody Diluent Time (min)	23.0	23.0	23.0	23.0	23.0	23.0	23.0	23.0
Primary Antibody Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Secondary Antibody Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Detection								

	1	2	3	4	5	6	7	8	9	10	11	12
A	Biot. Ladder						Sample					
B							Antibody Diluent					
C	Blocking						Primary Antibody					
D	Streptavidin ...						Secondary Antibody					
J												

Software Menus Active in the Assay Screen


You can use the following software menus:

- File
- Edit
- Instrument (when Compass is connected to Wes, Sally Sue or Peggy Sue)
- Window
- Help

The File and Edit menu options specific to the Assay screen are described next.

File Menu


The following File menu options are active:



- **New Assay** - Creates a new assay from a starter template.
- **Open Assay** - Opens an existing assay.
- **Save** - Saves the open assay.
- **Save As** - Saves the open assay under a different file name.
- **Import Protocol** - Imports a saved protocol file into an assay.
- **Import Template** - Imports a saved template file into an assay.
- **Export Protocol** - Exports the current protocol file for future use.
- **Export Template** - Exports the current template file for future use.
- **Print** - Prints the information in the Protocol or Template panes.
- **Exit** - Closes Compass.

Edit Menu

The following Edit menu options are active:




- **Copy** - Copies the information in the Protocol or Template panes into other documents.
- **Default Analysis** - Displays the default settings that will be used to analyze the run data generated with an assay.
- **Analysis** - Not active in this screen.
- **Preferences** - Set and save your preferences for data export, plot colors in the graph and Twitter settings. See Chapter 10, "Setting Your Preferences" for more information.

Reagent Color Coding

Immunoassays

The Assay screen uses color coding to identify various assay reagents in all panes. The Layout pane is shown in the following example:




- **Orange** - Samples and Ladder
- **Magenta** - Blocking reagent (Antibody Diluent)
- **Light Teal** - Primary antibody
- **Teal** - Secondary HRP conjugate
- **Gold** - Luminol/Peroxide mix

- **No color coding** - Separation Matrix (clearly designated)
- **No color coding** - Stacking Matrix (clearly designated)
- **Light Blue** - Water dispensed around Separation and Stacking Matrices

Total Protein Assays

The Assay screen uses color coding to identify various assay reagents in all panes. The Layout pane is shown in the following example:




- **Orange** - Samples and Ladder
- **Magenta** - Labeling reagent
- **Light Teal** - Antibody Diluent
- **Teal** - Total Protein Streptavidin-HRP
- **Gold** - Luminol/Peroxide mix
- **No color coding** - Separation Matrix (clearly designated)
- **No color coding** - Stacking Matrix (clearly designated)
- **Light Blue** - Water dispensed around Separation and Stacking Matrices

Opening an Assay

To open an existing assay:

1. Select **File** in the main menu and click **Open Assay**.




2. A list of the last five assays opened will display. Select one of these assays or click **Browse** to open the Assay folder and select a different assay.
3. To make changes to the assay, follow the instructions under "Immunoassays: Creating a New Assay" on page 25 for changing the assay parameters. Select **File** from the main menu and click **Save**.

Immunoassays: Creating a New Assay



To create a new assay, we recommend using one of the template assays and making any necessary modifications from there.

Step 1 - Open a Template Assay


1. Select **File** in the main menu and click **New Assay**:



2. Select the template assay for your instrument by selecting **Size** as your assay type, **Size Range** and **Cartridge Type** if running Wes, or choose **Open Assay** to select from the menu of saved assays.



If you're running a new size assay with split Running Buffer, the window has a Split Running Buffer check box to indicate the use of split buffer. This box is checked by default.




NOTES:

Compass v2.7 or higher and the latest embedded are needed to run the split Running Buffer assays.

If you make changes to the assay, such as defining the wells in the plate, be sure to select **File** and **Save** before proceeding.


Split Running Buffer will be noted above the protocol steps and plate layout if that assay was selected for Sally Sue or Peggy Sue.




Step 2 - Assign Assay Plate Reagents (Optional)

1. Click on the Layout tab. Displayed are the default row locations where each reagent should be placed on the assay plate.

NOTE: For Sally Sue and Peggy Sue, up to 96 different samples or conditions (lysates, proteins, antibodies or incubation times) can be assayed using the full eight cycles to accommodate experimental needs.

Wes plate layout:

- **Row A** - Biotinylated Ladder (A1: 12-230 kDa or 66-440 kDa) and Sample (A2-25)
- **Row B** - Blocking (Antibody Diluent)
- **Row C** - Antibody Diluent (C1) and Primary antibody (C2-C25)
- **Row D** - Streptavidin-HRP (D1) and Secondary HRP conjugate (D2-D25)
- **Row E** - Luminol-S/Peroxide mix
- **Row G** - Wash buffer
- **Row H** - Wash buffer
- **Row I** - Wash buffer
- **Row M-P** - Pre-filled Separation Reagents


Sally Sue/Peggy Sue plate layout:

- **Row A** - Biotinylated Ladder (A1: 12-230 kDa or 66-440 kDa) and Sample (A2-12)
- **Row B** - Blocking (Antibody Diluent)
- **Row C** - Antibody Diluent (C1) and Primary antibody (C2-C12)
- **Row D** - Streptavidin-HRP (D1) and Secondary HRP conjugate (D2-D12)
- **Row J** - Luminol-S/Peroxide mix


- **Row M** - Water (M5-M20)
- **Row N** - Water (N5-N6 and N19-N20) and Separation Matrix (N7-N18)
- **Row O** - Water (O5-O6 and O19-O20) and Stacking Matrix (O7-O18)
- **Row P** - Water (P5-P20)

NOTE: For details on sample, reagent and assay plate preparation, please refer to the product insert provided with the Simple Western kits.


2. If needed, well assignments can be modified as described below. Any row assignments changed in the **Layout** pane are updated in the Protocol pane automatically.
 - **To move a reagent row to another location** - Click the row in the **Layout** pane, then drag and drop it on the new location. The row from which it was moved will be reassigned as empty.




- **To insert a sample row** - Click an empty row or the row below where the new sample row should be inserted, then click the **S** icon (Insert a sample row) in the **Layout** pane toolbar. A new sample row will be added in the empty row or inserted above the selected row.



- **To insert a detection row** - Click an empty row or the row below where the new detection row should be inserted, then click the **D** icon (Insert a detection row) in the **Layout** pane toolbar. A new detection row will be added in the empty row or inserted above the selected row.



- **To insert a fourth or fifth incubation reagent** - Click an empty row or the row below where the new incubation reagent should be inserted, then click the **4** icon (Insert a quaternary row) in the **Layout** pane toolbar. A new incubation reagent row will be added in the empty row or inserted above the selected row.




A fifth incubation reagent can now be added by repeating the above and clicking the **5** icon (Insert a Quinary row) instead.

NOTES:

When using the fourth and fifth incubation steps, only five cycles can typically be performed per run as more water is needed per cycle for the additional wash steps.

Row J is the last row assignment that can be used on the assay plate. New, sample, incubation or luminol rows cannot be inserted if row J already has an assigned reagent. To insert a row in this case, you must first move the contents of row J to another area on the plate.

- **To delete a row** - Click the row to be deleted, then click the red **X** icon (Delete a row) in the toolbar. Only empty rows and fourth incubation rows can be deleted. Rows for required assay reagents cannot be deleted.



NOTES:

For Sally Sue and Peggy Sue, you can put samples, antibodies and blocking buffer in Rows A-J and in columns 1-12 or 13-24.

We recommend keeping two rows empty between the substrate and the closest dispensed HRP conjugate to avoid any interaction with the substrate.


Rows K and L are purposely left empty, don't put anything into these wells.

Separation and Stacking matrices can only be pipetted into the designated wells.

Water is used in the wells to reduce the potential of evaporation of the matrices.

Step 3 - Modifying the Assay Protocol (Optional for All Instruments)

1. Click on the Protocol tab. This pane shows the individual steps of the assay protocol and allows you to change parameters. When creating a new assay, the default protocol automatically assigns all reagent locations for the assay. The Sally Sue assay is shown as an example:



	Cycle 1	Cycle 2	Cycle 3	Cycle 4	Cycle 5	Cycle 6	Cycle 7	Cycle 8
▸ Separation Matrix								
▸ Stacking Matrix								
▸ Sample								
Separation Time (min)	40.0	40.0	40.0	40.0	40.0	40.0	40.0	40.0
▸ Separation Voltage (volts)	250	250	250	250	250	250	250	250
▸ Matrix Removal								
▸ Antibody Diluent Time (min)	23.0	23.0	23.0	23.0	23.0	23.0	23.0	23.0
▸ Primary Antibody Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
▸ Secondary Antibody Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
▸ Quaternary Antibody Time (min)	120.0	120.0	120.0	120.0	120.0	120.0	120.0	120.0
▸ Detection								

Each row contains an assay protocol step. Each step can contain a reagent assay plate row assignment and one or more parameter settings. To view the details for a step, click on the white arrow next to the step name. We recommend using the default protocol settings for the assays. An expanded list of the default protocol step parameters is shown here for a protocol with split Running Buffer:

Split Running Buffer	Cycle 1	Cycle 2	Cycle 3	Cycle 4	Cycle 5	Cycle 6	Cycle 7	Cycle 8
▾ Separation Matrix								
Well Row	N7							
Load Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
▾ Stacking Matrix								
Well Row	O7							
Load Time (sec)	15.0	15.0	15.0	15.0	15.0	15.0	15.0	15.0
▾ Sample								
Well Row	A1							
Load Time (sec)	9.0	9.0	9.0	9.0	9.0	9.0	9.0	9.0
Separation Time (min)	40.0	40.0	40.0	40.0	40.0	40.0	40.0	40.0
▾ Separation Voltage (volts)	250	250	250	250	250	250	250	250
Standards Exposure (sec)	5.0	5.0	5.0	5.0	5.0	5.0	5.0	5.0
EE Immobilization Time (sec)	250.0	250.0	250.0	250.0	250.0	250.0	250.0	250.0
▾ Matrix Removal								
Matrix Removal Time (sec)	230.0	230.0	230.0	230.0	230.0	230.0	230.0	230.0
Matrix Washes	3	3	3	3	3	3	3	3
Matrix Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
▾ Antibody Diluent Time (min)	23.0	23.0	23.0	23.0	23.0	23.0	23.0	23.0
Well Row	B1							
Washes	2	2	2	2	2	2	2	2
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
▾ Primary Antibody Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Well Row	C1							
Washes	2	2	2	2	2	2	2	2
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
▾ Secondary Antibody Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Well Row	D1							
Washes	2	2	2	2	2	2	2	2
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
▾ Detection								
Well Row	J1							
Detection Profile	7 Exposures							

- You can change the primary or secondary antibody incubation time. Click the cell in the value column next to Primary Ab Time (min) or Secondary Ab Time (min) and enter a new value in minutes:


Split Running Buffer	Cycle 1	Cycle 2	Cycle 3	Cycle 4	Cycle 5	Cycle 6	Cycle 7	Cycle 8
▸ Separation Matrix								
▸ Stacking Matrix								
▸ Sample								
Separation Time (min)	40.0	40.0	40.0	40.0	40.0	40.0	40.0	40.0
▸ Separation Voltage (volts)	250	250	250	250	250	250	250	250
▸ Matrix Removal								
▾ Antibody Diluent Time (min)	23.0	23.0	23.0	23.0	23.0	23.0	23.0	23.0
Well Row	B1							
Washes	2	2	2	2	2	2	2	2
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
▾ Primary Antibody Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Well Row	C1							
Washes	2	2	2	2	2	2	2	2
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
▾ Secondary Antibody Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Well Row	D1							
Washes	2	2	2	2	2	2	2	2
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
▸ Detection								

3. You can change the primary and secondary antibody incubation reagent row location. Click the cell in the value column next to Well Row and select a different row on the assay plate:

	Cycle 1	Cycle 2	Cycle 3	Cycle 4	Cycle 5	Cycle 6	Cycle 7	Cycle 8
Split Running Buffer								
Separation Matrix								
Stacking Matrix								
Sample								
Separation Time (min)	40.0	40.0	40.0	40.0	40.0	40.0	40.0	40.0
Separation Voltage (volts)	250	250	250	250	250	250	250	250
Matrix Removal								
Antibody Diluent Time (min)	23.0	23.0	23.0	23.0	23.0	23.0	23.0	23.0
Well Row	B1							
Washes	2	2	2	2	2	2	2	2
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
Primary Antibody Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Well Row	C1							
Washes	2	2	2	2	2	2	2	2
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
Secondary Antibody Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Well Row	D1							
Washes	2	2	2	2	2	2	2	2
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
Detection								

NOTE: Only rows you've designated as primary antibody in the Layout tab can be selected in the Well Row drop-down menu.

4. Additional exposures can be collected in the assay if desired. To do this, click the white arrow next to Detection to expand the row. Click the cell in the exposure column next to Detection Profile to open the Detection Profile window. Additional times can be added to the protocol by clicking the **Add** button, entering the values and selecting **OK**.



- You can modify any other protocol parameters as needed.


NOTES:

For more information on changing protocol step parameters other than incubation times, contact ProteinSimple Technical Support or your local Field Application Specialist.

When a protocol is being edited, an asterisk will appear in the Protocol tab to indicate changes have been made and should be saved.

Step 4 - Add Assay Notes (Optional)


- Click on the Notes tab.
- Click in the pane and type any assay or protocol notes. This information will be stored with the assay and run data.



NOTE: When notes are being edited, an asterisk will appear in the Notes tab to indicate changes have been made and should be saved.

Step 5 - Select a Schedule (Optional for Sally Sue and Peggy Sue)

Sally Sue and Peggy Sue can execute cycles serially or in parallel. To choose a schedule option, select **Edit** and click **Schedule**.



- **Serial** - Each cycle is executed in sequential order, and each cycle completes before the next one starts. Selecting this option result in the longest run times.
- **Overlapping** - Cycle steps are executed in parallel if timing permits shared use of the separation tray and robot. Each cycle will run through to completion once started. This option offers a reduction in overall run time without affecting the individual run time for each cycle.
- **Overlapping with Hold (default)** - Each cycle executes through the separate step in sequence and then enters a variable length hold step where capillaries remain in the incubation tray. More sample stability is provided with this option as proteins are 'locked' to the capillary wall and therefore more stable than samples on the assay plate. This option also results in the most significant reduction in overall run time.

*NOTE: We recommend using the **Overlapping with hold** option. For more information on when to use other schedule options, please contact ProteinSimple Technical Support or your local Field Application Specialist.*

Step 6 - Add Assay Plate Annotations (Optional)

Custom reagent names, notes and annotations can be entered for individual rows and wells on the assay plate. This information will be stored with assay and run data. During post-run analysis, this information will be used to apply custom analysis settings to specific sample names, blocking reagent names, primary antibody names or sample and primary antibody attributes. This will be explained in more detail in "Compass Analysis Settings Overview" on page 232.

NOTE: Template pane information can also be added or updated after a run is complete.


To enter annotations:

1. Click on the Template tab. The default annotations for reagent rows and individual wells on the assay plate will display:


	1	2	3	4	5	6	7	8	9	10	11	12
A	Biot. Ladder						Sample					
B							Antibody Diluent					
C	Blocking						Primary Antibody					
D	Streptavidin ...						Secondary Antibody					
J							Detection					

2. Change or add row and well annotations as needed. To do this:
 - a. **To enter annotations for a specific well** - Right click the well and select **Edit** or click **Edit** in the


upper right corner of the pane or double click the selected well. The following box will display:



Enter the reagent name, notes or specific attributes (for example, sample concentration or dilution ratio):




Click **OK**. The new information will display in the selected well:



- b. **To enter annotations for multiple wells or a row** - To select individual wells, press and hold the Control key and then select individual wells. To select a sequential set of wells or a full row, press and


hold the Shift key, then select the first well and last well. Next, right click and select **Edit** or click **Edit** in the upper right corner of the pane. The following box will display:



Enter the reagent name, notes and specific attributes (for example, sample concentration or dilution ratio):




Click **OK**. The new information will display in the selected wells:



*NOTE: When annotations are being edited, an asterisk will appear in the **Template** tab to indicate changes have been made and should be saved. Only Name and Attribute will be used by Compass to annotate the data.*

Step 7 - Save the Assay

1. Select **File** from the main menu and click **Save As**. Enter the assay name and click **Save**.




NOTE: New assays are saved in the Compass Assays directory.

Step 8 - Modify Default Analysis Parameters (Optional)

You can preset the parameters used to analyze run data generated with the assay.

1. Select **Edit** from the main menu and click **Default Analysis**. The following screen will display:




2. We recommend using the default parameters for Simple Western assays. However, you can modify any you want, then click **OK**. For detailed information on analysis parameters, please refer to "Compass Analysis Settings Overview" on page 232.

Total Protein Assays: Creating a New Assay



To create a new assay, we recommend using one of the template assays and making any necessary modifications from there.

Step 1 - Open a Template Assay


1. Select **File** in the main menu and click **New Assay**:



2. Select the template assay for your instrument by selecting **Total Protein Size** as your assay type, **Size Range** and **Cartridge Type** if running Wes, or choose **Open Assay** to select from the menu of saved assays.



If you're running a new Total Protein assay with split Running Buffer on Sally Sue or Peggy Sue, the window has a Split Running Buffer check box to indicate the use of split buffer. This box is checked by default.



NOTES:

Compass v2.7 or higher and the latest embedded are needed to run the split Running Buffer assays.

If you make changes to the assay, such as defining the wells in the plate, be sure to select **File** and **Save** before proceeding.

Split Running Buffer will be noted above the protocol steps and plate layout if that assay was selected for Sally Sue or Peggy Sue.

The screenshot shows the Compass software interface for a 'Sally Sue Total Protein High MW Split Running Buffer' assay. The top part displays a protocol table with columns for Cycle 1 through Cycle 8. The bottom part shows a plate layout grid with 12 columns and 10 rows (A-J).

	Cycle 1	Cycle 2	Cycle 3	Cycle 4	Cycle 5	Cycle 6	Cycle 7	Cycle 8
Split Running Buffer								
Separation Matrix								
Stacking Matrix								
Sample								
Separation Time (min)	50.0	50.0	50.0	50.0	50.0	50.0	50.0	50.0
Separation Voltage (volts)	275	275	275	275	275	275	275	275
Matrix Removal								
Biotin Labeling Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Antibody Diluent Time (min)	26.0	26.0	26.0	26.0	26.0	26.0	26.0	26.0
Total Protein HRP Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Detection								


	1	2	3	4	5	6	7	8	9	10	11	12
A	Biot. Ladder						Sample					
B	Antibody Diluent						Total Protein Biotin Labeling Reagent					
C							Antibody Diluent					
D							Total Protein Streptavidin HRP					
J							Detection					

Step 2 - Assign Assay Plate Reagents (Optional)

1. Click on the Layout tab. Displayed are the default row locations where each reagent should be placed on the assay plate.


NOTE: For Sally Sue and Peggy Sue, up to 96 different samples or conditions (lysates, proteins, antibodies or incubation times) can be assayed using the full eight cycles to accommodate experimental needs.

Wes plate layout:



- **Row A** - Biotinylated Ladder (A1: 12-230 kDa or 66-440 kDa) and Sample (A2-A25)
- **Row B** - Antibody Diluent (B1) and Labeling Reagent (B2-25)
- **Row C** - Antibody Diluent
- **Row D** - Total Protein Streptavidin-HRP
- **Row E** - Luminol-S/Peroxide mix
- **Row G** - Wash buffer
- **Row H** - Wash buffer
- **Row I** - Wash buffer
- **Row M-P** - Pre-filled Separation Reagents


Sally Sue/Peggy Sue plate layout:




- **Row A** - Biotinylated Ladder (A1: 12-230 kDa or 66-440 kDa) and Sample (A2-A12)
- **Row B** - Antibody Diluent (B1) and Total Protein Labeling Reagent (B2-12)
- **Row C** - Antibody Diluent
- **Row D** - Total Protein Streptavidin-HRP
- **Row J** - Luminol-S/Peroxide mix
- **Row M** - Water (M5-M20)
- **Row N** - Water (N5-N6 and N19-N20) and Separation Matrix (N7-N18)
- **Row O** - Water (O5-O6 and O19-O20) and Stacking Matrix (O7-O18)
- **Row P** - Water (P5-P20)

NOTE: For details on sample, reagent and assay plate preparation, please refer to the product insert provided with the Simple Western kits.


- If needed, well assignments can be modified. Any row assignments changed in the **Layout** pane are updated in the Protocol pane automatically.
 - **To move a reagent row to another location** - Click the row in the **Layout** pane, then drag and drop it on the new location. The row from which it was moved will be reassigned as empty.




- **To insert a sample row** - Click an empty row or the row below where the new sample row should be inserted, then click the **S** icon (Insert a sample row) in the **Layout** pane toolbar. A new sample row will be added in the empty row or inserted above the selected row.



- **To insert a detection row** - Click an empty row or the row below where the new detection row should be inserted, then click the **D** icon (Insert a detection row) in the **Layout** pane toolbar. A new detection row will be added in the empty row or inserted above the selected row.



- **To insert a fourth or fifth incubation reagent** - Click an empty row or the row below where the new incubation reagent should be inserted, then click the **4** icon (Insert a quaternary row) in the **Layout** pane toolbar. A new incubation reagent row will be added in the empty row or inserted above the selected row.




A fifth incubation reagent can now be added by repeating the above and clicking the **5** icon (Insert a Quinary row) instead.

NOTES:

When using the fourth and fifth incubation steps, only five cycles can typically be performed per run as more water is needed per cycle for the additional wash steps.

Row J is the last row assignment that can be used on the assay plate. New, sample, incubation or luminol rows cannot be inserted if row J already has an assigned reagent. To insert a row in this case, you must first move the contents of row J to another area on the plate.

- **To delete a row** - Click the row to be deleted, then click the red **X** icon (Delete a row) in the toolbar. Only empty rows and fourth incubation rows can be deleted. Rows for required assay reagents cannot be deleted.



NOTES:

We recommend keeping two rows empty between the substrate and the closest dispensed HRP conjugate to avoid any interaction with the substrate.

Rows K and L are purposely left empty, don't put anything into these wells.

Separation and Stacking matrices can only be pipetted into the designated wells.

Water is used in the wells to reduce the potential of evaporation of the matrices.

Step 3 - Modifying the Assay Protocol (Optional for All Instruments)

1. Click on the Protocol tab. This pane shows the individual steps of the assay protocol and allows you to change parameters. When creating a new assay, the default protocol automatically assigns all reagent locations for the assay. The Sally Sue assay is shown as an example:

	Cycle 1	Cycle 2	Cycle 3	Cycle 4	Cycle 5	Cycle 6	Cycle 7	Cycle 8
Split Running Buffer								
Separation Matrix								
Stacking Matrix								
Sample								
Separation Time (min)	40.0	40.0	40.0	40.0	40.0	40.0	40.0	40.0
Separation Voltage (volts)	250	250	250	250	250	250	250	250
Matrix Removal								
Biotin Labeling Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Antibody Diluent Time (min)	26.0	26.0	26.0	26.0	26.0	26.0	26.0	26.0
Total Protein HRP Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Detection								

Each row contains an assay protocol step. Each step can contain a reagent assay plate row assignment and one or more parameter settings. To view the details for a step, click on the white arrow next to the step name. We recommend using the default protocol settings for the assays. An expanded list of the default protocol step parameters is shown:

Protocol History Notes								
Split Running Buffer	Cycle 1	Cycle 2	Cycle 3	Cycle 4	Cycle 5	Cycle 6	Cycle 7	Cycle 8
Stacking Matrix								
Well Row	O7							
Load Time (sec)	15.0	15.0	15.0	15.0	15.0	15.0	15.0	15.0
Sample								
Well Row	A1							
Load Time (sec)	9.0	9.0	9.0	9.0	9.0	9.0	9.0	9.0
Separation Time (min)	40.0	40.0	40.0	40.0	40.0	40.0	40.0	40.0
Separation Voltage (volts)	250	250	250	250	250	250	250	250
Standards Exposure (sec)	5.0	5.0	5.0	5.0	5.0	5.0	5.0	5.0
EE Immobilization Time (sec)	250.0	250.0	250.0	250.0	250.0	250.0	250.0	250.0
Matrix Removal								
Matrix Removal Time (sec)	230.0	230.0	230.0	230.0	230.0	230.0	230.0	230.0
Matrix Washes	3	3	3	3	3	3	3	3
Matrix Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
Biotin Labeling Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Well Row	B1							
Washes	0	0	0	0	0	0	0	0
Wash Soak Time (sec)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Antibody Diluent Time (min)	26.0	26.0	26.0	26.0	26.0	26.0	26.0	26.0
Well Row	C1							
Washes	2	2	2	2	2	2	2	2
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
Total Protein HRP Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Well Row	D1							
Washes	2	2	2	2	2	2	2	2
Wash Soak Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
Detection								

- You can change the biotin labeling and Total Protein streptavidin-HRP incubation time. Click the cell in the value column next to Biotin Labeling Time (min) or Total Protein HRP Time (min) and enter a new value in minutes:


Protocol History Notes								
Split Running Buffer	Cycle 1	Cycle 2	Cycle 3	Cycle 4	Cycle 5	Cycle 6	Cycle 7	Cycle 8
Separation Matrix								
Well Row	N7							
Load Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
Stacking Matrix								
Sample								
Separation Time (min)	40.0	40.0	40.0	40.0	40.0	40.0	40.0	40.0
Separation Voltage (volts)	250	250	250	250	250	250	250	250
Matrix Removal								
Biotin Labeling Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Well Row	B1							
Washes	0	0	0	0	0	0	0	0
Wash Soak Time (sec)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Antibody Diluent Time (min)	26.0	26.0	26.0	26.0	26.0	26.0	26.0	26.0
Total Protein HRP Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Detection								

3. You can also change the biotin labeling and Total Protein streptavidin-HRP row location. Click the cell in the value column next to Well Row and select a different row on the assay plate:

Split Running Buffer	Cycle 1	Cycle 2	Cycle 3	Cycle 4	Cycle 5	Cycle 6	Cycle 7	Cycle 8
▲ Separation Matrix								
Well Row	N7							
Load Time (sec)	150.0	150.0	150.0	150.0	150.0	150.0	150.0	150.0
▷ Stacking Matrix								
▷ Sample								
Separation Time (min)	40.0	40.0	40.0	40.0	40.0	40.0	40.0	40.0
▷ Separation Voltage (volts)	250	250	250	250	250	250	250	250
▷ Matrix Removal								
▲ Biotin Labeling Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
Well Row	B1							
Washes	0	0	0	0	0	0	0	0
Wash Soak Time (sec)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
▷ Antibody Diluent Time (min)	26.0	26.0	26.0	26.0	26.0	26.0	26.0	26.0
▷ Total Protein HRP Time (min)	30.0	30.0	30.0	30.0	30.0	30.0	30.0	30.0
▷ Detection								

NOTE: Only rows you've designated as labeling reagent in the Layout tab can be selected in the Well Row drop-down menu.

4. Additional exposures can be collected in the assay if desired. To do this, click the white arrow next to Detection to expand the row. Click the cell in the exposure column next to Detection Profile to open the Detection Profile window. Additional times can be added to the protocol by clicking the **Add** button, entering the values and selecting **OK**.



5. You can modify any other protocol parameters as needed.

NOTES:

For more information on changing protocol step parameters other than incubation times, contact ProteinSimple Technical Support or your local Field Application Specialist.


When a protocol is being edited, an asterisk will appear in the Protocol tab to indicate changes have been made and should be saved.

Steps 4 - 8

Steps 4 through 8 for creating a Total Protein Assay are the same as when you're creating an Immunoassay. Please go to "Step 4 - Add Assay Notes (Optional)," starting on page 34 to continue.

Making Changes to an Existing Assay

1. Select **File** in the main menu and click **Open Assay**.




2. A list of the last assays opened will display. Select one of these assays or click **Browse** to open the assay folder and select a different assay.
3. Follow the steps in "Immunoassays: Creating a New Assay" on page 25 to make changes and save the assay.

Switching Between Open Assays

If more than one run file is open, you can switch between viewing the assays for each run. To do this:

1. Click the down arrow in the run box.




2. Select the run for the assay you want to view from the drop down list.


Creating a Template Assay

You can create and save template assays that can be used as a starting point for creating new assays. To do this:

1. Select **File** in the main menu. Click **Open Assay** to open an existing assay or click **New Assay** to open an existing template assay.
2. Follow the steps in “Immunoassays: Creating a New Assay” on page 25 to make changes to the assay.
3. When changes are complete, select **File** in the main menu and click **Save As**. Select the New Assays folder:



4. Type the name for the new template assay and click **Save**.
5. Select **File** in the main menu and click **New Assay**. The new template assay will now be available in the drop down list:



Viewing and Changing the Detection Exposures


To view the current detection profile, roll the cursor over the exposures cell in Protocol pane:

▲ Total Protein HRP Time (min)	30.0	30.0	30.0	30.0
Well Row	D1	D1	D1	D1
Washes	2	2	2	2
Wash Soak Time (sec)	150.0	150.0	150.0	150.0
▲ Detection				
Well Row	J1	J1	J1	J1
Detection Profile	7 Expos	7 Expos	7 Exposures	7 Exposures

signal, 5.0 sec
 signal, 15.0 sec
 signal, 30.0 sec
 signal, 60.0 sec
 signal, 120.0 sec
 signal, 240.0 sec
 signal, 480.0 sec


While we recommend using the default assay detection profile, you can change the profile, exposure times and exposure sequences if needed. To do this:

1. Select the exposures cell in the Protocol pane and click the button, or double-click in the cell. The following screen will display:



Each row represents an individual exposure that will be taken during the run.

- a. **To change an existing exposure time** - Click in the exposure cell and enter a new time in seconds:



- b. **To delete an existing exposure** - Select a type or exposure cell and click **Remove**.
 - c. **To add a new exposure** - Select **Add**. A new exposure will be added to the end of the list. Click in the exposure cell and enter an exposure time in seconds.
2. Click **OK** to save and exit.

Copying Protocols and Templates

The steps and parameters in the Protocol pane can be copied and d into other documents, as can the graphic image of the annotations in the Template pane.

Copying an Assay Protocol

1. Click on the Protocol tab.
2. Select **Edit** in the main menu and click **Copy**.
3. Open a document (Microsoft® Word®, Excel®, etc.). Right click in the document and select **Paste**. All assay protocol steps and parameters will be copied into the document as a list in the same format shown in the Protocol pane.

Copying an Assay Template


1. Click on the Template tab.
2. Select **Edit** in the main menu and click **Copy**.
3. Open a document (Microsoft® Word®, Excel®, etc.). Right click in the document and select **Paste**. The assay template will be copied into the document as an image, exactly as it is shown in the Template pane.

Printing Protocols and Templates

The information in the Protocol pane can be printed, as can a graphic image of the information in the Template pane.

Printing an Assay Protocol


1. Click on the Protocol tab.
2. Select **File** in the main menu, click **Print**, and then click **Print Protocol**.



All assay protocol steps and parameters will be printed as a list in the same format shown in the Protocol pane.

Printing an Assay Template

1. Click on the Template tab.
2. Select **File** in the main menu, click **Print**, and then click **Print Template**.



The assay template will print as an image, exactly as it is shown in the Template pane.

Importing and Exporting Protocols and Templates

The assay protocol in an open assay or run file can be exported as a separate file, as can the assay template annotation information. This allows the same assay protocol and template information to be imported into another assay at a later time, rather than having to re-enter assay protocol and template information manually.

Importing an Assay Protocol


NOTE: Importing an assay protocol imports information into the Protocol pane only.

1. Open the assay you want to import the assay protocol in to.
2. Select **File** in the main menu and click **Import Protocol**.
3. Select a protocol file (*.protocol) and click **OK**. The imported information will display in the Protocol pane.

Exporting an Assay Protocol

NOTE: Exporting an assay protocol exports information in the Protocol pane only.

1. Open the assay you want to export the assay protocol from.
2. Select **File** in the main menu and click **Export Protocol**. The following window displays:



3. The default directory is Compass/Assays. Change the directory if needed.
4. Enter a protocol name and click **Save**. The protocol will be saved as a *.protocol file.

Template Export and Import

The Template information including the sample and reagent names can be exported to a file. As an example, a default Wes Template is displayed in the Assay view.

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
A	Bio...																								
B																									
C	Blo...																								
D	Str...																								
E																									


There are three file format options:

- A comma separated CSV file that is best opened in a spreadsheet.
- An XML template file that has full control over the sample and reagent names and layout.
- A tab delimited TXT file that can be opened in a spreadsheet even when the decimal separator is set to a comma.


Exporting the Template to a CSV file

Follow the steps below to easily import a template into Compass, populating it with reagent names (sample, antibody, etc.) and their corresponding attributes (concentration, dilution factor, etc.). This example uses a Wes Template, but the same steps can be followed for a Sally Sue or Peggy Sue plate template.

1. Select **File** in the main menu and click **Export Template**. The following window displays:



Make sure the Save as type is set to CSV.



2. Enter a protocol name and click **Save**. The protocol will be saved as a *.csv file.
3. Open the CSV file in a spreadsheet program like Microsoft® Excel®.

	A	B	C	D	E	F
1	Biot. Ladder	Sample	Sample	Sample	Sample	Sample
2						
3	Antibody Diluent	Antibody Diluent	Antibody Diluent	Antibody Diluent	Antibody Diluent	Antibody Diluent
4						
5	Blocking	Primary Antibody				
6						
7	Streptavidin HRP	Secondary Antibody HRP				
8						
9	Luminol/Peroxide	Luminol/Peroxide	Luminol/Peroxide	Luminol/Peroxide	Luminol/Peroxide	Luminol/Peroxide

The names in the spreadsheet are arranged in the same order as the Compass Template, and the rows alternate between names and attributes.

NOTE: The default assay has no attributes so these rows will be empty.

4. Edit the names and add attributes, then save the spreadsheet as a CSV file.

NOTE: Make sure not to edit the first column of the spreadsheet, this corresponds to the ladder wells.

5. To import the edited CSV file into Compass, select **File** in the main menu, click **Import Template** and then browse to the .csv file you just saved. Once imported, the edited CSV file displays the edited Sample names and Primary Antibody with attributes in the Template.

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
A	Bio...	Patient A		Patient B		Patient C		Patient D		Patient E		Patient F		Patient G		Patient H									
B	Antibody Diluent																								
C	Bio...	1mg/mL	10mg/L	40mg/L	1mg/mL	10mg/L	40mg/L	1mg/mL	10mg/L	40mg/L	1mg/mL	10mg/L	40mg/L	1mg/mL	10mg/L	40mg/L	1mg/mL	10mg/L	40mg/L	1mg/mL	10mg/L	40mg/L	1mg/mL	10mg/L	40mg/L
D	Str...	Secondary Antibody HRP																							
E	Luminol/Peroxide																								

Template Cut and Paste

The names and attributes in the Template can be copied and pasted within Compass, between two copies of Compass, and between Compass and a spreadsheet like Microsoft® Excel®.

All you need to do is enter names in a spreadsheet row, then copy to the clipboard. Next, select a well in the Template and paste from the clipboard. The new names will be pasted into the row at the selected column and columns to the right. Here’s a quick example of how to do it:

1. Copy these three rows from a spreadsheet:

	A	B	C
1	S1	S2	S3

2. Select a well in the Template, A2 is selected for this example:

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	
A	Bio...																									
B																										
C	Blo...																									
D	Str...																									
E																										

3. Paste the names from the clipboard.

Note: Make sure not to select any of the ladder wells (A1, B1, C1, etc.)

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	
A	Bio...	S1	S2	S3																						
B																										
C	Blo...																									
D	Str...																									
E																										

4. To include attributes, add them in the next row:

		A	B	C
1	S1	S2	S3	
2		100	200	300

5. Pasting into the Template will update the names and attributes.

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	
A	Bio...	S1	S2	S3																						
B																										
C	Blo...																									
D	Str...																									
E																										

Multiple rows can also be copied and pasted with names and attributes on alternate rows:

	A	B	C
1	S1	S2	S3
2	100	200	300
3	AD1	AD2	AD3
4	101	201	301
5	AB1	AB2	AB3
6	102	202	302
7	HRP1	HRP2	HRP3
8	103	203	303
9	LUM1	LUM2	LUM3
10	104	204	304

Only wells in the same row can be selected, but multiple rows will be pasted into the rows below.

Note: Make sure not to select any of the ladder wells (A1, B1, C1, etc.)

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
A	Bio...	S1 100	S2 200	S3 300	Sample																				
B	Ant...	AD1 101	AD2 201	AD3 301	Antibody Diluent																				
C	Blo...	AB1 102	AB2 202	AB3 302	Primary Antibody																				
D	Str...	HRP1 103	HRP2 203	HRP3 303	Secondary Antibody HRP																				
E	Lum...	LUM1 104	LUM2 204	LUM3 304	Luminol/Peroxide																				

Chapter 3:

Running a Size Assay

Chapter Overview

- Starting a Run
- Stopping a Run

Starting a Run

Step 1 - Get Ready

1. Create or open the an assay file in Compass.
2. Prepare the assay plate using the information provided in the product insert.

IMPORTANT

To prevent well evaporation and ensure best results, keep a lid on the assay plate until ready to use.

3. Prepare the instrument following the procedure described in the User Guide.

IMPORTANT


Capillaries are light sensitive. Keep them covered until you are ready to transfer the capillary cartridge or box to the instrument.

4. Place assay plate into the sample tray of the instrument.

Step 2 - Start the Run

You can start a run in one of two ways depending on whether you want to run an assay using existing parameters or set up a new assay.


1. New run of an existing assay:
 - a. Select **File** in the main menu and click **Open Assay**.



- b. A list of the last five assays opened will display. Select one of these assays or click **Browse** to open the Assay folder and select a different assay.
2. Alternatively, choose **New Assay** and select one of the size assays to get the default assay conditions for using Wes, Sally Sue or Peggy Sue.
 - a. The **Start** button will display. This indicates that an assay has been loaded.




- a. Go to the Assay screen and verify this is the assay you want to use. If not, select **File** in the main menu, click **Open Assay**, and select another assay.
3. Click **Start** to begin the run.
 - a. **For Wes:** Instrument status will change to **Running**, the **Start** button will change to **Stop** and the **Run Progress Bar** will be displayed:



- b. **For Sally Sue/Peggy Sue:** This will launch the **Start Run Wizard** (continue on to the next step).


*NOTE: If you didn't clean the manifold in Sally Sue or Peggy Sue prior to starting the run, a message will appear. If you see this message, click **Yes** to cancel the run and perform the manifold cleaning.*

4. **Check Water and Waste.** The fluid levels in the accessory module bottles will be checked by the software. If the levels in both bottles will allow Sally Sue an Peggy Sue to complete the run, the wizard screen will display **Water Level OK** and **Waste Level OK** messages. Click **Next** to proceed.





*NOTE: If the waste level is too high or the water level is too low to complete the run, messages to indicate one or both will be presented in this screen. If this occurs, fill or empty each bottle as indicated using the procedures outlined in the Sally Sue/Peggy Sue User Guide. When this is complete, the error status will be automatically updated and allow the **Start Run Wizard** to proceed.*

- 5. Replace Sponge.** You should use a new sponge each time a new experimental run is started. Discard the old sponge and put a new sponge in the manifold wash station. Click **Next** to proceed.




- 6. Start.** The resource tray will automatically open. Fill the Wash Buffer, Running Buffer, and Matrix Removal Buffer cups and insert the capillary box in the primary capillary box position. Add a secondary capillary box if more capillaries are required. Click **Next** when finished. The resource tray will close.


NOTE: Make sure enough capillaries are loaded to complete the run. A capillary box should be inserted in the primary position first and if additional capillaries are needed, use the secondary position. When the primary box is depleted, Sally Sue or Peggy Sue will automatically move to the secondary box. Discard left-over Running Buffer, Wash Buffer and Matrix Removal Buffer prior to refilling bottles for new run.




NOTE: You can also refer to the labels on the resource tray for proper insertion of reagents.



7. **Load Sample Plate.** The sample tray will automatically open. Place the lidded 384-well plate on the cooling block in the tray, ensuring that the A1 well position is aligned with the upper left corner of the cold block. When this is complete, click **Next**. The sample tray will close.



A1 Position




NOTES:


Plate lids be used on sample plates for Sally Sue and Peggy Sue. If a lid is not detected, a message will be displayed in the Start Run Wizard. Compass software will reopen the sample tray to allow you to insert a lid.

When inserting the sample plate, ensure that the plate is firmly seated and level on the cold block. Plates that are not level can interfere with the movement of the sample tray.

8. **Begin the Automated Run.** The data file name will automatically default to the assay name appended with the current date and time. To change the file name, begin typing in the text box. To change the directory where the data file will be stored, click **Browse**:



Click **Start** to begin the run. Instrument status will change to **Running**, the **Start** button will change to **Stop** and the **Run Progress bar** will display:



The run will continue until complete (~14-19 hours depending on the assay).

Step 3 - Post-Run Procedures

1. Remove the capillary cartridge (Wes) or empty the capillary discard tray (Sally Sue/Peggy Sue).
2. Remove the assay plate.
3. Dispose of the assay plate and capillaries. Disposal will depend on the samples that have been assayed. If sample origins are unknown, we recommend that used capillaries cartridges, capillaries and plates be disposed of in biohazard waste.

!WARNING! SHARPS HAZARD

The capillaries may present a potential sharps hazard. Dispose of used capillaries in biomedical waste sharps containers.




!WARNING! BIOHAZARD

Samples and waste bottle contents should be handled by procedures recommended in the CDC/NIH manual: Biosafety in Microbiological and Biomedical Laboratories (BMBL). The manual is available from the U.S. Government Printing Office or online at <http://www.cdc.gov/biosafety/publications/bmbl5/>.

Depending on the samples used, waste bottle contents may constitute a biohazard. Use precaution when emptying the waste bottle. Dispose of waste bottle contents in accordance with good laboratory practices and local, state/provincial, or national environmental and health regulations. Read and understand the Material Safety Data Sheets (MSDSs) provided by the manufacturers of the chemicals in the waste vial before you store, handle, or dispose of chemical waste.

Stopping a Run

1. To stop a run, click **Stop**.
 - a. **Wes:** When the run stops, Wes's status will go to **Not Ready** and he'll start a process that plugs the capillaries in the cartridge. Once that's done, you can then remove the capillary cartridge and plate and discard them.
 - b. **Sally Sue/Peggy Sue:** When the run stops, Sally Sue's or Peggy Sue's status will go to **Not Ready** and a **Cleanup** button displays:



NOTE: If a run is stopped prior to completion, Sally Sue and Peggy Sue must perform a cleaning protocol to remove used capillaries from system trays. This is required to prepare the system for the next run.

2. Click **Cleanup**.



Allow Sally Sue and Peggy Sue to complete the cleaning protocol which takes about ten minutes. When complete, instrument status will change to **Ready** and a new run can be started.

Chapter 4:

Charge Assays

Chapter Overview

- Assay Screen Overview
- Reagent Color Coding
- Opening an Assay
- Creating a New Assay
- Making Changes to an Existing Assay
- Switching Between Open Assays
- Creating a Template Assay
- Viewing and Changing the Detection Exposures
- Copying Protocols and Templates
- Printing Protocols and Templates
- Importing and Exporting Protocols and Templates

Assay Screen Overview

The Assay screen is used to create, view, and edit assays. To access this screen, click **Assay** in the screen tab:



Assay Screen Panes

The Assay screen has four panes:

- **Layout** - Displays a map of the assay plate and shows where assay reagents will be located.
- **Protocol** - Lists individual assay protocol steps and parameters that Peggy Sue will execute for each of the 12 capillaries simultaneously.
- **Notes** - Lets you enter specific assay information that is saved with the assay and can be used for future reference.
- **Template** - Enter annotations for the individual well and row reagents in the assay plate.

The screenshot shows the 'Assay: Peggy Sue Charge' window. The 'Layout' pane displays a 12x12 grid representing the assay plate. The 'Protocol' pane shows a table with columns for Cycle 1 through Cycle 8 and rows for Sample, Separation, Immobilization Time (sec), Primary Antibody Time (min), Secondary Antibody Time (min), and Detection. The 'Template' pane shows a 12x4 grid with rows labeled A, B, C, and J, and columns labeled 1 through 12. The grid is color-coded: Row A is orange (Sample), Row B is light blue (Primary Antibody), Row C is teal (Secondary Antibody), and Row J is yellow (Luminol/Peroxide).

	Cycle 1	Cycle 2	Cycle 3	Cycle 4	Cycle 5	Cycle 6	Cycle 7	Cycle 8
Sample								
Separation								
Immobilization Time (sec)	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0
Primary Antibody Time (min)	120.0	120.0	120.0	120.0	120.0	120.0	120.0	120.0
Secondary Antibody Time (min)	60.0	60.0	60.0	60.0	60.0	60.0	60.0	60.0
Detection								

	1	2	3	4	5	6	7	8	9	10	11	12
A	Sample											
B	Primary Antibody											
C	Secondary Antibody											
J	Luminol/Peroxide											

Software Menus Active in the Assay Screen


The following software menus are available:

- File
- Edit
- Instrument (when Compass is connected to Peggy Sue)
- Window
- Help

The File and Edit menu options specific to the Assay screen are described next.

File Menu


The following File menu options are active:



- **New Assay** - Creates a new assay from a starter template.
- **Open Assay** - Opens an existing assay.
- **Save** - Saves the open assay.
- **Save As** - Saves the open assay under a different file name.
- **Import Protocol** - Imports a saved protocol file into an assay.
- **Import Template** - Imports a saved template file into an assay.
- **Export Protocol** - Exports the current protocol file for future use.
- **Export Template** - Exports the current template file for future use.
- **Print** - Prints the information in the Protocol or Template panes.
- **Exit** - Closes Compass.

Edit Menu


The following Edit menu options are active:



- **Copy** - Copies the information in the Protocol or Template panes into other documents.
- **Default Analysis** - Displays the default settings that will be used to analyze the run data generated with an assay.
- **Analysis** - Not active in this screen.
- **Preferences** - Set and save custom preferences for data export, plot colors in the graph and Peggy Sue's Twitter settings. See Chapter 10, "Setting Your Preferences" for more information.

Reagent Color Coding

The Assay screen uses color coding to identify various assay reagents in all panes. The Layout pane is shown in the following example:




- **Orange** - Samples
- **Light Blue** - Primary antibody
- **Teal** - Tertiary antibody
- **Yellow** - Secondary HRP conjugate
- **Gold** - Luminol/Peroxide mix

Opening an Assay

To open an existing assay:

1. Select **File** in the main menu and click **Open Assay**.




2. A list of the last five assays opened will display. Select one of these assays or click **Browse** to open the Assay folder and select a different assay.

Creating a New Assay


To create a new assay, we recommend using the Peggy Sue template assay and modifying from there as needed.

Step 1 - Open a Template Assay

1. Select **File** in the main menu and click **New Assay**:



2. A list of template assays that can be used as a starting point for new assays will display. Click **Peggy Sue Charge**. The Peggy Sue template assay and default settings will display in the Assay screen:




Step 2 - Assign Assay Plate Reagents (Optional)

1. Click on the Layout tab. This pane displays the default row locations of where each reagent should be placed on the assay plate.

NOTE: Up to 96 different samples or conditions (lysates, proteins, antibodies or incubation times) can be assayed using the full eight cycles to accommodate experimental needs.


- **Row A** - Samples




- **Row B** - Primary Antibody
- **Row C** - Secondary Antibody
- **Row D** - Tertiary Antibody (optional)
- **Row E** - Luminol-S/Peroxide mix

NOTE: For details on sample, reagent and assay plate preparation, please refer to the product insert provided with the Simple Western kits.


2. If needed, well assignments can be modified as described below. Any row assignments changed in the **Layout** pane are updated in the Protocol pane automatically.
 - **To move a reagent row to another location** - Click the row in the **Layout** pane, then drag and drop it on the new location. The row from which it was moved will be reassigned as empty.




- **To insert a new row** - Click the row below where the new one should be inserted, then click **Insert an empty row** (circle icon) in the **Layout** pane toolbar. A new row will be inserted above the selected row.




- **To insert a sample row** - Click an empty row or the row below where the new sample row should be inserted, then click **Insert a sample row** (S icon) in the **Layout** pane toolbar. A new sample row will be added in the empty row or inserted above the selected row.



- **To insert a detection row** - Click an empty row or the row below where the new detection row should be inserted, then click **Insert a detection row** (D icon) in the **Layout** pane toolbar. A new detection row will be added in the empty row or inserted above the selected row.



- **To insert a fourth or fifth incubation reagent** - Click an empty row or the row below where the new incubation reagent should be inserted, then click **Insert a Quaternary Row** (4 icon) in the **Layout** pane toolbar. A new incubation reagent row will be added in the empty row or inserted above the selected row.




A fifth incubation reagent can now be added by repeating the above and clicking **Insert a Quinary Row** (5 icon) instead.

NOTES:

When using the fourth and fifth incubation steps, only five cycles can typically be performed per run as more water is needed per cycle for the additional wash steps.

Row J is the last row assignment that can be used on the assay plate. New, sample, incubation or luminol rows cannot be inserted if row J already has an assigned reagent. To insert a row in this case, you must first move the contents of row J to another area on the plate.


- **To delete a row** - Click the row to be deleted, then click **Delete** (red x icon) in the toolbar. Only empty rows and fourth incubation rows can be deleted. Rows for required assay reagents cannot be deleted.



NOTE: Samples, antibodies and blocking buffer can be dispensed in Rows A-J and in columns 1-12 or 13-24. Rows K-P cannot be used for assay reagents.

Step 3 - Modifying the Assay Protocol (Optional)

1. Click on the Protocol tab. This pane displays the individual steps of the assay protocol and allows you to change parameters as needed. When creating a new assay, a default protocol will display which automatically assigns all reagent locations for Cycle 1:



The screenshot shows a software window with two tabs: "Protocol" (active) and "Notes". The "Protocol" tab contains a table with the following data:


	Cycle 1	
▶ Sample		
▶ Separation		
▶ Immobilization Time (sec)	100.0	
▶ Primary Ab Time (min)	120.0	
▶ Secondary Ab Time (min)	60.0	
▶ Tertiary Ab Time (min)	10.0	
▶ Detection		

Each row contains an assay protocol step. Each step can contain a reagent assay plate row assignment and one or more parameter settings. To view the details for a step, click on the white arrow next to the step name. We recommend using the default protocol settings for Peggy Sue assays. An expanded list of the default protocol step parameters is shown:

*Protocol		Add ▼ Remove
		Cycle 1
▲ Sample		
Well Row	A1	
Load Time (sec)	25.0	
▲ Separation		
Separation Profile	Power 1 Step	
Standards Exposure (sec)	3.0	
▲ Immobilization Time (sec)	100.0	
Washes	2	
Wash Load Time (sec)	20.0	
Wash Soak Time (sec)	150.0	
▲ Primary Ab Time (min)	120.0	
Well Row	B1	
Load Time (sec)	2.0	
Washes	2	
Wash Load Time (sec)	20.0	
Wash Soak Time (sec)	150.0	
▲ Secondary Ab Time (min)	60.0	
Well Row	C1	
Load Time (sec)	2.0	
Washes	2	
Wash Load Time (sec)	20.0	
Wash Soak Time (sec)	150.0	
▲ Tertiary Ab Time (min)	10.0	
Well Row	D1	
Load Time (sec)	2.0	
Washes	2	
Wash Load Time (sec)	20.0	
Wash Soak Time (sec)	150.0	
▲ Detection		
Well Row	F1	
Wash Load Time (sec)	2.0	
Detection Profile	5 Exposures	

- Five incubation steps are allowed per protocol. You can select the type of incubation for each step. The available incubation types and their default Simple Western use is as follows:
 - **First incubation** - Primary antibody
 - **Second incubation** - Secondary antibody
 - **Third incubation** - User defined (tertiary antibody)
 - **Fourth incubation** - User defined (quaternary antibody)
 - **Fifth incubation** - User defined (custom)

To change the type, click the incubation step name and select an option from the drop down list.



- If needed, change the primary incubation time. To do this, click the cell in the value column next to Primary Ab Time (min) and enter a new value in minutes:


Cycle 1	
Sample	
Separation	
Immobilization Time (sec)	100.0
Primary Ab Time (min)	120.0
Well Row	B1
Load Time (sec)	2.0
Washes	2
Wash Load Time (sec)	20.0
Wash Soak Time (sec)	150.0
Secondary Ab Time (min)	60.0
Tertiary Ab Time (min)	10.0
Detection	

4. If needed, change the primary incubation reagent row location. To do this, click the cell in the value column next to Well Row and select a different row on the assay plate:

Cycle 1	
Sample	
Separation	
Immobilization Time (sec)	100.0
Primary Ab Time (min)	120.0
Well Row	B1
Load Time (sec)	B1
Washes	C1
Wash Load Time (sec)	20.0
Wash Soak Time (sec)	150.0
Secondary Ab Time (min)	60.0
Tertiary Ab Time (min)	10.0
Detection	


NOTE: Only rows designated as primary antibody in the Layout tab can be selected in the Well Row drop-down menu.

5. If needed, change the secondary incubation time. To do this, click the cell in the value column next to Secondary Ab Time (min) and enter a new value in minutes:



	Cycle 1
▶ Sample	
▶ Separation	
▶ Immobilization Time (sec)	100.0
▶ Primary Ab Time (min)	120.0
▲ Secondary Ab Time (min)	60.0
Well Row	C1
Load Time (sec)	2.0
Washes	2
Wash Load Time (sec)	20.0
Wash Soak Time (sec)	150.0
▶ Tertiary Ab Time (min)	10.0
▶ Detection	


6. If needed, change the secondary incubation reagent row location. To do this, click the cell in the value column next to Well Row and select a different row on the assay plate:




	Cycle 1
Sample	
Separation	
Immobilization Time (sec)	100.0
Primary Ab Time (min)	120.0
Secondary Ab Time (min)	60.0
Well Row	C1
Load Time (sec)	2.0
Washes	2
Wash Load Time (sec)	20.0
Wash Soak Time (sec)	150.0
Tertiary Ab Time (min)	10.0
Detection	

NOTE: Only rows designated as secondary antibody in the Layout tab can be selected in the Well Row drop-down menu.

7. Additional exposures can be collected in the assay if desired. To do this, click the white arrow next to Detection to expand the row. Click the cell in the exposure column next to Detection Profile to open the Detection Profile window. Additional times can be added to the protocol by clicking the **Add** button, entering the values and selecting **OK**.




8. Modify any other protocol parameters as needed.
9. To add additional cycles to the assay protocol, click in any cell with a value in a cycle column.
 - To add one cycle, either click **Add** or click the down arrow next to Add and select **1 Cycle**
 - Select **4 Cycles** from the drop down menu to add four additional cycles
 - Select **All Cycles** from the drop down menu to add the number of cycles needed to reach the maximum of eight



New cycle columns will display using the same parameters used for cycle 1.

	Cycle 1	Cycle 2	Cycle 3	Cycle 4	Cycle 5	Cycle 6	Cycle 7	Cycle 8
▶ Sample								
▶ Separation								
▶ Immobilization Time (sec)	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0
▶ Primary Ab Time (min)	120.0	120.0	120.0	120.0	120.0	120.0	120.0	120.0
▶ Secondary Ab Time (min)	60.0	60.0	60.0	60.0	60.0	60.0	60.0	60.0
▶ Tertiary Ab Time (min)	10.0	10.0	10.0	10.0	10.0	10.0	10.0	10.0
▶ Detection								

Repeat steps 1-8 to change parameters for the added cycles. As cycles are added and reagent locations are selected, cycle number assignments will update in the Layout pane assay plate map:




NOTES:

For more information on changing protocol step parameters other than incubation times, contact ProteinSimple Technical Support.

When a protocol is being edited, an asterisk will appear in the Protocol tab to indicate changes have been made and should be saved.

Step 4 - Add Assay Notes (Optional)


1. Click on the Notes tab.
2. Click in the pane and type any assay or protocol notes. This information will be stored with the assay and run data.



NOTE: When notes are being edited, an asterisk will appear in the Notes tab to indicate changes have been made and should be saved.

Step 5 - Select a Schedule (Optional)

Peggy Sue can execute cycles serially or in parallel. To choose an option, select **Edit** and click **Schedule**.



- **Serial** - Each cycle is executed in sequential order, and each cycle completes before the next one starts. Selecting this option result in the longest run times.
- **Overlapping** - Cycle steps are executed in parallel if timing permits shared use of the separation tray and robot. Each cycle will run through to completion once started. This option offers a reduction in overall run time without affecting the individual run time for each cycle.
- **Overlapping with Hold (default)** - Each cycle executes through the separate step in sequence and then enters a variable length hold step where capillaries remain in the incubation tray. More sample stability is provided with this option as proteins are 'locked' to the capillary wall and therefore more stable than samples on the assay plate. This option also results in the most significant reduction in overall run time.

NOTE: We recommend using the overlapping with hold option. For more information on when to use other schedule options, please contact ProteinSimple Technical Support.


Step 6 - Add Assay Plate Annotations (Optional)

Custom reagent names, notes and annotations can be entered for individual rows and wells on the assay plate. This information will be stored with assay and run data. During post-run analysis, this information will be used to apply custom analysis settings to specific sample names, primary antibody names or sample and primary antibody attributes. This will be explained in more detail in “Compass Analysis Settings Overview” on page 232.


NOTE: Template pane information can also be added or updated after a run is complete.

To enter annotations:


1. Click on the **Template** tab. The default annotations for reagent rows and individual wells on the assay plate will display:




2. Change or add row and well annotations as needed. To do this:
 - a. **To enter annotations for a specific well** - Right click the well and select **Edit** or click **Edit** in the upper right corner of the pane or double click the selected well. The following box will display:




Enter the reagent name, notes or specific attributes (for example, sample concentration or dilution ratio):




Click **OK**. The new information will display in the selected well:



- b. **To enter annotations for multiple wells or a row** - To select individual wells, press and hold the Control key and then select individual wells. To select a sequential set of wells or a full row, press and hold the Shift key, then select the first well and last well. Next, right click and select **Edit** or click **Edit** in the upper right corner of the pane. The following box will display:




Enter the reagent name, notes and specific attributes (for example, sample concentration or dilution ratio):



A dialog box titled "Well Content" with a close button (X) in the top right corner. It contains two input fields: "Name:" with the text "HeLa" and "Attribute:" with the text "1 mg/mL". At the bottom, there are two buttons: "OK" and "Cancel".

Click **OK**. The new information will display in the selected wells:




A screenshot of an assay template grid with 12 columns (1-12) and 6 rows (A-J). The grid is color-coded by row: Row A is orange, B is purple, C is light blue, D is teal, and J is yellow. Row A contains "Biot. Ladder" in column 1 and "HeLa" in column 7, with "1 mg/mL" in columns 2-6 and 8-12. Row B contains "Antibody Diluent" in columns 7-12. Row C contains "Blocking" in column 1 and "Primary Antibody" in columns 7-12. Row D contains "Streptavidin ..." in column 1 and "Secondary Antibody" in columns 7-12. Row J contains "Detection" in columns 7-12.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Biot. Ladder	1 mg/mL	HeLa	1 mg/mL								
B							Antibody Diluent					
C	Blocking						Primary Antibody					
D	Streptavidin ...						Secondary Antibody					
J							Detection					

*NOTE: When annotations are being edited, an asterisk will appear in the **Template** tab to indicate changes have been made and should be saved. Only Name and Attribute will be used by Compass to annotate the data.*

Step 7 - Save the Assay

1. Select **File** from the main menu and click **Save As**. Enter the assay name and click **Save**.




NOTE: New assays are saved in the Compass Assays directory.

Step 8 - Modify Default Analysis Parameters (Optional)

You can preset the parameters used to analyze run data generated with the assay.


1. Select **Edit** from the main menu and click **Default Analysis**. The following screen will display:



2. We recommend using the default parameters for Simple Western assays. However, you can modify any parameters as needed, then click **OK**. For detailed information on analysis parameters, please refer to "Compass Analysis Settings Overview" on page 232.

Making Changes to an Existing Assay

1. Select **File** in the main menu and click **Open Assay**.




2. A list of the last assays opened will display. Select one of these assays or click **Browse** to open the assay folder and select a different assay.
3. Follow the steps in "Creating a New Assay" on page 79 to make changes and save the assay.

Switching Between Open Assays

If more than one run file is open, you can switch between viewing the assays for each run. To do this:

1. Click the down arrow in the run box.




2. Select the run for the assay you want to view from the drop down list.


Creating a Template Assay

You can create and save template assays that can be used as a starting point for creating new assays. To do this:

1. Select **File** in the main menu. Click **Open Assay** to open an existing assay or click **New Assay** to open an existing template assay.
2. Follow the steps in “Creating a New Assay” on page 79 to make changes to the assay.
3. When changes are complete, select **File** in the main menu and click **Save As**. Select the New Assays folder:




4. Type the name for the new template assay and click **Save**.
5. Select **File** in the main menu and click **New Assay**. The new template assay will now be available in the drop down list:




Viewing and Changing the Detection Exposures

To view the current detection profile, roll the cursor over the exposures cell in Protocol pane:




While we recommend using the default assay detection profile, you can change the profile, exposure times and exposure sequences if needed. To do this:

1. Select the exposures cell in the Protocol pane and click the button, or double-click in the cell. The following screen will display:



Each row represents an individual exposure that will be taken during the run.

- a. **To change an existing exposure time** - Click in the exposure cell and enter a new time in seconds:



- b. **To delete an existing exposure** - Select a type or exposure cell and click **Remove**.
 - c. **To add a new exposure** - Select **Add**. A new exposure will be added to the end of the list. Click in the exposure cell and enter an exposure time in seconds.
2. Click **OK** to save and exit.

Copying Protocols and Templates

The steps and parameters in the Protocol pane can be copied and d into other documents, as can the graphic image of the annotations in the Template pane.

Copying an Assay Protocol

1. Click on the Protocol tab.
2. Select **Edit** in the main menu and click **Copy**.
3. Open a document (Microsoft® Word®, Excel®, etc.). Right click in the document and select. All assay protocol steps and parameters will be copied into the document as a list in the same format shown in the Protocol pane.

Copying an Assay Template


1. Click on the Template tab.
2. Select **Edit** in the main menu and click **Copy**.
3. Open a document (Microsoft® Word®, Excel®, etc.). Right click in the document and select. The assay template will be copied into the document as an image, exactly as it is shown in the Template pane.

Printing Protocols and Templates

The information in the Protocol pane can be printed, as can a graphic image of the information in the Template pane.

Printing an Assay Protocol


1. Click on the Protocol tab.
2. Select **File** in the main menu, click **Print**, and then click **Print Protocol**.



All assay protocol steps and parameters will be printed as a list in the same format shown in the Protocol pane.

Printing an Assay Template

1. Click on the Template tab.
2. Select **File** in the main menu, click **Print**, and then click **Print Template**.



The assay template will print as an image, exactly as it is shown in the Template pane.

Importing and Exporting Protocols and Templates

The assay protocol in an open assay or run file can be exported as a separate file, as can the assay template annotation information. This allows the same assay protocol and template information to be imported into another assay at a later time, rather than having to re-enter assay protocol and template information manually.

Importing an Assay Protocol


NOTE: Importing an assay protocol imports information into the Protocol pane only.

1. Open the assay you want to import the assay protocol in to.
2. Select **File** in the main menu and click **Import Protocol**.
3. Select a protocol file (*.protocol) and click **OK**. The imported information will display in the Protocol pane.

Exporting an Assay Protocol

NOTE: Exporting an assay protocol exports information in the Protocol pane only.

1. Open the assay you want to export the assay protocol from.
2. Select **File** in the main menu and click **Export Protocol**. The following window displays:



3. The default directory is Compass/Assays. Change the directory if needed.
4. Enter a protocol name and click **Save**. The protocol will be saved as a *.protocol file.

Importing an Assay Template


NOTE: Importing an assay template imports information into the Template pane only.

1. Open the assay you want to import the assay template in to.
2. Select **File** in the main menu and click **Import Template**.
3. Select a template file (*.template) and click **OK**. The imported information will display in the Template pane.

Exporting an Assay Template

NOTE: Exporting an assay template exports information in the Template pane only.

1. Open the assay you want to export the assay template from.
2. Select **File** in the main menu and click **Export Template**. The following window displays:



3. The default directory will be Compass/Assays. Change the directory if needed.
4. Enter a template name and click **Save**. The template will be saved as a *.template file.

Chapter 5:

Running a Charge Assay on Peggy Sue

Chapter Overview

- Starting a Run
- Stopping a Run

Starting a Run

Step 1 - Get Ready

1. Open Compass software.
2. Prepare instrument: empty waste, refill water and add a new manifold sponge.
3. Create or open desired assay file.
4. Prepare assay plate following the procedure described in the product insert.

IMPORTANT

To prevent well evaporation and ensure best results, keep a lid on the assay plate until ready to use.

5. While plate is spinning, add Wash Buffer, Anolyte and Catholyte to resource tray cups. Place capillary box in the designated resource tray position.


IMPORTANT

Capillaries are light sensitive. Keep the cover on the box until you are ready to transfer the capillary box to the resource tray.

6. Place assay plate into the sample tray of the instrument and press **Start**.

Step 2 - Start the Run

1. New run of an existing assay:
 - a. Select **File** in the main menu and click **Open Assay**.




- b. A list of the last five assays opened will display. Select one of these assays or click **Browse** to open the Assay folder and select a different assay.
2. Alternatively, choose **New Assay** and select **Peggy Sue Charge** to get the default Peggy Sue assay conditions.
- c. The **Start** button will display. This indicates than an assay has been loaded.



- a. Go to the Assay screen and verify this is the assay you want to use. If not, select **File** in the main menu, click **Open Assay**, and select another assay.
3. Click **Start**. This will launch the **Start Run Wizard**.


*NOTE: If the manifold was not cleaned prior to starting the run, a message indicating this will display. If this occurs, click **Yes** to cancel the run and perform the manifold cleaning.*

4. **Check Water and Waste.** The fluid levels in the accessory module bottles will be checked by the software. If the levels in both bottles will allow Peggy Sue to complete the run, the wizard screen will display **Water Level OK** and **Waste Level OK** messages. Click **Next** to proceed.



*NOTE: If the waste level is too high or the water level is too low to complete the run, messages to indicate one or both will be presented in this screen. If this occurs, fill or empty each bottle as indicated using the procedures outlined in the Peggy Sue User Guide. When this is complete, the error status will be automatically updated and allow the **Start Run Wizard** to proceed.*

- 5. Replace Sponge.** A new sponge should be used each time a new experimental run is started. Discard the old sponge and put a new sponge in the manifold wash station.




- 6. Start.** The resource tray will automatically open. Fill the Wash Buffer, Anolyte and Catholyte cups and insert the capillary box in the primary capillary box position. Add a secondary capillary box if more capillaries are required. Click **Next** when finished. The resource tray will close.


NOTE: Make sure enough capillaries are loaded to complete the run. A capillary box should be inserted in the primary position first and if additional capillaries are needed, use the secondary position. When the primary box is depleted, Peggy Sue will automatically move to the secondary box. Discard leftover Running Buffer, Wash Buffer and Matrix Removal Buffer prior to refilling bottles for new run.




NOTE: You an also refer to the labels on the resource tray for proper insertion of reagents.



7. **Load Sample Plate.** The sample tray will automatically open. Place the lidded 384-well plate on the cooling block in the tray, ensuring that the A1 well position is aligned with the upper left corner of the cold block. When this is complete, click **Next**. The sample tray will close.



A1 Position




NOTES:

Peggy Sue requires that plate lids be used on sample plates. If a lid is not detected, a message will be displayed in the Start Run Wizard. Compass software will reopen the sample tray to allow you to insert a lid.

When inserting the sample plate, ensure that the plate is firmly seated and level on the cold block. Plates that are not level can interfere with the movement of the sample tray.

8. **Begin the Automated Run.** The data file name will automatically default to the assay name appended with the current date and time. To change the file name, begin typing in the text box. To change the directory where the data file will be stored, click **Browse**:



Click **Start** to begin the run. Instrument status will change to running, and the stop button and progress bar will display:



The run will continue until complete (~14-19 hours depending on the assay).

Step 3 - Post-Run Procedures

1. Empty the capillary discard tray.
2. Remove the assay plate.
3. Dispose of the assay plate and capillaries. Disposal will depend on the samples that have been assayed. If sample origins are unknown, we recommend that used capillaries and plates be disposed of in bio-hazard waste.

!WARNING! SHARPS HAZARD

The capillaries may present a potential sharps hazard. Dispose of used capillaries in biomedical waste sharps containers.




!WARNING! BIOHAZARD

Samples and waste bottle contents should be handled by procedures recommended in the CDC/NIH manual: Biosafety in Microbiological and Biomedical Laboratories (BMBL). The manual is available from the U.S. Government Printing Office or online at <http://www.cdc.gov/biosafety/publications/bmbl5/>.

Depending on the samples used, waste bottle contents may constitute a biohazard. Use precaution when emptying the waste bottle. Dispose of waste bottle contents in accordance with good laboratory practices and local, state/provincial, or national environmental and health regulations. Read and understand the Material Safety Data Sheets (MSDSs) provided by the manufacturers of the chemicals in the waste vial before you store, handle, or dispose of chemical waste.


Stopping a Run

1. To stop a run, click **Stop**. When the run stops, instrument status will go to Not Ready and a Cleanup button displays:



NOTE: If a run is stopped prior to completion, Peggy Sue must perform a cleaning protocol to remove used capillaries from system trays. This is required to prepare the system for the next run.

2. Click **Cleanup**.



Allow Peggy Sue to complete the cleaning protocol which takes about ten minutes. When complete, instrument status will change to Ready and a new run can be started.

Chapter 6:


Run Status

Chapter Overview

- Run Summary Screen Overview
- Opening Run Files
- Viewing File and Run Status Information
- Watching Standards Separation Movies
- Viewing Current and Voltage Plots
- Switching Between Open Run Files
- Closing Run Files

Run Summary Screen Overview

The Run Summary screen is used to monitor run progress, watch movies of the fluorescent standards separation, and view current and voltage plots for a run. To access this screen, click Run Summary in the screen tab:



Run Summary Screen Panes

The Run Summary screen has three panes:

- **Status** - Displays run file information and current status of a run in progress.
- **Separation** - Lets you view a movie of the fluorescent standards separation for each cycle of the experimental run.
- **IV Plot** - Lets you view plots of the total current and voltage measured during separation for all capillaries for each cycle of the experimental run.

Run: 2015-07-17_14-03-48_Peggy Sue Size Split Running Buffer 8-cycle

Status | History

Path: C:\Users\ppiatt\Desktop
 Assay: Peggy Sue Size Split Running Buffer
 Kit Info: Regular: 12-230 kDa, Split Running Buffer
 Schedule: Overlapping with hold
 Instrument: Peggy-Sue : Peggy Sue SW0105

Started: Fri 2:09 PM Jul 17, 2015 PDT
 Completed: Sat 5:57 AM Jul 18, 2015 PDT

Cycle	Sample	Sep	Hold	B	1*	Detect	Results
1	2:09 PM	2:17 PM	3:33 PM	1:22 AM	1:44 AM	2:23 AM	2:42 AM
2	3:33 PM	3:41 PM	4:57 PM	1:49 AM	2:11 AM	2:49 AM	3:08 AM
3	4:57 PM	5:05 PM	6:21 PM	2:15 AM	2:37 AM	3:16 AM	3:34 AM
4							

The right pane shows a video of 'Cycle 1' with a play button at the bottom.

Software Menus Active in the Run Summary Screen


The following software menus are available:

- File
- Edit
- Instrument (when Compass is connected to an instrument)
- Window
- Help

The File and Edit menu options specific to the Run Summary screen are described next.

File Menu


The following **File** menu options are active:



- **Open Run** - Opens a run file.
- **Add Run** - Open and view other run files in addition to the one that is already open.
- **Close** - Closes the run file currently being viewed.
- **Close All** - Closes all open run files.
- **Exit** - Closes Compass.

Edit Menu

The following Edit menu options are active:




- **Preferences** - Set and save custom preferences for data export, plot colors in the graph and Twitter settings. See Chapter 10, *Setting Your Preferences* for more information.

Opening Run Files

You can open one run file or multiple run files at a time to compare information between runs.

Opening One Run File


1. Select **File** in the main menu and click **Open Run**.




2. A list of the last 10 runs opened will display. Select one of these runs or click **Browse** to open the Runs folder and select a different file.

Opening Multiple Run Files

1. To open the first run file, select **File** in the main menu and click **Open Run**.




2. A list of the last 10 runs opened will display. Select one of these runs or click **Browse** to open the Runs folder and select a different file.
3. To open another run file, select **File** in the main menu and click **Add Run**.



4. A list of the last 10 runs opened will display. Select one of these runs or click **Browse** to open the Runs folder and select a different file.
5. Repeat the last two steps to open additional runs.

Viewing File and Run Status Information

Information specific to each run file is shown in the Status pane:



The run file name, path (directory location) and assay used is displayed along with instrument serial number and the run start/complete date and time.

- **To go to the run file directory location** - Double click the path hyperlink, or right-click and select **Open Directory**.
- **To copy the path** - Right-click on the path hyperlink and click **Copy**. The path can then be copied into documents. The path can also be copied into the Windows Explorer address bar to launch Compass and open the run file automatically.
- **Kit info** - Compass v2.7 and higher displays the type of kit used to run the assay (regular for immunoassays, total protein or charge), the molecular range and whether or not the split Running Buffer was used.
- **Plate S/N** - For Wes, the plate serial number (S/N) information is captured and displayed for all assays.

Status	History
Run	2015-07-31_3A7_3B1_C6138_029_K562_6
Path	C:\Users\fdeng\Documents\Compass\Run
Kit Info	High Molecular Weight: 66-440 kDa
Instrument	Wes : Wes WS2284 - WS2284
Plate S/N	7670330039 - Split Running Buffer
Started	Fri 1:21 PM Jul 31, 2015 PDT
Completed	Fri 4:06 PM Jul 31, 2015 PDT

Assay Steps: Size-based Assays

Each assay step that is executed during the run is represented by an icon. The time the individual step is scheduled to start is shown under each icon. Details on what happens during each assay step is as follows:

Step	Description
 Sample 12:54 PM	Sample Loading Step - Capillaries are moved to the assay plate in the sample tray where Separation Matrix, Stacking Matrix, biotinylated ladder and samples are aspirated. Capillaries are then transferred to the separation tray.

Step	Description
Sep  12:56 PM	<p>Separation Step - Samples and fluorescent standards are separated in the capillaries. Capillaries are then exposed to UV light, which immobilizes the separated proteins to the walls of each capillary. During the separation, voltage and current are monitored and displayed in the IV Plot. Separation of the fluorescent standards is recorded and a movie of the separation is compiled for viewing when separation is complete in the Separation pane. After separation, Matrix Removal Buffer is aspirated.</p>
Block  2:00 PM	<p>Blocking Step - Capillaries are moved to the assay plate in the sample tray and blocking reagent (Antibody Diluent) is aspirated. Capillaries are then transferred to an incubator tray. When incubation is complete, Wash Buffer is aspirated.</p>
1°  2:16 PM	<p>Primary Antibody or Total Protein Labeling Reagent (1°) Step - Capillaries are moved to the assay plate in the sample tray and primary antibody or labeling reagent is aspirated. Capillaries are then transferred to an incubator tray. When incubation is complete, Wash Buffer is aspirated.</p>
2°  4:22 PM	<p>Secondary Antibody or Total Protein Streptavidin-HRP (2°) Step - Capillaries are moved to the assay plate in the sample tray and secondary HRP conjugate and Streptavidin-HRP (Immunoassays), or Total Protein Streptavidin-HRP (Total Protein Assays) is aspirated. Capillaries are then transferred to an incubator tray. When incubation is complete, Wash Buffer is aspirated.</p>
Detect  5:28 PM	<p>Detect Step - Capillaries are moved to the assay plate in the sample tray and Luminol-Peroxide solution is aspirated. Capillaries are then transferred to the separation tray where the emitted chemiluminescent light is detected with the CCD camera.</p>
Results  6:00 PM	<p>Results Step - Results are available in the Analysis screen.</p>
Hold  1:10 PM	<p>Hold Step - This step is used for overlapping or overlapping with hold schedule options. A cycle that has completed a specific step will go into a hold step while the remaining cycles execute the same step in parallel. For example, once cycle 1 completes the separate step, it may go into a hold step while cycle 2 executes the separate step. This execute and hold pattern will continue until all cycles have completed the separate step.</p>

When a run is in progress, icons for steps that have not executed will be grey (inactive). In the following example, the Secondary Antibody (2°) step is executing and the Detect and Results steps have not started:



Assay Steps: Charge-based Assays

Each assay step that is executed during the run is represented by an icon. The time the individual step is scheduled to start is shown under each icon. Details on what happens during each assay step is as follows:

Step	Description
 12:54 PM	Sample Loading Step - Capillaries are moved to the assay plate in the sample tray and samples are aspirated. Capillaries are then transferred to the separation tray.
 12:56 PM	Separation Step - Samples, ampholyte mix and fluorescent pl standards are separated in the capillaries. Capillaries are then exposed to UV light, which immobilizes the separated proteins to the walls of each capillary. During the separation, voltage and current are monitored and displayed in the IV Plot. Separation of the fluorescent standards is recorded and a movie of the separation is compiled for viewing in the Separation pane after separation is complete.
 2:16 PM	Primary Antibody (1°) Step - Capillaries are moved to the assay plate in the sample tray and primary antibody is aspirated. Capillaries are then transferred to an incubator tray. When incubation is complete, Wash Buffer is aspirated.
 4:22 PM	Secondary Antibody (2°) Step - Capillaries are moved to the assay plate in the sample tray and secondary HRP-conjugated antibody is aspirated. Capillaries are then transferred to an incubator tray. When incubation is complete, Wash Buffer is aspirated.
 5:28 PM	Detect Step - Capillaries are moved to the assay plate in the sample tray and Luminol-Peroxide solution is aspirated. Capillaries are then transferred to the separation tray where the emitted chemiluminescent light is detected with the CCD camera.
 6:00 PM	Results Step - Results are available in the Analysis screen.

Step	Description
Hold  1:10 PM	Hold Step - This step is used for overlapping or overlapping with hold schedule options. A cycle that has completed a specific step will go into a hold step while the remaining cycles execute the same step in parallel. For example, once cycle1 completes the separate step, it may go into a hold step while cycle 2 executes the separate step. This execute and hold pattern will continue until all cycles have completed the separate step.


When a run is in progress, icons for steps that have not executed will be grey (inactive). In the following example, the Secondary Antibody (2°) step is executing and the Detect and Results steps have not started:



Watching Standards Separation Movies

You can view a movie of the fluorescent standards separation in all 12 capillaries. To do this:


1. Click the **Separation** tab.



2. The player control panel has play/pause, rewind and fast forward buttons, and a slider bar that allows you to scroll through the movie manually:




Click **Play** (button on far left) to view the movie. In the examples below, standards for a size assay are on the left, and standards for a charge assay are on the right:



NOTE: Complete separation movies of the fluorescent standards are not available until the separation step has finished executing. If the movie is played while the separation step is executing, the movie will only show separation progress up to the current point in time.

Viewing Current and Voltage Plots

You can view plots of the total current and voltage measured during separation for all 12 capillaries. To do this, click the **IV Plot** tab. In the examples below, the IV plot for a size assay is on the left, and the IV plot for a charge assay is on the right:



The blue Y-axis and plot shows the run voltage in volts (V), and the red Y-axis and plot shows the run current in micro amps (μA). The X-axis displays time in seconds.


- **To zoom in on an area of the plot** - Hold the mouse button down and draw a box around the area with the mouse.
- **To zoom out** - Click **Zoom Out** in the upper right corner of the pane.

NOTE: The IV plot for a run in progress will not be available until the separation step starts executing. The plot is then displayed in real-time.

Switching Between Open Run Files

If more than one run file is open, you can switch between viewing the run information in each. To do this:

1. Click the down arrow in the run box.



2. Select the run you want to view from the drop down list.

Closing Run Files

If more than one run file is open, you can close just one file or all the open files at the same time.

- **To close the run file being viewed** - Select **File** from the main menu and click **Close**.
- **To close all open run files** - Select **File** from the main menu and click **Close All**.

Chapter 7:


Controlling Wes, Sally Sue and Peggy Sue

Chapter Overview

- Instrument Control
- Self Test
- Viewing and Changing System Properties
- Viewing Log Files
- Status Modes

Instrument Control

The Instrument menu allows you to control Wes, Sally Sue and Peggy Sue.




NOTE: Instrument menu options are active only when a computer with Compass software is connected directly to Wes, Sally Sue or Peggy Sue.

Starting a New Run

To start a new run, select **Instrument** in the main menu and click **New Run**. Then follow the steps described in “Step 2 - Start the Run” on page 64 for size assays or “Step 2 - Start the Run” on page 106 for charge assays.


Opening Trays (Sally Sue and Peggy Sue)

To open any of the five trays, select **Instrument** and click **Open Trays**. The tray control window will appear:



Open a tray by clicking on its button. The button will become highlighted indicating the tray is open.

NOTE: Only one tray can be open at a time.



To close a tray, click the corresponding tray button again.

NOTE: If the tray control window is closed when a tray is open, the tray will close automatically.

Cleaning

Two cleaning options are available for Sally Sue and Peggy Sue.

Manual Clean


This option is used for general manual cleaning and cleaning the manifold head. To do a manual cleaning, select **Instrument** and click **Manual Clean**. The manifold head will move to a safe position for easy access and the vacuum will turn on.

NOTE: Please contact Protein Simple Technical Support if you have any questions regarding the manifold cleaning procedure.


Cleanup

This option is a fully automated cleaning step. The manifold head is flushed, the separation tray troughs are aspirated and washed, and any capillaries left in the trays or gripper are picked up and discarded. This option should be selected when the instrument has not been used for more than a week or if a run error occurs. Cleaning takes about eight minutes to complete.

To start the protocol, select **Instrument** and click **Cleanup**. A window will appear with instructions:



Sally Sue's or Peggy Sue's status will change to cleaning, and the stop button and the cleaning progress bar display. The **Assay** screen provides cleaning status details:




When cleaning is complete, instrument status will change to **Ready**.



Cleaning After a Run Error

Additional cleaning steps are required if an error occurs that stops the run. When this happens, the red Error status light on Sally Sue's or Peggy Sue's front panel will come on.

Click on the **Reset** button displayed in Compass software. The following instructions will appear:




If the troughs in the separation tray are empty, click on **OK** and proceed with "Cleanup" on page 131.

If Running Buffer is present in the separation tray, click on **Cancel** and manually remove the buffer. Evaporation of the Running Buffer will result in a highly viscous residue which the automatic cleaning feature cannot remove.

To remove the Running Buffer:


1. Select **Instrument** and click **Open Trays**.



2. Click **Separation** to open the separation tray.
3. Add 800 μL of deionized water to the troughs in the separation tray and soak for 20 minutes.
4. Remove the water by either aspirating with a pipette or with the vacuum wand located on the inside of Sally Sue's or Peggy Sue's left door.
5. Repeat the steps above until the Running Buffer or residues are completely removed.
6. To complete the cleaning process, select **Instrument** and click **Cleanup**.

Self Test

Wes, Sally Sue and Peggy Sue can perform a series of self tests to check for proper instrument performance. To start the test, select **Instrument** and click **Self Test**. The test takes approximately two minutes.




NOTE: We recommend performing the self test prior to starting a run.

To view the test log at completion of the test, select **Instrument**, click **Properties** and click **View Test Log**. See "Self Test Logs" on page 138 for more information.

Viewing and Changing System Properties

Select **Instrument** and click **Properties** to display system properties which include:

- Name
- Location
- Type
- Serial number
- Instrument software version (firmware)
- Network name and address
- Date and time of the instrument clock




- **To change system name or location** - click in the name or location boxes and enter the new information.
- **To sync the instrument clock with the computer** - click **Set to PC time**.

Viewing Log Files


Error Logs

1. Select **Instrument** in the main menu and click **Properties** to display system properties.
2. Click **View Error Log**. A list of system logs will display:



Name	Date	Size
embedded.log	2012-10-30 13:59	4729499
temperature.log	2012-10-30 13:58	24198
temperature.log.2012-10-29	2012-10-30 09:39	134780
lamp.log	2012-10-30 03:28	153166
error.log	2012-10-29 17:06	873831
temperature.log.2012-10-28	2012-10-29 09:39	134811
temperature.log.2012-10-27	2012-10-28 09:38	134807
temperature.log.2012-10-26	2012-10-27 09:37	134725
temperature.log.2012-10-25	2012-10-26 09:37	532087
stdout.log	2012-10-24 06:54	244772
temperature.log.2012-10-21	2012-10-22 10:44	134793
temperature.log.2012-10-20	2012-10-21 10:43	134782
embedded.log.1	2012-10-02 21:50	8388603
embedded.log.2	2012-07-23 19:26	8388540
embedded.log.3	2012-05-24 23:12	8388565
embedded.log.4	2012-03-31 07:00	8388547


3. Select a log file and click **View**. The log details will display:



4. Click **Save File As** to save a copy of the log file.

Self Test Logs


1. Select **Instrument** in the main menu and click **Properties** to display system properties.
2. Click **View Test Log**. A list of self test logs will display:



The screenshot shows a window titled "PL0003 Self Tests" with a close button in the top right corner. The window contains a table with two columns: "Test Date" and "Status". The table lists ten test entries with their respective dates and outcomes. At the bottom of the window, there are two buttons: "View" and "Done".

Test Date	Status
2012-10-23_14.33.20	PASSED
2012-10-11_13.11.50	PASSED
2012-10-10_16.22.29	FAILED
2012-10-10_16.15.12	FAILED
2012-10-10_14.25.48	STOPPED
2012-10-10_14.24.41	FAILED
2012-10-05_14.05.02	PASSED
2012-10-05_14.02.32	PASSED
2012-10-05_13.59.20	PASSED
2012-06-19_16.05.22	PASSED

3. Select a log file and click **View**. The individual test details will display:



Name	Start	Duration (sec)	Result	Failure Reason
Dark Masters	12:00:21	0.024	PASSED	
Camera	12:00:10	11.720	PASSED	
Fan Temp Sensors	12:00:08	0.056	PASSED	
Level	12:00:08	0.161	PASSED	
Filter Not In	12:00:07	0.143	PASSED	
Filter Not Out	12:00:06	0.163	PASSED	
Filter Out	12:00:04	0.593	PASSED	
Filter In	12:00:03	0.665	PASSED	
LED Low Not On	12:00:02	0.025	PASSED	
LED Low On	12:00:00	0.017	PASSED	
LED Not On	11:59:59	0.016	PASSED	
LED On	11:59:58	0.016	PASSED	
Voltage On	11:59:48	8.965	PASSED	
Lid Remover Not Down	11:59:47	0.174	PASSED	
Lid Remover Not Up	11:59:46	0.147	PASSED	
Lid Remover Up	11:59:44	0.468	PASSED	
Lid Remover Down	11:59:43	0.309	PASSED	
Pipet Not Down	11:59:41	0.052	PASSED	
Pipet Not Up	11:59:40	0.053	PASSED	
Pipet Up	11:59:39	0.248	PASSED	
Pipet Down	11:59:38	0.321	PASSED	

4. Click **Save File As** to save a copy of the log file.

Status Modes

The instrument status bar displays status, buttons and progress bars depending on what Wes, Sally Sue or Peggy Sue is doing.

- **Ready/Start button** - The instrument is ready and an assay is loaded. Click **Start** to begin a run.
- **Not Ready/Clean button** - The instrument is not ready and must perform system cleaning. Click **Clean** to start the cleaning protocol.
- **Not Ready/Reset button** - The instrument is not ready and must reinitialize. Click **Reset** to start the initialization protocol.
- **Running/Stop button** - The instrument is running an assay. The run name, time the run started and when it will complete display in the run progress bar. Click **Stop** to stop the run.
- **Cleaning/button not active** - The instrument is running a cleaning protocol. The time the cleaning protocol started and when it will complete display in the run progress bar.
- **Error/Reset button** - An error has occurred. Go to the **Status** window in the **Run Summary** screen to view details. When the source of the error is corrected, click **Reset**.

Chapter 8:

Size Assay Data Analysis


Chapter Overview

- Analysis Screen Overview
- Opening Run Files
- How Run Data is Displayed in the Analysis Screen
- Viewing Run Data
- Compass Run Data Notifications and Warnings
- Checking Your Results
- Group Statistics
- Copying Data Views and Results Tables
- Exporting Run Files
- Changing Sample Protein Identification
- Changing the Virtual Blot View
- Changing the Electropherogram View
- Closing Run Files
- Compass Analysis Settings Overview
- Advanced Analysis Settings
- Images Analysis Settings
- Peak Fit Analysis Settings
- Peak Names Settings
- Standard Curve Settings
- System or Loading Control Settings

- Standard Curve Settings
- Importing and Exporting Analysis Settings

Analysis Screen Overview


The Analysis screen is used to view run data including electropherograms, capillary images, lane view data and tabulated results. Any post-run analysis is also done in this screen. To access this screen, click **Analysis** in the screen tab:



Analysis Screen Panes

The Analysis screen has six panes, each displays the following data for up to 96 capillaries per experimental run:

- **Experiment** - Lists the assay protocol steps and assay template information.
- **Graph** - Displays electropherogram data for sample proteins, fluorescent standards or capillary registrations.
- **Image** - Displays a 12-capillary image of the separated sample proteins, fluorescent standards or capillary registrations.
- **Lane** - Displays data for sample proteins, fluorescent standards or capillary registrations as bands in individual lanes. This virtual blot-like image is similar to traditional blot results.
- **Peaks** - Lists the tabulated results for sample proteins and information for identified fluorescent standards and capillary registrations.
- **Capillaries** - Displays a list of the sample proteins Compass named automatically using the user-defined peak name analysis parameters.



NOTE: The reported molecular weight for sample proteins in Compass may vary slightly from predicted molecular weights based on sample and assay conditions.

Software Menus Active in the Analysis Screen


The following software menus are available:

- File
- Edit
- View
- Instrument (when Compass is connected to Wes, Sally Sue or Peggy Sue)
- Window
- Help

The File, Edit and View menu options specific to the Analysis screen are described next.

File Menu


The following File menu options are active:



- **Open Run** - Opens a run file.
- **Add Run** - Opens and views other run files in addition to those that are already open.
- **Close** - Closes the run file currently being viewed.
- **Close All** - Closes all open run files.
- **Save** - Saves changes to the open run file.
- **Save As** - Saves changes to the open run file under a different file name.
- **Export Tables** - Exports the results for all capillaries in the run in .txt format.
- **Export Spectra** - Exports the raw data traces for each capillary in the run in .txt or .cdf format.
- **Exit** - Closes Compass.

Edit Menu


The following Edit menu options are active:



- **Copy** - Lets you copy data shown in the graph, lane, peaks or capillaries panes. See “Copying Data Views and Results Tables” on page 194 for more information.
- **Analysis** - Displays the analysis settings used to analyze the run data and lets you change them as needed. See “Compass Analysis Settings Overview” on page 232 for more information.
- **Preferences** - Lets you set and save custom preferences for data export, plot colors in the graph and Twitter settings. See Chapter 10, “Setting Your Preferences” for more information.

View Menu

The following View menu options are active:



- **Single View** - Displays data in a per capillary (single) view format.
- **Multiple View** - Displays data in a per 12- or 25-capillary (multiple) view format.
- **Standards** - Lets you change the data view to show only the fluorescent standards.
- **Registration** - Lets you change the data view to show only the capillary registrations (Sally Sue and Peggy Sue only).
- **Samples** - Lets you change the data view to show sample proteins.
- **Filter** - Lets you display data only for specific capillaries or named proteins.
- **View Region** - Lets you change the molecular weight (x-axis) range of the data displayed.

- **Show Hidden**- Shows capillaries that are hidden from the data view.


Opening Run Files

You can open one or multiple run files at a time to compare data between runs.

Opening One Run File

NOTE: If you need to open a run file when another run is executing, launch another instance of Compass first, then open the file.

1. Select **File** in the main menu and click **Open Run**.




2. A list of the last 10 runs opened will display. Select one of these runs or click **Browse** to open the Runs folder and select a different file.


Opening Multiple Run Files

NOTE: If you need to open a run file when another run is executing, launch another instance of Compass first, then open the file.


1. To open the first run file, select **File** in the main menu and click **Open Run**.



2. A list of the last 10 runs opened will display. Select one of these runs or click **Browse** to open the Runs folder and select a different file.
3. To open another run file, select **File** in the main menu and click **Add Run**.



4. A list of the last 10 runs opened will display. Select one of these runs or click **Browse** to open the Runs folder and select a different file. When a run is added, its data will append to the open run file, and will display as a second set of up to 25 capillaries (Wes) or 96 capillaries (Sally Sue/Peggy Sue) in all screen panes. The second run file name will also appear in the Compass title bar:




5. Repeat the last two steps to add additional runs.

How Run Data is Displayed in the Analysis Screen

Data in the run file is compartmentalized for easy review in one of six panes in the Analysis screen.

Experiment Pane: Assay and Capillary Information

The experiment pane displays assay and capillary information for each of the 25 capillaries (Wes) or 96 capillaries (Sally Sue/Peggy Sue) in the run. A maximized view of the experiment pane is shown below.



Sample	Primary	Cycle	Cap	S	1	2	3
✓ Biotinylated Ladder	Blocking	1	1	C1	D1	E1	F1
✓ K562	anti-ERK1/2	1	2	C2	D2	E2	F2
K562	anti-ERK1/2	1	3	C3	D3	E3	F3
✓ K562	anti-ERK1/2	1	4	C4	D4	E4	F4
✓ K562	anti-ERK1/2	1	5	C5	D5	E5	F5
✓ K562	anti-ERK1/2	1	6	C6	D6	E6	F6
K562	anti-ERK1/2	1	7	C7	D7	E7	F7
K562	anti-ERK1/2	1	8	C8	D8	E8	F8
K562	anti-ERK1/2	1	9	C9	D9	E9	F9
K562	anti-ERK1/2	1	10	C10	D10	E10	F10
K562	anti-ERK1/2	1	11	C11	D11	E11	F11
RTU K562	anti-ERK1/2	1	12	C12	D12	E12	F12

- **To view all columns** - Click the **Experiment** tab, then use the scroll bar or click **Maximize** in the upper right corner.
- **To resize columns** - Click the **Experiment** tab. Roll the mouse over a column border until the sizing arrow appears, then click and drag to resize.

Experiment pane column descriptions for the default assay are as follows:

- **Sample** - Sample name. If sample names were entered in the assay template (Assay screen), those names will display here. Otherwise, Sample (default name) will display.

NOTE: Data notification icons will display in the sample column if Compass detected a potential analysis issue or data was manually modified by the user. For more information see "Compass Run Data Notifications and Warnings" on page 176.


- **Primary** - Primary antibody name. If primary antibody names were entered in the assay template (Assay screen), those names will display here. Otherwise, Primary (default name) will display.
- **Cycle** - Run cycle number. There are 25 capillaries in one cycle for Wes and 12 capillaries in one cycle for Sally Sue and Peggy Sue.

NOTE: Sally Sue and Peggy Sue can run up to eight cycles (12 capillaries per cycle) in an experiment. Each cycle is designated in the Experiment summary.

- **Cap** - Capillary number.
- **S** - Well on the assay plate used for sample.
- **1** - Well on the assay plate used for primary antibody or Total Protein labeling reagent.
- **2** - Well on the assay plate used for secondary HRP-conjugate or Total Protein Streptavidin-HRP.

Graph Pane: Electropherogram Data


Click the **Graph** tab to view data for sample proteins, fluorescent standards or capillary registrations. Data for samples is shown in the following example, and proteins are displayed as peaks:



More Graph view options will be described in more detail in "Changing the Electropherogram View" on page 209.



Image Pane: Capillary Separation Image Data

Click the **Image** tab to view final separation images of sample proteins, fluorescent standards or capillary registrations for up to 96 capillaries per experimental run. Image data for samples is shown in the following example:



Lane Pane: Virtual Blot-Like Image Data

Click the **Lane** tab to view data for sample proteins, fluorescent standards or capillary registrations as bands in individual lanes. This virtual blot-like image is similar to traditional blot results. Data for samples in the lane view for an Immunoassay (immunodetected proteins) and from a Total Protein Assay are displayed as bands in the examples below.




To view information for a band, roll the mouse over a band until the info box appears.

NOTE: The reported molecular weight for sample proteins detected by an Immunoassay or a Total Protein Assay in Compass may vary slightly from predicted molecular weights based on sample and assay conditions.

Lane data displayed in the virtual blot is automatically aligned by Compass. To view raw, unaligned lane data and learn more about virtual blot viewing options, see “Changing the Virtual Blot View” on page 203.

Std Curve Pane: Standard Curve Fit Data

Click the **Std Curve** tab to view a linear or 4-parameter curve fit of your standard curve proteins. Data for the samples in the std curve view is shown in the following example where proteins are displayed as dots.



Concentration of the sample proteins is automatically determined by Compass based on defined values of the standard curve.

Peaks Pane: Calculated Results

Click the **Peaks** tab to view tabulated results for sample proteins, fluorescent standards or capillary registrations. Each row in the table shows the individual results for each peak detected in each capillary. Data for samples in the peaks table for both an Immunoassay and a Total Protein Assay are shown in the following examples:

Sample	Primary	Cap	Peak	Name	Position	MW (kDa)	Height	Area	% Area	Width	S/N
Biotinylat...	Blocking	C1:1	5	Ldr 116	761	116	3175.2	59215		17.5	207.0
Biotinylat...	Blocking	C1:1	6	Ldr 180	865	180	3775.6	78654		19.6	229.9
K562	anti-E...	C1:2	1	ERK2	522	43	1380.3	21641	79.6	14.7	289.7
K562	anti-E...	C1:2	2	ERK1	549	47	353.2	5552	20.4	14.8	35.9
K562	anti-E...	C1:3	1	ERK2	528	43	1400.1	21546	78.8	14.5	205.6
K562	anti-E...	C1:3	2	ERK1	555	47	342.7	5792	21.2	15.9	25.7
K562	anti-E...	C1:4	1	ERK2	518	43	1373.1	21542	79.1	14.7	304.9
K562	anti-F...	C1:4	2	FRK1	545	47	359.6	5700	20.9	14.9	42.0

Sample	Primary	Cap	Peak	Name	Position	MW (kDa)	Height	Area	% Area	% Total	Width	S/N	Baseline
Biot. Lad...	Antibo...	C1:1	5	Ldr 180	632	180	3547.2	66440		16.56	29.0	148.9	1034.6
Biot. Lad...	Antibo...	C1:1	6	Ldr 230	657	230	3272.5	60980		15.20	35.0	138.2	1034.5
HeLa	Antibo...	C1:2	8	DNAK	536	73	1621.6	41485	100.0	14.11	31.0	18.8	1264.3
HeLa + 3...	Antibo...	C1:3	7	DNAK	536	74	7017.3	131139	100.0	33.17	41.0	370.6	1408.6
HeLa + 1...	Antibo...	C1:4	10	DNAK	536	75	5663.3	122330	100.0	23.36	41.0	219.0	1088.8
HeLa + 7...	Antibo...	C1:5	8	DNAK	536	74	3467.4	77588	100.0	21.14	37.0	115.8	1310.0
HeLa + 3...	Antibo...	C1:6	10	DNAK	536	74	2452.4	58625	100.0	17.03	35.0	71.7	1313.9
HeLa	Antibo...	C1:7	7	DNAK	533	72	1379.5	36572	100.0	14.21	32.0	17.3	1356.2
HeLa + 3...	Antibo...	C1:8	8	DNAK	537	75	6249.9	113059	100.0	29.93	39.0	306.9	1322.4
HeLa + 1...	Antibo...	C1:9	9	DNAK	537	75	5472.6	104042	100.0	21.65	36.0	162.8	1134.9

NOTES:

Peaks that Compass names automatically using the user-defined peak name analysis parameters are color-coded.

The reported molecular weight for sample proteins detected using an Immunoassay or a Total Protein Assay in Compass may vary slightly from predicted molecular weights based on sample and assay conditions.

- **To view all rows** - Click the **Peaks** tab, then use the scroll bar or click **Maximize** in the upper right corner.
- **To resize columns** - Click the **Peaks** tab. Roll the mouse over a column border until the sizing arrow appears, then click and drag to resize.

Peak table column descriptions are as follows:

- **Sample** - Sample name. If sample names were entered in the assay template (Assay screen), those names will display here. Otherwise, Sample (default name) will display.
- **Primary** - Primary antibody name. For Immunoassays, if primary antibody names were entered in the assay template (Assay screen), those names will display here. Otherwise, Primary (default name) will display. For Total Protein Assays, Antibody Diluent will display as the default name.
- **Cap** - Cycle and capillary number. For example, C1:3 indicates cycle 1, capillary 3.

NOTE: Sally Sue and Peggy Sue run up to eight cycles in an experiment, 12 capillaries at a time. Capillary and cycle numbers are displayed in the Experiment tab.

- **Peak** - Peak number. Peaks are numbered in order of detection.
- **Name** - Peak name. Displays peaks that Compass named automatically using the user-defined peak name analysis parameters. Cells will be blank if Compass was not able to name the peak or if naming parameters were not entered.
- **Position** - Displays the pixel position of a peak in the image.
- **MW (kDa)** - Displays the calculated molecular weight in kDa for the peak (shown for sample data only).
- **Height** - Displays the calculated peak height.
- **Area** - Displays the calculated peak area (shown for sample data only).
- **% Area** - Reported when area is calculated using the Gaussian method (default for Immunoassays, see "Peak Find Settings" on page 245 for more information). Displays the calculated percent area for the named peak compared to all named peaks. This value results from dividing the individual peak area by the sum of all named peak areas for the capillary and multiplying by 100 (shown for named peak sample data only).
- **% Total** - Reported when area is calculated using the Dropped Line method (default for Total Protein Assays, see "Peak Find Settings" on page 245 for more information). Displays the calculated percent area for the named peak compared to the total area measured in the capillary. This value results from dividing the individual peak area by the entire area under the curve for the capillary and multiplying by 100.
- **Conc (Concentration)** - Displays the calculated concentration of protein for the named peak. This column will not appear unless a standard curve is defined.
- **Corr. Area** - Displays the calculated area for the named peak compared to the area for the loading control. This column will not appear unless a loading control is defined.
- **Width** - Displays the calculated peak width (shown for sample data only).
- **S/N** - Displays the calculated signal to noise ratio for the peak (shown for sample data only).

Capillaries Pane: User-Specified Peak Names

Click the **Capillaries** tab to view tabulated results for sample proteins associated with specific primary antibodies or total protein detection in the run data. Compass labels these sample peaks automatically using user-defined peak name settings. Each row in the table shows the individual results for the named peaks detected in each capillary. Data for samples in the capillaries table is shown in the following example for an Immunoassay with Gaussian fit.

Sample	Primary	Capillary	System Co...	ERK1
HeLa	ERK1+...	C1:2	44704	55899
HeLa	ERK1+...	C1:3	42886	63846
HeLa	ERK1+...	C1:4	50291	65303
HeLa	ERK1+...	C1:5	45954	57780
HeLa	ERK1+...	C1:6	45887	54196
HeLa	ERK1+...	C1:7	48500	68243
HeLa	ERK1+...	C1:8	47257	56120
HeLa	ERK1+...	C1:9	49193	59797
HeLa	ERK1+...	C1:10	47426	57923
HeLa	ERK1+...	C1:11	42898	45543
HeLa	ERK1+...	C1:12	43542	49701

For Total Protein Assays where area calculation is handled via Dropped Lines by default, the Capillaries tab is shown below. Compass labels the sample peaks automatically using user-defined peak name settings. Each row in the table shows the individual results for the named peaks detected in each capillary.

Sample	Primary	Capillary	Total Area	DNAK
HeLa	Antibo...	C1:2	294092	14.11
HeLa + 3...	Antibo...	C1:3	395361	33.17
HeLa + 1...	Antibo...	C1:4	523599	23.36
HeLa + 7...	Antibo...	C1:5	367106	21.14
HeLa + 3...	Antibo...	C1:6	344206	17.03
HeLa	Antibo...	C1:7	257414	14.21
HeLa + 3...	Antibo...	C1:8	377685	29.93
HeLa + 1...	Antibo...	C1:9	480506	21.65
HeLa + 7...	Antibo...	C1:10	483433	20.36

NOTES:

Peaks that Compass names automatically with user-defined peak name settings are color-coded.

Information displayed for fluorescent standards and capillary registration data will be for identified standards or registration peaks.


- **To view all rows** - Click the **Capillaries** tab, then use the scroll bar or click **Maximize** in the upper right corner.
- **To resize columns** - Click the **Capillaries** tab. Roll the mouse over a column border until the sizing arrow appears, then click and drag to resize.

Capillaries table column descriptions are as follows:

- **Sample** - Sample name. If sample names were entered in the assay template (Assay screen), those names will display here. Otherwise, Sample (default name) will display.
 - **Primary** - Primary antibody name. If primary antibody names were entered in the assay template (Assay screen), those names will display here. Otherwise, Primary (default name) will display.
 - **Capillary** - Cycle and capillary number. For example, C1:3 indicates cycle 1, capillary 3.
-


NOTE: Sally Sue and Peggy Sue run 12 capillaries at a time in a cycle and are able to run up to eight cycles in an experiment. The information on cycle and capillary number is displayed in the Experiment tab.

- **Peak Name Columns** - An individual column per peak name will display for every peak identified by name in the run data. Cells for capillaries in these columns will be blank if Compass did not find peaks automatically using the user-defined peak name analysis parameters (or none were entered).
 - **To view % area in the peak name columns** - the Gaussian Fit area calculation setting must first be selected. To do this, select **Analysis** from the **Edit** menu, then select **Peak Fit** page. In the Peak Find box, select **Gaussian Fit** for the Area Calculation setting (default for Immunoassays). Next, select **Area %** in the upper right corner of the Capillaries pane. This displays the calculated percent area for the named peak compared to all named peaks. This value results from dividing the individual peak area by the sum of all named peak areas for the capillary and multiplying by 100.



Sample	Primary	Capillary	System Co...	Area	Area %	Corr. Area	Conc
HeLa	ERK1+...	C1:2	44704	55899			
HeLa	ERK1+...	C1:3	42886	63846			
HeLa	ERK1+...	C1:4	50291	65303			
HeLa	ERK1+...	C1:5	45954	57780			
HeLa	ERK1+...	C1:6	45887	54196			
HeLa	ERK1+...	C1:7	48500	68243			
HeLa	ERK1+...	C1:8	47257	56120			
HeLa	ERK1+...	C1:9	49193	59797			
HeLa	ERK1+...	C1:10	47426	57923			
HeLa	ERK1+...	C1:11	42898	45543			
HeLa	ERK1+...	C1:12	43542	49701			

- **To view peak area in the peak name columns (default)** - Select **Area** in the upper right corner of the pane. This displays calculated peak area for the individual peak only.
- **To view corrected peak area in the peak name columns** - Select **Corr. Area** in the upper right corner of the pane. This displays the corrected peak area for the named peak compared with the loading control peak.
- **To view concentration in the peak name columns** - Select **Conc.** in the upper right corner of the pane. This displays the calculated concentration for the named peak compared with the standard curve.
- **To view % total in the peak name columns** - The Dropped Lines area calculation setting must first be selected. To do this, select **Analysis** from the **Edit** menu, then select **Peak Fit** page. In the Peak Find box, select **Dropped Lines** for the Area Calculation setting (default for Total Protein Assays). Next, select **% Total** in the upper right corner of the Capillaries pane. This displays the calculated percent total for the named peak compared to the total area measured in the capillary. This value results from dividing the individual peak area by the entire area under the curve for the capillary and multiplying by 100.



Sample	Primary	Capillary	Total Area	DNAK	% Total	Area	% Area
HeLa	Antibo...	C1:2	294092	14.11			
HeLa + 3...	Antibo...	C1:3	395361	33.17			
HeLa + 1...	Antibo...	C1:4	523599	23.36			
HeLa + 7....	Antibo...	C1:5	367106	21.14			
HeLa + 3....	Antibo...	C1:6	344206	17.03			
HeLa	Antibo...	C1:7	257414	14.21			
HeLa + 3...	Antibo...	C1:8	377685	29.93			
HeLa + 1...	Antibo...	C1:9	480506	21.65			
HeLa + 7....	Antibo...	C1:10	483433	20.36			

Viewing Run Data


Each run file contains the following data for up to 96 capillaries:

- **Sample data** - For the proteins in the sample.
- **Standards data** - For the fluorescent standards run with each sample.
- **Registration data (Sally Sue and Peggy Sue only)** - For tracking capillaries as they are moved for various assay steps.


Data can be viewed for one, all, or individually selected capillaries.

Switching Between Sample, Standards and Registration Data Views


You can switch between viewing sample, standards and registration data in a run file using the View bar or View menu:




Data buttons in the View bar:

 Standards : Show Standards

 Reg. : Show Registrations

 Samples : Show Samples

- **To view sample data** - Click **Show Samples** in the View bar or select **View** in the main menu and click **Samples**:




- Data in this view is for sample proteins only.
- Graph view data displays electropherograms in chemiluminescence units (y-axis) and molecular weight in kDa (x-axis).
- Lane view data displays sample proteins only.
- Image view data displays sample proteins only.
- Results for each protein are shown in the peaks and capillaries tables.

NOTE: The reported molecular weight for sample proteins detected by an Immunoassay or a Total Protein Assay in Compass may vary slightly from predicted molecular weights based on sample and assay conditions.

For information on checking and identifying sample peaks, see “Step 4 – Checking the Ladder” on page 183 or “Step 5 – Checking Samples” on page 185.


- **To view standards data** - Click **Show Standards** in the View bar or select **View** in the main menu and click **Standards**:



- Data in this view is for the fluorescent standards only. These standards are premixed with each sample prior to analysis.
- Graph view data displays electropherograms in fluorescence (y-axis) and position (x-axis).
- Lane view data displays fluorescent standards only.
- Image view data displays fluorescent standards only.
- Standards are highlighted in both the peaks and capillaries tables and marked with an **S**.

For information on checking and identifying standards peaks, see “Step 2 – Checking Fluorescent Sizing Standards” on page 178.

- **To view registration data (Sally Sue and Peggy Sue only) - Click Show Registrations in the View bar or select View in the main menu and click Registration:**



- Data in this view shows capillary registration only.
- Graph view data displays electropherograms in fluorescence (y-axis) and position (x-axis).
- Lane view data displays capillary registration only.
- Image view data displays capillary registration only.
- Registrations are highlighted in both the peaks and capillaries tables and marked with an R.


Because capillaries in Sally Sue and Peggy Sue must be moved multiple times during the run, their positions are tracked to account for potential changes in imaging position as they are moved to and from the separation block. To do this, one of the standards is used as a capillary registration. An image of the standard used for registration is taken at the end of the separation step prior to chemiluminescent imaging. The position of this peak in the registration image is then compared to the position of the same standard peak in the standards image. The shift value listed in the peaks table is the change in pixel position between the two images. Compass data is then aligned to account for any changes in capillary position.

For information on checking and identifying registration peaks, see “Step 3 – Checking Capillary Registrations (Sally Sue and Peggy Sue Only)” on page 182.


Selecting and Displaying Capillary Data

You can choose to view data from one, multiple, or all capillaries at once.

- **To look at data for one capillary** - Click a row in the experiment pane. Data for just the row selected will display in all data views and tables. The following example shows sample data displayed in the lane and peaks panes when one sample is selected:



- **To look at data for multiple non-sequential capillaries** - Hold the **Ctrl** key and select the rows you want to view in the experiment pane. Data for only the rows selected will display in all data views and results tables. The following example shows sample data displayed in the graph (in multiple view) and capillaries panes when multiple rows are selected:




- **To look at data for multiple sequential capillaries** - Hold the **Shift** key and select the first and last rows you want to view in the experiment pane. Data for only the selected rows will display in all data views and results tables. The following example shows standards data displayed in the image and peaks panes when multiple sequential rows are selected:

The screenshot shows the Compass software interface with the following components:

- Experiment Pane:** A list of experiment rows. Rows 3 through 8 are selected (highlighted in blue).
- Graph Pane:** A dark image showing four horizontal peaks, each highlighted with a yellow rectangular box.
- Peaks Table:** A table displaying peak data for the selected rows. The table has columns: Sample, Primary, Cap, Peak, Name, Position, MW (kDa), Height, Area, % Area, Width, S/N, and B.


Sample	Primary	Cap	Peak	Name	Position	MW (kDa)	Height	Area	% Area	Width	S/N	B
HeLa	PBK Cl...	C3:7	1	PBK	786	109	1750.5	17322	100.0	9.3	239.8	
HeLa	PBK Cl...	C3:8	1	PBK	776	108	5690.1	56555	100.0	9.3	677.7	
HeLa	PBK Cl...	C3:9	1	PBK	792	108	11855.3	127118	100.0	10.1	1108.4	
HeLa	PBK Cl...	C3:10	1		724	61	222.3	3618		15.3	30.5	
HeLa	PBK Cl...	C3:10	2	PBK	791	107	16793.0	214744	100.0	12.0	1654.2	

- **To look at data for all capillaries** - Hold the **Shift** key and select the first and last rows in the experiment pane. Data for all rows will display in all data views and results tables. The following example shows standards data displayed in the lane and peaks panes when all rows are selected. The pane scroll bars can be used to view all 8 cycles:



Switching Between Single and Multiple Views of the Capillaries

You can switch between displaying run data in the graph, image and lane panes in a per capillary (single) format or up to 96 capillary format. This view is selected in the View bar or the View menu.



Capillary view buttons in the View bar:




Single View




Multiple View

- **To view data in a per capillary format** - Click **Single View** in the View bar or select **View** in the main menu and click **Single View**:




Data for the row(s) selected in the experiment pane will display as follows:


- Electropherograms in the graph pane display overlaid or stacked (depending on the graph option chosen).
- Results for only the selected row(s) will display in the peaks and capillaries tables.
- The selected row(s) are highlighted in the image pane:



- Lanes for only the selected row(s) are displayed in the lane pane:




- **To view data in a multiple capillary format** - Click **Multiple View** in the View bar or select **View** in the main menu and click **Multiple View**:




Data for the row(s) selected in the experiment pane will display as follows:

- Electropherograms corresponding to the selected row(s) are highlighted and displayed in the graph pane.
- Results for selected capillaries display in the peaks and capillaries tables and are highlighted.
- The selected row(s) are highlighted in the image pane:



- All selected lanes display in the lane pane, and lanes corresponding to the selected row(s) are highlighted.




Hiding Capillary Data

Capillary data can be hidden from view if needed. To do this:


1. Click the **Experiment** tab.

- Select the rows you want to hide, then right click on one of the selected rows and click **Hide**:



Data for selected rows will be hidden in all data views and results tables (except for the image pane).

- To view hidden rows** - Select **View** in the main menu and click **Show Hidden**. Hidden rows will become visible again in all panes, and hidden rows will be marked with an X in the experiment pane:




- To unhide rows** - Select the hidden row(s). Right click on one of the selected rows and click **Unhide**.

Setting Run Data Display Filters

The filter lets you auto-hide specific capillaries in the run file data, or only show data for peaks that Compass identified automatically using the user-defined peak name analysis parameters.

NOTE: When more than one run file is open, filter settings will be applied to all files.

- To filter data to show specific capillaries only** - Select **View** in the main menu and click **Filter**. Uncheck the boxes for the capillaries you do not want shown, then click **OK**.



Only data for the checked capillaries will display in data views and results tables, unchecked capillaries will be hidden.

- **To filter data to show named peaks only** - Select **View** in the main menu and click **Filter**. Select **Show named peaks only** then click **OK**. Only data for peaks that Compass identified automatically using the user-defined peak name analysis parameters will display in data views and results tables, all other peak data will be hidden.

Compass Run Data Notifications and Warnings

If Compass detects a potential data issue, a notification or warning will display next to the row in the experiment pane. A list of warnings, notifications and corrective actions are as follows:

- 
Manual correction of sample data notification - Indicates the user modified sample data manually, such as adding or removing a sample peak. Rolling the mouse over the icon displays the type of data modification made.

Sample	Primary	Cycle	
✓ High ...	ERK1/2	1	1
Low p...	ERK1/2	1	1
High ...	ERK1/2		Baseline Manual

- 
Standards warning - Indicates that one of the standards may not be identified properly. This can be resolved by manually identifying the standards. Please refer to “Step 2 – Checking Fluorescent Sizing Standards” on page 178 for details. Rolling the mouse over the icon displays warning details.

Sample	Primary	5	2
Sample	Primary	5	3

Standards Warning: Low Confidence

- 
Manual correction of standards data notification - Indicates the user modified the standards data manually. Rolling the mouse over the icon displays the type of data modification made.

- 
Registrations warning - Indicates that a capillary registration was not identified properly. This can be resolved by manually identifying the registration peak. Please refer to “Step 3 – Checking Capillary Registrations (Sally Sue and Peggy Sue Only)” on page 182 for details. Rolling the mouse over the icon displays warning details.

ERK Hi...	ERK1/2	1	
ERK Hi...	ERK1/2	1	

Registration Warning: Large Registration Shift

-  • **Manual correction of registrations notification** - Indicates the user modified the capillary registration data manually. Rolling the mouse over the icon displays the type of data modification made.
-  • **Peak fit warning** - Indicates that a peak cannot be fit properly. This is typically caused either by a very small peak in the electropherogram or a peak that is very close to the end of the molecular weight range. This can be resolved by removing the peak. Please refer to “Step 4 – Checking the Ladder” on page 183 or “Step 5 – Checking Samples” on page 185 for details. Rolling the mouse over the icon displays warning details.

	Kit low-pho...	anti-H...	2	4
	Kit low-pho...	anti-E...	2	5
Peak Fit Warning: Too many iterations				


Checking Your Results

Compass detects proteins, fluorescent standards and capillary registrations and reports results automatically with no user-interaction. However, we recommend you review the data as outlined in this section as good general practice to ensure accurate results. If a data warning is presented in the experiment pane, the following steps will also help you identify and correct any issues.

Step 1 – Review the Fluorescent Sizing Standards Movie

All capillaries should have three fluorescent sizing standards. To verify the standards separated in the capillary correctly:

1. When the run has completed, click the **Run Summary** screen tab.
2. Click the **Separation** tab and play the movie (this will be the fluorescent standards separation).



- For each capillary, verify that the standards visibly separated. Each fluorescent point in the capillary represents one of the three sizing standards.


Step 2 – Checking Fluorescent Sizing Standards


To verify the standards are identified correctly:

- Click the **Analysis** screen tab.
- Click **Show Standards** in the View bar. Verification that the standards have been correctly identified can be done in either the graph or lane panes.

Graph Pane:

- Click **Single View** in the View bar.
- Click on the first row in the experiment pane, then click the **Graph** tab. Check that the electropherogram has three standard peaks labeled Std 1, Std 29 and Std 230 (for 12-230 kDa size assays) or Std 57 and Std 280 (for 66-440 kDa size assays). They will also be identified with a green **S** in the peaks table.





If standards are not identified correctly, they can be manually corrected as follows:

- **If an incorrect peak is identified as a standard** - Right click the peak in the electropherogram or peaks table and select **Not a Standard**. Compass should correctly reassign the remaining peaks as standards, and update the peaks table.
- **To set an unidentified peak as a standard** - Right click the peak in the electropherogram or peaks table and select **Force Standard**. Compass will assign the peak as a standard, and correctly reassign the remaining standard peaks. A lock icon indicating the standard was set manually will display next to the peak in the peaks table.


*NOTE: To remove standards peak assignments that were made manually, right click on the peak in the electropherogram or peaks table and click **Clear**.*

- Repeat the previous steps for the remaining rows in the experiment pane to make sure all standards are identified correctly.

Lane Pane:

- Click **Multiple View** in the View bar.
- Click on the first row in the experiment pane, then click the **Lane** tab. Standards will be bands and identified with a green outline. Check that all standard bands are labeled: Std 1, Std 12 and Std 230

(for 12-230 kDa size assays) or Std 57 and Std 280 (for 66-440 kDa size assays). They will also be identified with a green **S** in the peaks table. To view band labels, roll the mouse over the individual bands.




If standards are not identified correctly, they can be manually corrected as follows:

- **If an incorrect band is identified as a standard** - Right click the band in the lane or peaks table and select **Not a Standard**. Compass should correctly reassign the remaining bands as standards.
 - **To set an unidentified band as a standard** - Right click the band in the lane or peaks table and select **Force Standard**. Compass will assign the band as a standard, and correctly reassign the remaining standard bands.
- c. Repeat the previous steps for the remaining rows in the experiment pane to make sure all standards are identified correctly.

Step 3 – Checking Capillary Registrations (Sally Sue and Peggy Sue Only)

All capillaries should have a single capillary registration peak. To verify the registrations are identified correctly:

1. Click the **Analysis** screen tab.
2. Click **Show Registrations** and **Single View** in the View bar.
3. Click the **Graph** tab.
4. Click on the first row in the experiment pane. The registration peak will be the first and largest peak in the electropherogram. Check that the registration peak is identified and labeled Reg 1 in the electropherogram. It will also be identified with a purple **R** in the peaks table.



5. If the registration peak is not identified correctly, right click the correct registration peak in the electropherogram or peaks table and select **Registration Force**. Compass will assign the new peak as the capillary registration. A lock icon indicating the registration was set manually will display next to the peak in the peaks table.

*NOTE: To remove registration peak assignments that were made manually, right click on the peak in the electropherogram or peaks table and click **Registration Clear**.*

6. Repeat the previous steps for the remaining rows in the experiment pane to make sure all registrations are identified correctly.


Step 4 – Checking the Ladder

Depending on the size assay you are running, the biotinylated ladder should have the following sizing standards: 12, 40, 66, 116, 180 and 230 kDa (for 12-230 kDa size assays), or 66, 116, 200, 280 and 440 kDa (for 66-440 kDa assays). To verify the ladder standards are identified correctly:

1. Click the **Analysis** screen tab.
2. Click **Show Samples** in the View bar. Verification that the ladder standards have been correctly identified can be done in either the graph or lane panes, but manual corrections must be done in the graph pane.

Graph Pane:

- a. Click **Single View** in the View bar.
- b. Click on the row in the experiment pane that contains the ladder (typically row 1), then click the **Graph** tab. Check that the electropherogram has either six ladder peaks (for 12-230 kDa size assays) or five ladder peaks (for 66-440 kDa size assays). In the example below, the electropherogram has six peaks labeled Ldr 12, Ldr 40, Ldr 66, Ldr 116, Ldr 180 and Ldr 230.



If ladder peaks are not identified correctly, they can be manually corrected as follows:


- **If an incorrect peak is identified as a ladder peak** - Right click the peak in the electropherogram or peaks table and select **Remove peak**. Compass should correctly reassign the remaining peaks as ladder standards.
- **To set an unidentified peak as a ladder peak** - Right click the peak in the electropherogram or peaks table and select **Add Peak**. Compass will assign the peak as a ladder standard, and correctly reassign the remaining ladder standards peaks.

*NOTE: To remove ladder peak assignments that were made manually and go back to the original view of the data, right click in the electropherogram and click **Clear All**.*

Lane Pane:

- Click either **Single View** or **Multiple View** in the View bar.

- b. Click on the row in the experiment pane that contains the ladder (typically row 1), then click the **Lane** tab. Check that the lane has either six ladder bands (for 12–230 kDa size assays) or five ladder bands (for 66–440 kDa size assays). In the example below, the lane has six peaks labeled Ldr 12, Ldr 40, Ldr 66, Ldr 116, Ldr 180 and Ldr 230. To view band labels, roll the cursor over the individual bands. If ladder bands are not identified correctly, they must be corrected in the graph pane as described in the previous section.



Step 5 – Checking Samples

All sample proteins in the graph and lane panes will be labeled automatically with the calculated protein size.


NOTE: The reported molecular weight for immunodetected and total protein sample proteins in Compass may vary slightly from predicted molecular weights based on sample and assay conditions.

To verify that sample proteins are identified correctly:

1. Click the **Analysis** screen tab.
2. Click **Show Samples** in the View bar. Verification that sample proteins have been correctly identified can be done in either the graph or lane panes, but manual corrections must be done in the graph pane.

Graph Pane:

- a. Click **Single View** in the View bar.
- b. Click on the row in the experiment pane that contains the sample you wish to check, then click the **Graph** tab.



If sample peaks are not identified correctly, they can be manually corrected as follows:

- **If an incorrect peak is identified as a sample peak** - Right click the peak in the electropherogram or peaks table and select **Remove peak**. Compass will no longer identify it as a sample peak in the electropherogram and the peak data will be removed in the results tables.


- **To set an unidentified peak as a sample peak** - Right click the peak in the electropherogram or peaks table and select **Add Peak**. Compass will calculate and display the results for the peak in the results tables and identify the peak in the electropherogram.

*NOTE: To remove sample peak assignments that were made manually and go back to the original view of the data, right click in the electropherogram and click **Clear All**.*

- c. Repeat the previous steps for the remaining sample rows in the experiment pane to make sure all sample proteins are identified correctly.

Lane Pane:

- a. Click either **Single View** or **Multiple View** in the View bar.
- b. Click on the row in the experiment pane that contains the sample you wish to check, then click the **Lane** tab. To view sample protein band labels, roll the cursor over the individual bands. If sample bands are not identified correctly, they must be corrected in the graph pane as described in the previous section.



Step 6 – Assigning Peak Names (Optional)

Compass can identify and automatically name sample proteins associated with specific primary antibodies in the run data using user-specified peak name settings. For more information on how to do this see “Peak Names Settings” on page 252.

Group Statistics

The Grouping View is used to analyze replicates by calculating the mean, standard deviation and CV of named proteins (see “Peak Names Settings” on page 252 for detailed information on entering named proteins). Statistics for each protein are also plotted for easy comparison between samples, antibodies, and proteins.

Using Groups


1. A group is automatically created for capillaries with the same sample and primary antibody name, so to use this feature, this information must first be entered into the Assay template. To do this:
 - a. Click the **Assay** tab and go to the Template pane.
 - b. Enter sample names and primary antibody names as described in “Step 6 - Add Assay Plate Annotations (Optional)” on page 36. Be sure to enter the same sample and/or primary antibody names for the groups of samples you want to calculate statistics for.

In this example there are two sample types, **Sample A** and **Sample B** which were run with two different antibodies, **Primary 1** and **Primary 2**.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Biot. La...	SampleA	SampleB	SampleA	SampleB	SampleA	SampleB	SampleA	SampleB		Control	
B	Antibody Diluent											
C	Antibody...	Primary 1	Primary 2	Primary 1	Primary 2	Primary 1	Primary 2	Primary 1	Primary 2	Primary 1	Primary 2	Primary 1
D	Streptav...	Secondary Antibody										
J	Detection											

Each of the two samples were run with each of the two antibodies twice in every cycle. This creates four groups for the combination of two samples and two antibodies.

2. Select the **Analysis** tab. Group your results and view the associated statistics by selecting **View** and clicking on **Grouping**. In the box, click **Enable Grouping**.



3. Select a grouping option by clicking the box next to the option. These options allow you to group capillaries in multiple ways:



- **Group across runs** - Groups capillaries from multiple runs. This option creates fewer groups and generates statistics across multiple runs.
- **Group across cycles** - Groups capillaries run in different cycles.
- **No option selected** - When only one run is open, groups will only contain capillaries from the same cycle. When more than one run is open, groups are created for each run.

Viewing Statistics

Peak and Capillary Groups

The Peak Groups pane reports statistics for each named protein in every group. Each group shows the statistics for named proteins which includes average area, standard deviation and %CV. The number in parenthesis after the sample name indicates the number of capillaries in the group.

Sample	Primary	Cap	Name	Corr Area	Std.Dev.	% CV	SEM
▶ HeLa (24)	ERK1+System Control		ERK1	127918	12886	10.1	2630
▶ HeLa (24)	ERK1+System Control		System Control	100000	0.0000	0.0	0.0000

Clicking the arrow next to a group lists the individual capillaries in the group and reported data for each capillary:

Sample	Primary	Cap	Name	Corr Area	Std.Dev.	% CV	SEM
HeLa (24)	ERK1+System Control		ERK1	127918	12886	10.1	2630
HeLa	ERK1+System Control	C1:2	ERK1	125041			
HeLa	ERK1+System Control	C1:3	ERK1	148874			
HeLa	ERK1+System Control	C1:4	ERK1	129850			
HeLa	ERK1+System Control	C1:5	ERK1	125733			
HeLa	ERK1+System Control	C1:6	ERK1	118107			
HeLa	ERK1+System Control	C1:7	ERK1	140706			
HeLa	ERK1+System Control	C1:8	ERK1	118756			

The **Capillary Groups** pane pivots the **Peak Groups** results to show statistics for named protein peaks in individual columns.

Sample	Primary	Capillary	Syste...	Std.Dev.	%CV	SEM	ERK1	Std.Dev.	%CV	SEM
HeLa (24)	ERK1+System Control		100000	0.0000	0.0	0.0000	127918	12886	10.1	2630

Group Plots


The mean values for named peaks in each group are plotted in bar graphs with error bars showing the standard deviation. The plots compare different antibodies for the same sample and different samples for the same antibody to allow a choice of presentation. The y-axis on the plot will be concentration if a standard curve is defined.



Hiding or Removing Capillaries in Group Analysis

Hidden capillaries are not included in groups. However, hiding capillaries provides an easy way to reject individual capillaries from the statistical analysis. See “Hiding Capillary Data” on page 173 for details on how to do this.

Larger groups of capillaries can also be excluded from analysis. To do this, select **View** and click **Filter**.




Uncheck the box next to the cycles or capillaries you wish to remove and click **OK**. This data will now be removed from the grouped view statistics.

Copying Data Views and Results Tables

You can copy and paste data and results tables into other documents, or save a data view as a graphic file.

Copying Data Views

1. Click in the graph or lane pane.
2. Select **Edit** in the main menu and click **Copy**, or right click and select **Copy**.
3. If you selected copy from the graph pane the following window will display, click **Copy**. If you selected copy from the lane pane skip to the next step.




4. Open a document (Microsoft® Word®, Excel®, PowerPoint®, etc.). Right click in the document and select **Paste**. A graphic of the copied data view will be pasted into the document.

Copying Results Tables

1. Click in the peaks or capillaries pane.
2. Select one or multiple rows.
3. Select **Edit** in the main menu and click **Copy**, or right click and select **Copy**.
4. Open a document (Microsoft® Word®, Excel®, PowerPoint®, etc.). Right click in the document and select **Paste**. Data for the rows selected will be pasted into the document.

Saving the Graph View as an Image File

1. Click in the graph pane.
2. Select **Edit** in the main menu and click **Copy**, or right click and select **Copy**.
3. Select an image option (EMF, PNG or PDF) in the pop-up window, then click **Save**.



4. Select a directory to save the file to and enter a file name, then click **OK**.

Exporting Run Files

Results tables and raw plot data can be exported for use in other applications.

Exporting Results Tables

To export the information in the peaks and capillaries tables:


1. Click **File** in the main menu and click **Export Tables**.
2. Select a directory to save the files to and click **OK**. Data will be exported in .txt format.

NOTE: To exclude export of standards data or export results table data in .csv format, see "Setting Data Export Options" on page 413.

Exporting Raw Sample Electropherogram Data

To export raw sample plot data:

1. Click **File** in the main menu and click **Export Spectra**.




- **To export data in .txt format** - Select **Text Format**. Plots will be exported in one file for all capillaries.
 - **To export data in .cdf format** - Select **Andi Format**. Plots will be exported in one file per capillary.
2. Select a directory to save the files to and click **OK**. Data will be exported in the selected format.

Changing Sample Protein Identification

Compass allows you to customize what sample proteins are reported in the results tables by making manual adjustments in the electropherogram or peaks table.

Adding or Removing Sample Data

1. Click **Show Samples** in the View bar.
2. Click **Single View** in the View bar.
3. Click on the row in the experiment pane that contains the sample you wish to correct, then click the **Graph** tab.
 - **To remove a peak from the data** - Right click the peak in the electropherogram or peaks table and select **Remove peak**. Compass will no longer identify it as a sample peak in the electropherogram and the peak data will be removed in the results tables.



- **To add an unidentified peak to the data** - Right click the peak in the electropherogram or peaks table and select **Add Peak**. Compass will calculate and display the results for the peak in the results tables and identify the peak in the electropherogram.

NOTES:


To remove sample peak assignments that were made manually and go back to the original view of the data, right click in the electropherogram and click **Clear All**.

Virtual blot data in the lane pane will also update to reflect changes made in the graph pane.


Hiding Sample Data

You can hide the results for a sample protein in the results tables without completely removing it from the reported results. To do this:

1. Click **Show Samples** in the View bar.
2. Click **Single View** in the View bar.
3. Click on the row in the experiment pane that contains the sample you wish to correct, then click the **Graph** tab.
4. Right click the peak in the electropherogram or peaks table and select **Hide peak**. Compass will hide the peak data in the results tables.



- To view hidden peak data, click **View** in the main menu and click **Show Hidden**. Hidden peak data will display in the results table and be marked with an **X**.



- To unhide a peak, right click on the peak in the electropherogram or peaks table and select **Unhide Peak**.

Changing Peak Names for Sample Data

If Compass did not automatically name a sample protein associated with a specific primary antibody, this can be adjusted manually. To do this:

- Click **Show Samples** in the View bar.
- Click **Single View** in the View bar.
- Click on the row in the experiment pane that contains the sample you wish to correct, then click the **Graph** tab.

- Right click the peak in the electropherogram or peaks table and click **Name**, then click a name in the list. Compass will change the peak name in the electropherogram and results tables, and adjust peak names for other sample proteins accordingly.



NOTES:


For details on how to specify peak name settings, see "Peak Names Settings" on page 252.

Virtual blot data in the lane pane will also update to reflect changes made in the graph pane.


Displaying Sample Data for Named Peaks Only

You can adjust the sample data to display results only for user-specified named peaks. To do this:

1. Click **Show Samples** in the View bar.
2. Click **View** in the main menu and click **Filter**.
3. Check the **Show Named Peaks only** box and click **OK**.



Compass will hide peak data in the electropherogram and results tables for all unnamed peaks, and instead only display data for named peaks in the electropherogram and results tables.



Changing the Virtual Blot View

Options in the lane pane let you change the contrast or invert the virtual blot, remove baseline noise, change lane labels or overlay standards data on sample lanes.

The lane pane toolbar has the following options:

-  Contrast Adjustment
-  Invert
-  Edit Labels
-  Remove Baseline
-  Overlay Standards Data

Adjusting the Contrast


1. Click the **Contrast Adjustment** button. The contrast tool will display:



2. Click the bar and drag it up or down to adjust the contrast.
3. When finished, click **X** to close the tool.

Inverting the Virtual Blot

1. Click the **Invert** button. The virtual blot image will invert:



2. Click the **Invert** button again to return to the default view.

Selecting Lane Labels

The labels shown above the lanes in the virtual blot can be customized. To do this:

1. Click the **Edit Labels** button. The label box will display:




2. Check one or multiple label boxes, and uncheck those you don't want to display. To remove labels completely, uncheck all boxes. The following label options are available:
 - **Sample** - Sample name. If sample names were entered in the assay template (Assay screen), those names will display here. Otherwise, Sample (default name) will display.
 - **Blocking** - Blocking reagent name. If a name was entered in the assay template (Assay screen), that name will display here. Otherwise, Blocking (default name) will display.
 - **Primary** - Primary antibody name. If primary antibody names were entered in the assay template (Assay screen), those names will display here. Otherwise, Primary (default name) will display.
 - **Secondary Ab** - Secondary antibody name. If a secondary antibody name was entered in the assay template (Assay screen), those names will display here. Otherwise, Secondary (default name) will display.
 - **Capillary** - Cycle and capillary number. For example, C1:3 indicates cycle 1, capillary 3

NOTE: Sally Sue and Peggy Sue run up to eight cycles (12 capillaries at a time), so the cycle number displayed will always be C1-C8 depending on the number of cycles programmed.

- **Attributes** - Attribute text. If attribute information was entered for any of the rows or wells in the assay template (Assay screen), up to four can be selected and displayed as lane labels.

Viewing the Uncorrected Sample Baseline

1. Click the **Remove Baseline** button (active for sample data only). This will remove the automatic baseline correction.



2. Click the **Remove Baseline** button again to return to the default, baseline corrected view.

Overlaying Standards Data on Sample Lanes

Data for the standards can be overlaid on the sample data in the virtual blot so you can view the raw, unaligned (uncorrected) lane data. To do this:

1. Click the **Overlay Standards Data** button (active for sample data only). An overlay of the raw sample and standards data will display:



2. Click the **Overlay Standards Data** button again to return to the default, auto-aligned view.

Moving Lanes in the Virtual Blot View

Lanes in the virtual blot view can be moved to make it easier to compare specific lanes of data. To do this:

1. Click a lane in the virtual blot. Hold the mouse button down and drag the lane to a new position.



2. Release the mouse button. The lane will now be repositioned in the virtual blot view.

Changing the Electropherogram View

Options in the graph pane let you zoom and scale electropherograms, overlay or stack plots and change the peak and plot information displayed.

The graph pane toolbar has the following options:



Auto Scale



Graph Options




Stack the Plots




Overlay the Plots

Autoscaling the Electropherogram

Click the **Autoscale** button to scale the y-axis to the largest peak in the electropherogram.




Click the **Autoscale** button again to return to default scaling.



Stacking Multiple Electropherograms

Electropherogram data for multiple capillaries can be stacked vertically in the graph pane for comparison. To do this:

1. Click **Single View**.
2. Select multiple rows in the experiment pane.
3. Click the **Stack the Plots** button. The individual electropherograms for each row selected will stack in the graph pane.




You can customize the colors used for the stacked plot display. For more information see "Selecting Custom Plot Colors for Graph Overlay" on page 414.

Overlaying Multiple Electropherograms

The electropherogram data for multiple capillaries can also be overlaid on top of each other for comparison in the graph pane. To do this:


1. Click **Single View**.
2. Select multiple rows in the experiment pane.
3. Click the **Overlay the Plots** button. The individual electropherograms for each row selected will overlay in the graph pane.




You can customize the colors used for the overlay plot display. For more information see “Selecting Custom Plot Colors for Graph Overlay” on page 414.

Zooming

To zoom in on a specific area of the electropherogram, hold the mouse button down and draw a box around the area with the mouse:




To return to default scaling, right click in the electropherogram and click **Zoom Out**.




Customizing the Data Display

You can customize electropherogram peak labels, plot labels and display options. To do this, select the **Graph Options** button.




Peak Labels

You can customize the labels used to identify peaks in the electropherogram with these options:




Matching Peak Names - Checking this box will draw vertical lines through each named peak. Using this option with Stack the Plots or Overlay the Plots features is useful for visual comparison of named peaks across multiple capillaries.



- **Peak Names** - Checking this box will display peak name labels on all named peaks in the electropherogram.

NOTE: If more than one peak label option is selected, peak name labels will always be used for named peaks.




- **Peak Values** - Checking this box will display the molecular weight labels on all peaks in the electropherogram.

NOTES:

When viewing standards and registration data, this option is called Peak Positions. Labels displayed are peak positions rather than molecular weight.


If more than one peak label option is selected, peak name labels will always be used for named peaks.

The reported molecular weight for immunodetected and total protein sample proteins in Compass may vary slightly from predicted molecular weights based on sample and assay conditions.




Baseline and Grid Options

You can view the calculated baseline fit, peak integration and show grid lines with these options.




- **Fitted peaks** - Checking this box will display how the peaks were fit by the software.
 - For Immunoassays, the software uses Gaussian fit by default:

NOTE: This option is only available for sample data.




- For Total Protein Assays, the software uses Dropped Lines fit by default:




- **Baseline Fit** - Checking this box will display the calculated baseline for the peaks. Baseline points will also display for regions of the electropherogram considered to be at baseline.

NOTE: This option is only available for sample data.



- **Grid Lines** - Checking this box will display grid lines in the graph area.



Plot Labels


You can customize the plot labels displayed on the electropherogram with these options.

Plot Label	
<input checked="" type="checkbox"/> Sample	<input type="checkbox"/> Attribute
<input type="checkbox"/> Primary	<input type="checkbox"/> Attribute
<input checked="" type="checkbox"/> Capillary	<input type="checkbox"/> Exposure

Plot labels are shown on the right side of the graph pane.


- **Sample** - Checking this box will display the sample name. If sample names were entered in the assay template (Assay screen), those names will display here. Otherwise, Sample (default name) will display.
- **Primary** - Checking this box will display the primary antibody name. If primary antibody names were entered in the assay template (Assay screen), those names will display here. Otherwise, Primary (default name) will display.
- **Capillary** - Checking this box will display the cycle and capillary number. For example, C1:3 indicates cycle 1, capillary 3.
- **Attributes** - Checking this box will display attribute text. If attribute information was entered for samples or primary antibodies in the assay template (Assay screen), they can be selected as plot labels.

- **Exposure** - Checking this box will display the exposure time(s) used for the data.
The following example shows an electropherogram with all plot labels selected:



Selecting Data Viewing Options


The graph view menu provides you with multiple options for changing what type of electropherogram data is displayed. To access the view menu, click the down arrow next to the graph pane toolbar:




A check mark next to the menu option indicates it is currently selected. Multiple options can be selected at a time.

NOTE: Unless otherwise noted, graph view menu options are available for sample data only.

- **Sample** - Clicking this option will display raw, uncorrected sample data.




- **Sample Baseline Corrected** - Clicking this option will display sample data with the baseline subtracted (zeroed). This is the default view.




- **Baseline Fit** - Clicking this option will display the calculated baseline for the raw sample data. In the example that follows, both Baseline fit and Sample are selected.

NOTE: This option is selected automatically when Baseline Fit is selected in graph options.




- **Baseline Points** - Clicking this option will display regions of the electropherogram considered to be at baseline. In the example that follows, both Baseline Points and Sample are selected.


NOTE: This option is selected automatically when Baseline Fit is selected in graph options.



- **Fit** - Clicking this option will display the bounding envelope of the fitted peaks as calculated by the software for the raw sample data. In the example that follows, both Fit and Sample are selected.




- **Fit Baseline Corrected** - Clicking this option will display the fitted peaks as calculated by the software for the sample baseline corrected data. In the example that follows, both Fit Baseline Corrected and Sample are selected, the fit plot is on the bottom.



NOTE: When viewing standards or registration data, this option is called Standards Fluorescence or Registration Fluorescence, respectively.

- **Standards** - Checking this box aligns the molecular weight of the raw standards data to the sample data and overlays both electropherograms in the graph pane.




Adding and Removing Baseline Points

Points in the baseline can be added or removed as needed. To do this:

1. Click the **Graph Options** button in the graph pane toolbar and check **Baseline Fit**. This will display baseline points for the raw sample data.
2. Use the mouse to draw a box around the area you want to correct. This will zoom in on the area.


3. Right click a baseline point and click **Add Baseline Point** or **Remove Baseline Point**.




*NOTE: To clear the manual addition or removal of baseline points and go back to the original view of the data, right click in the electropherogram and click **Clear All**.*

Selecting the X-Axis Molecular Weight Range


The molecular weight range used for the x-axis can be changed. To do this, Select **View** in the main menu and click **View Region**. The following pop-up window will display:




- Selecting **Analysis** will match the range of the run data that was selected in the Peak Fit Analysis Settings in both the electropherogram and virtual blot view.



- Selecting **Full** will display the range of the run data for the entire capillary in both the electropherogram and virtual blot view. (Insert View region Full)



- Selecting **Custom** allows you to manually enter the Lower and Upper ranges of the run data to be displayed in both the electropherogram and virtual blot view.(Insert View region Custom)



NOTE: You can change the default x-axis range that Compass uses. For more information, see "Peak Fit Analysis Settings" on page 244.


Closing Run Files

If more than one run file is open, you can close just one file or all the open files at the same time.


- **To close one of multiple open run files** - In the experiment pane, click on one of the sample rows in the file. Next click **File** from the main menu and click **Close**.
- **To close all open run files** - Select **File** from the main menu and click **Close All**.

Compass Analysis Settings Overview

Compass has a variety of analysis features and settings that you can modify as needed to enhance run data. To access these settings, select **Edit** in the main menu and click **Analysis**. If more than one run file is open, select the run file you want to view analysis settings for from the list:



The following screen will display:




To move between analysis pages in this window, click on an option in the list on the left. The following analysis settings can be user-customized in Compass:

- **Advanced** - Lets you customize analysis settings for samples, standards and image data.
- **Images** - Lets you view the chemiluminescent exposures taken during the run and view data for different exposures in the Analysis screen.
- **Peak Fit** - Lets you customize peak fit settings for sample data.
- **Peak Names** - Lets you enter custom naming settings for sample proteins associated with specific blocking reagents, primary antibodies or attributes and have Compass automatically label the peaks in the run data.
- **Standards** - Lets you customize the molecular weight and positions Compass uses to identify ladder standards, fluorescent standards and registration peaks, as well as change the capillary used for the ladder.

Advanced Analysis Settings

The advanced analysis settings page lets you view and change analysis settings for samples, standards and image data. To access these settings, select **Edit** in the main menu and click **Analysis**, then click **Advanced** in the options list.

NOTE: Settings can be modified in an assay prior to starting a run, or in a run file once the run has finished executing. Analysis changes made to an executing run will not be saved to the final run file.



- Click **Import** to import an analysis settings file. This will be explained in more detail in “Importing Analysis Settings” on page 284.
- Click **Export** to export the current analysis settings file. This will be explained in more detail in “Exporting Analysis Settings” on page 284.
- Click **Restore Original** to restore Compass default settings.
- Click **OK** to save changes and exit.
- Click **Cancel** to exit without saving changes.

Standards Settings

- **Peak Width** - The approximate width (at full width half max) used to filter out fluorescence artifacts which improves recognition of standards.
- **Allowable Drift** - The distance in pixels that standards are expected to move compared to their entered positions in the standards analysis page. This setting assists in recognition of standards.

Sample Settings

- **Peak Fit Starting Width Ratio** - Focuses the peak fit towards the peak center and aids the overall peak fitting.

Image Settings

- **Median Filter Threshold Ratio** - Pixel ratio used to filter out camera artifacts.
- **Median Filter Threshold Limit** - Pixel threshold value used to filter out camera artifacts.

Advanced Analysis Settings Groups


Advanced analysis settings are saved as a group, and multiple settings groups can be created. Specific group settings can then be applied to individual capillaries, sample names or attributes in the run data.

NOTES:

We recommend using the Compass default values for advanced analysis settings. These settings are included in the default Advanced group.

Analysis settings are run-file specific. However, settings can be imported or exported for use with other run files. For more information see "Importing and Exporting Analysis Settings" on page 284.

Analysis groups are displayed in the analysis settings box:




The Advanced group shown contains the Compass default analysis settings. You can make changes to this group and create new groups.


To view settings for a group, click on the group name in the analysis settings box.

Creating a New Analysis Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Advanced** in the options list.
2. Click **Add** under the analysis settings box. A new group will be created:



3. Click on the new group and enter a new name.




4. Modify standards, sample or image parameters as needed.
5. To use the new group as the default analysis settings for the run file data, click the arrow in the drop down list next to Default, then click the new group from the list. Analysis settings in the new group will then be applied to the run data.



6. Click **OK** to save changes.

Changing the Default Analysis Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Advanced** in the options list.
2. Click the arrow in the drop down list next to Default, then click a new default group from the list.



3. Click **OK** to save changes. Analysis settings in the group selected will be applied to the run data.

Modifying an Analysis Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Advanced** in the options list.
2. Click on the group in the analysis settings box you want to modify.



3. Modify standards, sample or image parameters as needed.
4. Click **OK** to save changes. The new analysis settings will be applied to the run data.

Deleting an Analysis Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Advanced** in the options list.
2. Click on the group in the analysis settings box you want to delete and click **Remove**.




3. Click **OK** to save changes.

Applying Analysis Groups to Specific Run Data


1. Select **Edit** in the main menu and click **Analysis**, then click **Advanced** in the options list.
2. Click on the group in the analysis settings box you want to apply to specific run data.



3. Application of analysis groups to specific run data is done in the override box. Click **Add** under the override box. A default override data set will be created from sample information found in the run file.




4. Click the cell in the **Apply To** column, then click the down arrow.




5. Select an option from the drop down list. This will apply the settings group selected to specific run data as follows:
- **Sample names** - All sample names entered in the assay template (Assay screen) will display in the drop down list, otherwise the default name of Sample will be shown. Select a sample name to apply group settings to all capillaries that use this sample name in the run file.
 - **Attributes** - All sample attributes entered in the assay template (Assay screen) will display in the drop down list. Select an attribute to apply group settings to all capillaries that use this attribute in the run file.

- **Cycle 1-8** - When this option is selected, group settings will be applied to all capillaries in the cycle.
- **Custom settings** - Lets you choose specific capillaries to apply the group settings to. When this option is selected, the following pop-up box appears to let you enter a specific capillary number or range of capillaries:




6. If you need to change the analysis group used for a data set, click the cell in the **Settings** column and click the down arrow. Select a group from the drop down list.



7. Repeat the previous steps to apply other groups to specific run data.
8. To remove a data set, click on its cell in the **Apply To** column, then click **Remove**.
9. Click **OK** to save changes.

Images Analysis Settings

The Images analysis settings page lets you see what chemiluminescent exposures were taken during the run, and view data for different exposures in the Analysis screen. To access these settings, select **Edit** in the main menu and click **Analysis**, then click **Images** in the options list.



- Click **Import** to import an analysis settings file. This will be explained in more detail in “Importing Analysis Settings” on page 284.
- Click **Export** to export the current analysis settings file. This will be explained in more detail in “Exporting Analysis Settings” on page 284.
- Click **Restore Original** to restore Compass default settings.
- Click **OK** to save changes and exit.
- Click **Cancel** to exit without saving changes.

Exposure Settings


The exposure used for the sample data displayed in the Analysis screen is shown in the All box:

Cycle	Luminescence
All	Multi-Image Analysis
1	Multi-Image Analysis
2	Multi-Image Analysis
3	Multi-Image Analysis
4	Multi-Image Analysis
5	Multi-Image Analysis
6	Multi-Image Analysis
7	Multi-Image Analysis
8	Multi-Image Analysis

NOTE: Sally Sue and Peggy Sue run up to eight cycles (12 capillaries per cycle). The individual exposure setting selected for a specific cycle is applied for all 12 capillaries in that cycle unless selected for All cycles.

- **Multi-Image Analysis** - Sample data displayed in the Analysis screen is compiled from all exposures taken during the run and utilized to calculate the chemiluminescent signal output at time zero of the chemiluminescent reaction. This calculation represents the kinetics of the chemiluminescent reaction and helps to eliminate signal burn out that may occur with stronger signals and longer exposure times.
- **Exposure #** - Sample data displayed in the Analysis screen is for this specific exposure only.

To see the number of exposures and exposure times used for the run data, click the arrow in the drop down list next to **All**.




NOTE: The number of exposures taken and exposure times shown were specified in the assay protocol (Assay screen) prior to executing the run and cannot be modified.

Changing the Sample Data Exposure Displayed

The exposure used for the sample data displayed in the Analysis screen can be changed. To do this:

1. Select **Edit** in the main menu and click **Analysis**, then click **Images** in the options list.
2. Click the arrow in the drop down list next to **All** and select an exposure setting:




3. Click **OK** to save changes and exit. Sample data for the exposure selected will display in the Analysis screen.

Peak Fit Analysis Settings

The peak fit analysis settings page lets you view and change peak fit settings for sample data. To access these settings, select **Edit** in the main menu and click **Analysis**, then click **Peak Fit** in the options list.

NOTE: Settings can be modified in an assay prior to starting a run, or in a run file once the run has finished executing. Analysis changes made to an executing run will not be saved to the final run file.



- Click **Import** to import an analysis settings file. This will be explained in more detail in "Importing Analysis Settings" on page 284.

- Click **Export** to export the current analysis settings file. This will be explained in more detail in “Exporting Analysis Settings” on page 284.
- Click **Restore Original** to restore Compass default settings.
- Click **OK** to save changes and exit.
- Click **Cancel** to exit without saving changes.

Range Settings


- **Minimum** - The molecular weight value (in kDa) below which peaks will not be identified. This value will also be used as the default lower MW range for the data displayed in the electropherogram and virtual blot.
- **Maximum**: The molecular weight value (in kDa) above which peaks will not be identified. This value will also be used as the default upper MW range for the data displayed in the electropherogram and virtual blot.

Baseline Settings


- **Threshold** - The variance, or roughness, in a baseline data segment below which a point will be called part of the baseline.
- **Window** - How long baseline data segments are expected to be in pixels. Shorter segments will allow the baseline to follow plateau sections of the signal.
- **Stiffness** - The amount the baseline is allowed to vary from a straight line. Settings between 0.1 and 1.0 make the baseline fit closer to a straight line. Settings from 1.0 to 10.0 will make the baseline fit follow the data more closely.

Peak Find Settings

- **Threshold** - The minimum signal to noise ratio required for a peak to be identified. A setting of 1.0 will detect many peaks, a setting of 10.0 will detect fewer peaks.
- **Width** - The approximate peak width (at full width half max) in pixels used to detect peaks. The minimum value for this setting is 3.0. Larger widths help to eliminate the detection of shoulder peaks and noise peaks.
- **Area Calculation** - Two fits are used, either Gaussian Fit or Dropped Lines. These settings can be changed before or after the run is finished.
 - For Immunoassays, peak area is calculated using Gaussian distribution by default:



- For Total Protein Assays, peak area is calculated using Dropped Lines. This type of area calculation is also often called the perpendicular drop method. This method is preferred when peaks overlap or are close to each other as they are in Total Protein Assays. This method draws two vertical lines from the left and right bounds of the peak down to the x-axis and then measures the total area bounded by the signal curve, the x-axis (y=0 line), and the two vertical lines.



Peak Fit Analysis Settings Groups


Peak fit settings are saved as a group, and multiple settings groups can be created. Specific group settings can then be applied to individual capillaries, sample names or attributes in the run data.

NOTES:

We recommend using the Compass default values for peak fit analysis settings. These settings are included in the default Peak Fit group.

Analysis settings are run-file specific. However, settings can be imported or exported for use with other run files. For more information see "Importing and Exporting Analysis Settings" on page 284.

Peak fit groups are displayed in the analysis settings box:




The Peak Fit group shown contains the Compass default analysis settings. You can make changes to this group and create new groups.


To view settings for a group, click on the group name in the analysis settings box.

Creating a New Peak Fit Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Fit** in the options list.
2. Click **Add** under the analysis settings box. A new group will be created:



3. Click on the new group and enter a new name.




4. Modify range, baseline or peak find parameters as needed.
5. To use the new group as the default peak fit settings for the run file data, click the arrow in the drop down list next to Default, then click the new group from the list. Peak fit settings in the new group will then be applied to the run data.



6. Click **OK** to save changes.

Changing the Default Peak Fit Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Fit** in the options list.
2. Click the arrow in the drop down list next to Default, then click a new default group from the list.



3. Click **OK** to save changes. Peak fit settings in the group selected will be applied to the run data.

Modifying a Peak Fit Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Fit** in the options list.
2. Click on the group in the analysis settings box you want to modify.



3. Modify range, baseline or peak find parameters as needed.
4. Click **OK** to save changes. The new peak fit settings will be applied to the run data.

Deleting a Peak Fit Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Fit** in the options list.
2. Click on the group in the analysis settings box you want to delete and click **Remove**.




3. Click **OK** to save changes.

Applying Peak Fit Groups to Specific Run Data


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Fit** in the options list.
2. Click on the group in the analysis settings box you want to apply to specific run data.



- Application of peak fit groups to specific run data is done in the override box. Click **Add** under the override box. A default override data set will be created from sample information found in the run file.




- Click the cell in the **Apply To** column, then click the down arrow.




- Select an option from the drop down list. This will apply the settings group selected to specific run data as follows:
 - Sample names** - All sample names entered in the assay template (Assay screen) will display in the drop down list, otherwise the default name of Sample will be shown. Select a sample name to apply group settings to all capillaries that use this sample name in the run file.

- **Attributes** - All sample attributes entered in the assay template (Assay screen) will display in the drop down list. Select an attribute to apply group settings to all capillaries that use this attribute in the run file.
- **Cycle 1-8** - When this option is selected, group settings will be applied to all capillaries in the cycle.
- **Custom settings** - Lets you choose specific capillaries to apply the group settings to. When this option is selected, the following pop-up box appears to let you enter a specific capillary number or range of capillaries:



6. If you need to change the peak fit group used for a data set, click the cell in the **Settings** column and click the down arrow. Select a group from the drop down list.




7. Repeat the previous steps to apply other groups to specific run data.
8. To remove a data set, click on its cell in the **Apply To** column, then click **Remove**.
9. Click **OK** to save changes.

Peak Names Settings

The peak names analysis settings page lets you enter custom naming settings for sample proteins. Compass can automatically label sample peaks in the run data associated with specific blocking reagent names, primary antibody names or attributes. To access these settings, select **Edit** in the main menu and click **Analysis**, then click **Peak Names** in the options list.

NOTE: Settings can be modified in an assay prior to starting a run, or in a run file once the run has finished executing. Analysis changes made to an executing run will not be saved to the final run file.



- Click **Import** to import an analysis settings file. This will be explained in more detail in “Importing Analysis Settings” on page 284.


- Click **Export** to export the current analysis settings file. This will be explained in more detail in “Exporting Analysis Settings” on page 284.
- Click **Restore Original** to restore Compass default settings.
- Click **OK** to save changes and exit.
- Click **Cancel** to exit without saving changes.

Peak Names Analysis Settings Groups

Peak name settings are saved as a group, and multiple settings groups can be created. Specific group settings can then be applied to individual capillaries, blocking reagent names or primary antibody names and attributes in the run data.

NOTE: Analysis settings are run-file specific. However, settings can be imported or exported for use with other run files. For more information see “Importing and Exporting Analysis Settings” on page 284.

Peak name groups are displayed in the analysis settings box:




The Peak Names group shown is a Compass template group. You can make changes to this group and create new groups.


To view settings for a group, click on the group name in the analysis settings box.

Creating a Peak Names Group

1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Names** in the options list.
2. Click on the **Peak Names 1** template group in the analysis settings box.



3. Enter a new name for the group.



4. Click in the first cell in the **Name** column in the analysis settings peak table.

Analysis Settings: ERK1/2

Name	MW	Color	Range (%)	Control	Show
Peak1	100		10	<input type="checkbox"/>	<input checked="" type="checkbox"/>

5. Enter a sample protein name associated with the primary antibody used in the run.

Analysis Settings: ERK1/2

Name	MW	Color	Range (%)	Control	Show
ERK2	100		10	<input type="checkbox"/>	<input checked="" type="checkbox"/>

6. Click in the first cell in the **MW** column.

Analysis Settings: ERK1/2

Name	MW	Color	Range (%)	Control	Show
ERK2	100		10	<input type="checkbox"/>	<input checked="" type="checkbox"/>

7. Enter the molecular weight (in kDa) for the sample protein.

Analysis Settings: ERK1/2


Name	MW	Color	Range (%)	Control	Show
ERK2	42		10	<input type="checkbox"/>	<input checked="" type="checkbox"/>

8. Click in the first cell in the **Color** column, then click the button.

Analysis Settings: ERK1/2

Name	MW	Color	Range (%)	Control	Show
ERK2	42		10	<input type="checkbox"/>	<input checked="" type="checkbox"/>

The color selection box will display:



The color selected will be used to identify the sample protein peak in the peaks and capillaries table in the Analysis screen.

9. Click a color or define a custom color and click **OK**. The color selection will update in the table:

Analysis Settings: ERK1/2

Name	MW	Color	Range (%)	Show
ERK2	42		10	<input checked="" type="checkbox"/>

10. Click in the first cell in the **Range (%)** column.

Analysis Settings: ERK2/2

Name	MW	Color	Range (%)	Control	Show
ERK2	42		10	<input type="checkbox"/>	<input checked="" type="checkbox"/>

- Enter a range window for the MW entered. Compass will automatically name peaks found within this percent of the molecular weight. For example, if the molecular weight entered is 40 kDa and a 10% range is used, all peaks between 36 and 44 kDa will be identified with this peak name.

NOTE: The reported molecular weight for immunodetected and total protein sample proteins in Compass may vary slightly from predicted molecular weights based on sample and assay conditions.

- Select the checkbox in the first cell of the **Show** column. This will turn peak naming on for the sample protein.

Analysis Settings: Protein


Name	MW	Color	Range (%)	Control	Show
ERK2	42		10	<input type="checkbox"/>	<input checked="" type="checkbox"/>

To turn peak naming off for a particular sample protein, deselect the checkbox in the Show column.

- To add another sample protein, click **Add** under the analysis settings peak table:

Analysis Settings: Protein


Name	MW	Color	Range (%)	Control	Show
ERK2	42		10	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Peak2	100		10	<input type="checkbox"/>	<input checked="" type="checkbox"/>



4. Enter information in the analysis settings peak table as described in “Creating a Peak Names Group” on page 253.
5. Click **OK** to save changes.

Modifying a Peak Names Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Names** in the options list.
2. Click on the group in the analysis settings box you want to modify.



3. Change the information in the analysis settings peak table as described in “Creating a Peak Names Group” on page 253.
4. Click **OK** to save changes.

Deleting a Peak Names Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Names** in the options list.
2. Click on the group in the analysis settings box you want to delete and click **Remove**.




3. Click **OK** to save changes.

Applying Peak Names Groups to Run Data


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Names** in the options list.
2. Click on the group in the analysis settings box you want to apply to specific run data.




3. Application of peak names groups to specific run data is done in the override box. Click **Add** under the override box. A default override data set will be created.




4. Click the cell in the **Apply To** column, then click the down arrow.



5. Select an option from the drop down list. This will apply the settings group selected to specific run data as follows:
- **All** - When this option is selected, peak names group settings will be applied to all capillaries in the run file.
 - **Blocking reagent** - When this option is selected, group settings will be applied to all capillaries that use this reagent name in the run file.
 - **Primary antibody names** - All primary antibody names entered in the assay template (Assay screen) will display in the drop down list, otherwise the default name of Primary will be shown. Select a name to apply group settings to all capillaries that use the primary antibody name in the run file.
 - **Attributes** - All primary antibody attributes entered in the assay template (Assay screen) will display in the drop down list. Select an attribute to apply group settings to all capillaries that use this attribute in the run file.
 - **Cycle 1-8** - When this option is selected, group settings will be applied to all capillaries in the cycle.
 - **Custom settings** - Lets you choose specific capillaries to apply the group settings to. When this option is selected, the following pop-up box appears to let you enter a specific capillary number or range of capillaries:




6. If you need to change the peak names group used for a data set, click the cell in the **Settings** column and click the down arrow. Select a group from the drop down list.



7. Repeat the previous steps to apply other groups to specific run data.
8. To remove a data set, click on its cell in the **Apply To** column, then click **Remove**.


- Click **OK** to save changes. Named peaks will be identified with a peak name label in the electropherogram and virtual blot, and will be color coded in the peaks and capillaries tables:



Standard Curve Settings

To use a standard curve to quantitate the concentration of a target protein detected either by an Immunoassay or a Total Protein Assay, first create peak names groups as described earlier for your standard curve protein and the target protein. In the example below, GST-AKT1 is the standard curve protein and AKT1 is the target protein.

NOTE: Analysis settings are run file-specific. However, settings can be imported or exported for use with other run files. For more information see "Importing and Exporting Analysis Settings" on page 284.



To set up a standard curve:

1. Click the **Define...** button.

Standard Curve

Peak: GST-AKT

Fit: Linear


Capillary	Concentration
Jurkat	62.5
Jurkat	31.2

Add Remove

Units: pg/uL

OK Cancel

2. From the Peak drop down list, select the peak name for your standard curve protein.
3. Choose either a **Linear** or **4 Parameter (4PL)** curve fit from the **Fit** drop down list.
4. Select the checkbox in the first cell of the **Capillary** column. Sample names and attributes will be displayed in the list. Any number included in the name will auto-populate the Concentration column when selected.
5. To add another concentration, click **Add** under the Standard Curve table.
6. Repeat the previous steps to enter information for other concentrations. In the following example, six concentrations were entered:



The Standard Curve dialog box contains the following elements:

- Peak: GST-AKT
- Fit: Linear
- Table with columns Capillary and Concentration:

Capillary	Concentration
62.5 pg/uL	62.5
31.2 pg/uL	31.2
15.6 pg/uL	15.6
7.8 pg/uL	7.8
3.9 pg/uL	3.9
0 pg/uL	0.0

Buttons: Add, Remove, OK, Cancel


Units: pg/uL

To remove a concentration, select its row and click **Remove**.

7. Enter the concentration units in the box (for example, pg/ μ L).
8. Click **OK** to save changes.

Applying Peak Names Groups to Run Standard Curve


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Names** in the options list.
2. Click on the group in the analysis settings box you want to apply to specific run data.




The Analysis Settings dialog box contains the following elements:

- Analysis Settings
- List of peak names: AKT1, GST-AKT
- Buttons: Add, Remove

3. Application of peak names groups to specific run data is done in the override box. Click **Add** under the override box. A default override data set will be created.




4. Click the cell in the Apply To column, then click the down arrow.




5. Select an option from the drop down list. This will apply the settings group selected to specific run data as follows:
 - a. **All** - When this option is selected, peak names group settings will be applied to all capillaries in the run file.
 - b. **Blocking reagent** - When this option is selected, group settings will be applied to all capillaries that use this reagent name in the run file.
 - c. **Primary antibody names** - All primary antibody names entered in the assay template (Assay screen) will display in the drop down list, otherwise the default name of Primary will be shown. Select a name to apply group settings to all capillaries that use the primary antibody name in the run file.
 - d. **Attributes** - All primary antibody attributes entered in the assay template (Assay screen) will display in the drop down list. Select an attribute to apply group settings to all capillaries that use this attribute in the run file.
 - e. **Cycle 1-8** - When this option is selected, group settings will be applied to all capillaries in the cycle.


- f. **Custom settings** - Lets you choose specific capillaries to apply the group settings to. When this option is selected, the following pop-up box appears to let you enter a specific capillary number or range of capillaries:



6. If you need to change the peak names group used for a data set, click the cell in the Settings column and click the down arrow. Select a group from the drop down list.




7. Repeat the previous steps to apply other groups to specific run data.
8. To remove a data set, click on its cell in the **Apply To** column, then click **Remove**.
9. Click **OK** to save changes. The curve fit will be visible on the Std Curve tab and the concentration of the proteins will appear in the Peaks table:



System or Loading Control Settings


It is possible to use a system or loading control protein to normalize data between capillaries, between runs and between instruments. First create peak names groups as described earlier for your control protein and the target protein. In the example below, System Control is the control protein and ERK1 is the target protein.

NOTE: Analysis settings are run file-specific. However, settings can be imported or exported for use with other run files. For more information see "Importing and Exporting Analysis Settings" on page 284.



To identify the control protein:


1. Click the **Control** box.
2. Select a Reference Capillary from the pull-down menu. The peak area of this capillary will be used to normalize the peak area of the named peaks.
3. To manually assign the Control Area, click the box and enter the peak area of your control protein.
4. Click **OK** to save changes. Compass will automatically normalize the peak area of the target protein against the control protein. The corrected area appears in the Capillaries table:




Standards Settings

The standards analysis settings page lets you view and change the molecular weight and position for ladder and fluorescent standards, set the registration peak, and change the capillary used for the ladder. To access these settings, select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list. The default standard settings for both the 12-230 kDa and 66-440 kDa size ranges are shown in the following examples:

NOTE: Settings can be modified in an assay prior to starting a run, or in a run file once the run has finished executing. Analysis changes made to an executing run will not be saved to the final run file.





- Click **Import** to import an analysis settings file. This will be explained in more detail in “Importing Analysis Settings” on page 284.
- Click **Export** to export the current analysis settings file. This will be explained in more detail in “Exporting Analysis Settings” on page 284.
- Click **Restore Original** to restore Compass default settings.
- Click **OK** to save changes and exit.
- Click **Cancel** to exit without saving changes.

Standards Analysis Settings Groups


Standards settings are saved as a group, and multiple settings groups can be created. Specific group settings can then be applied to individual capillaries, sample names or attributes in the run data.

NOTES:

We recommend using the Compass default values for standards analysis settings. These settings are included in the default Standards group.

Analysis settings are run-file specific. However, settings can be imported or exported for use with other run files. For more information see "Importing and Exporting Analysis Settings" on page 284.

Standards groups are displayed in the analysis settings box:



The Standards group shown contains the Compass default settings. You can make changes to this group and create new groups.

To view settings for a group, click on the group name in the analysis settings box.

Changing the Capillary Used for the Ladder

Known ladder standards are used to calculate the molecular weights of unknown sample proteins. As noted in each of the Master Kit Product Inserts, we strongly recommend that you use capillary 1 for the ladder. However, you can change the ladder capillary as needed, or opt to not use a ladder at all.

NOTES:


When the ladder capillary is set to none, fluorescent standards information is used to calculate sample protein molecular weight instead of the ladder.

The reported molecular weight for immunodetected and total protein sample proteins in Compass may vary slightly from predicted molecular weights based on sample and assay conditions.

To change the ladder capillary:

1. Select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list.

- Click the arrow in the drop down list next to Ladder Capillary, then click a capillary number or none from the list.




Compass will use the data in the selected capillary to recalculate molecular weights for sample proteins in the run data using the information in the ladder table. If none is selected, Compass will instead use the information in the fluorescent peaks table (fluorescent standards) to calculate molecular weight for sample proteins.


NOTE: When the ladder capillary is set to none, the ladder table becomes inactive and cannot be modified.

Creating a New Standards Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list.
2. Click **Add** under the analysis settings box. A new group will be created:



3. Click on the new group and enter a new name.



4. Click in the first cell in the **MW** column in the Fluorescent Peaks table.




New Standards

Fluorescent Peaks


MW	Position	Fit	Registration
1	200	<input type="checkbox"/>	<input checked="" type="checkbox"/>
29	400	<input checked="" type="checkbox"/>	<input type="checkbox"/>
230	900	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Add Remove


5. Enter the molecular weight (in kDa) for the fluorescent standard.



6. Click in the first cell in the **Position** column.




7. Enter the position of the fluorescent standard peak.



NOTE: Standards peak positions are relative to each other. Only the difference in their position is used to help identify the standard peaks. When entering standard peak information for the first time, review the standards data in the Analysis screen to find the correct peak position.

8. Repeat the steps above for the remaining standards in the table.

- **To add another standard** - Click **Add** under the peak table, then modify the information in the new row.
 - **To remove a standard** - Select its row and click **Remove**.
9. Select which standard should be used for capillary registration by clicking the checkbox in the **Registration** column. The first standard is typically used for the registration.




The screenshot shows a dialog box titled "New Standards" with a sub-header "Fluorescent Peaks". It contains a table with the following data:

MW	Position	Fit	Registration
2	115	<input type="checkbox"/>	<input checked="" type="checkbox"/>
29	400	<input checked="" type="checkbox"/>	<input type="checkbox"/>
230	900	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Below the table are two buttons: "Add" and "Remove".

NOTE: In order for Compass to perform data analysis, at least one peak must be selected for registration.

10. Select which standards should be used for molecular weight determination of sample proteins by clicking the checkbox in the **Fit** column. The standards not used for registration are typically used for fit.




The screenshot shows a dialog box titled "New Standards" with a sub-header "Fluorescent Peaks". It contains a table with the following data:

MW	Position	Fit	Registration
2	115	<input type="checkbox"/>	<input checked="" type="checkbox"/>
29	400	<input checked="" type="checkbox"/>	<input type="checkbox"/>
230	900	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Below the table are two buttons: "Add" and "Remove".

11. Click the arrow in the drop down list next to Ladder Capillary, then click a capillary number or none from the list. Capillary 1 is typically used for the ladder.




Compass will use the data in the selected capillary to calculate the molecular weights for sample proteins using the information in the ladder table. If none is selected, Compass will instead use the information in the fluorescent peaks table (fluorescent standards) to calculate molecular weight for sample proteins.

NOTES:

When the ladder capillary is set to none, the ladder table becomes inactive and cannot be modified.

The reported molecular weight for immunodetected and total protein sample proteins in Compass may vary slightly from predicted molecular weights based on sample and assay conditions.

12. If a ladder capillary was selected, click in the first cell in the **MW** column in the ladder table.



Enter the molecular weight (in kDa) for the ladder standard.

Ladder Capillary 3

MW	Position	Fit
15	350	<input checked="" type="checkbox"/>
40	520	<input checked="" type="checkbox"/>
66	550	<input checked="" type="checkbox"/>
116	700	<input checked="" type="checkbox"/>
180	850	<input checked="" type="checkbox"/>
230	900	<input checked="" type="checkbox"/>

Add Remove

13. Click in the first cell in the **Position** column.

Ladder Capillary 3

MW	Position	Fit
15	200	<input checked="" type="checkbox"/>
40	520	<input checked="" type="checkbox"/>
66	550	<input checked="" type="checkbox"/>
116	700	<input checked="" type="checkbox"/>
180	850	<input checked="" type="checkbox"/>
230	900	<input checked="" type="checkbox"/>

Add Remove

Enter the position of the ladder standard peak.

Ladder Capillary 3

MW	Position	Fit
15	205	<input checked="" type="checkbox"/>
40	520	<input checked="" type="checkbox"/>
66	550	<input checked="" type="checkbox"/>
116	700	<input checked="" type="checkbox"/>
180	850	<input checked="" type="checkbox"/>
230	900	<input checked="" type="checkbox"/>

Add Remove

NOTE:

Ladder peak positions are relative to each other. Only the difference in their position is used to help identify the ladder peaks. When entering ladder peak information for the first time, review the ladder data in the Analysis screen to find the correct peak position.

14. Repeat the steps above for the remaining ladder standards in the table.

- **To add another ladder standard** - Click **Add** under the table, then modify the information in the new row.
 - **To remove a ladder standard** - Select its row and click **Remove**.
15. Select which ladder standards should be used for molecular weight determination of sample proteins by clicking the checkbox in the **Fit** column. All ladder standards are typically used for fit.

Ladder Capillary 3

MW	Position	Fit
15	205	<input checked="" type="checkbox"/>
45	515	<input checked="" type="checkbox"/>
70	615	<input checked="" type="checkbox"/>
100	715	<input checked="" type="checkbox"/>
120	815	<input checked="" type="checkbox"/>
195	920	<input checked="" type="checkbox"/>

Add Remove

16. To use the new group as the default analysis settings for the run file data, click the arrow in the drop down list next to Default, then click the new group from the list. Analysis settings in the new group will then be applied to the run data.

Analysis Settings

Standard

New Standards

Add Remove

Default

Standard


Standard

New Standards

17. Click **OK** to save changes.

Changing the Default Standards Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list.
2. Click the arrow in the drop down list next to Default, then click a new default group from the list.



3. Click **OK** to save changes. Analysis settings in the group selected will be applied to the run data.

Modifying a Standards Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list.
2. Click on the group in the analysis settings box you want to modify.



3. Modify fluorescent standards and ladder standards information as described in "Creating a New Standards Group" on page 275.
4. Click **OK** to save changes. The new analysis settings will be applied to the run data.

Deleting an Analysis Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list.
2. Click on the group in the analysis settings box you want to delete and click **Remove**.




3. Click **OK** to save changes.

Applying Analysis Groups to Specific Run Data


1. Select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list.
2. Click on the group in the analysis settings box you want to apply to specific run data.




3. Application of analysis groups to specific run data is done in the override box. Click **Add** under the override box. A default override data set will be created from sample information found in the run file.




4. Click the cell in the **Apply To** column, then click the down arrow.



5. Select an option from the drop down list. This will apply the settings group selected to specific run data as follows:
- **Sample names** - All sample names entered in the assay template (Assay screen) will display in the drop down list, otherwise the default name of Sample will be shown. Select a name to apply group settings to all capillaries that use this sample name in the run file.
 - **Attributes** - All sample attributes entered in the assay template (Assay screen) will display in the drop down list. Select an attribute to apply group settings to all capillaries that use this attribute in the run file.
 - **Cycle 1-8** - When this option is selected, group settings will be applied to all capillaries in the cycle.
 - **Custom settings** - Lets you choose specific capillaries to apply the group settings to. When this option is selected, the following pop-up box appears to let you enter a specific capillary number or range of capillaries:



6. If you need to change the analysis group used for a data set, click the cell in the **Settings** column and click the down arrow. Select a group from the drop down list.



7. Repeat the previous steps to apply other groups to specific run data.
8. To remove a data set, click on its cell in the **Apply To** column, then click **Remove**.
9. Click **OK** to save changes.

Importing and Exporting Analysis Settings

The analysis settings in a run file can be exported as a separate file. This allows the same analysis settings to be imported into other assays or run files at a later time, rather than having to re-enter settings manually.

Importing Analysis Settings


NOTE: Importing an analysis settings file populates the settings in all Analysis pages.

1. Open the run file or assay you want to import analysis settings to.
2. Select **Edit** in the main menu and click **Default Analysis** (Assay screen) or **Analysis** (Analysis screen).
3. Click **Import** on any page.
4. Select a settings file (*.settings) and click **OK**. The imported settings will display in all analysis pages.

Exporting Analysis Settings

NOTE: Exporting an analysis settings file exports the settings in all Analysis pages.

1. Open the run file or assay you want to export analysis settings from.
2. Select **Edit** in the main menu and click **Default Analysis** (Assay screen) or **Analysis** (Analysis screen).
3. Click **Export** on any page. The following window displays:



4. The default directory is Compass/Assays. Change the directory if needed.
5. Enter a file name and click **Save**. The settings will be saved as a *.settings file.

Chapter 9:

Charge Assay Data Analysis

Chapter Overview

- Analysis Screen Overview
- Opening Run Files
- How Run Data is Displayed in the Analysis Screen
- Viewing Run Data
- Compass Run Data Notifications and Warnings
- Checking Your Results
- Group Statistics
- Copying Data Views and Results Tables
- Exporting Run Files
- Changing Sample Protein Identification
- Changing the Virtual Blot View
- Changing the Electropherogram View
- Closing Run Files
- Compass Analysis Settings Overview
- Advanced Analysis Settings
- Images Analysis Settings
- Peak Fit Analysis Settings
- Peak Names Settings
- Standards Settings
- Importing and Exporting Analysis Settings

Analysis Screen Overview


The Analysis screen is used to view run data including electropherograms, capillary images, lane view data and tabulated results. Any post-run analysis is also done in this screen. To access this screen, click **Analysis** in the screen tab:



Analysis Screen Panes

The Analysis screen has six panes, each displays the following data for up to 96 capillaries per experimental run:

- **Experiment** - Lists the assay protocol steps and assay template information.
- **Graph** - Displays electropherogram data for immunodetected sample proteins, fluorescent standards or capillary registrations.
- **Image** - Displays a 12-capillary image of the separated immunodetected sample proteins, fluorescent standards or capillary registrations.
- **Lane** - Displays data for immunodetected sample proteins, fluorescent standards or capillary registrations as bands in individual lanes. This virtual blot-like image is similar to traditional blot results.
- **Peaks** - Lists the tabulated results for immunodetected sample proteins and information for identified fluorescent standards and capillary registrations.
- **Capillaries** - Displays a list of the immunodetected sample proteins Compass named automatically using the user-defined peak name analysis parameters.



NOTE: The reported pI for immunodetected sample proteins in Compass may vary slightly from predicted pIs based on sample, buffer and assay conditions.

Software Menus Active in the Analysis Screen


The following software menus are available:

- File
- Edit
- View
- Instrument (when Compass is connected to Peggy Sue)
- Window
- Help

The File, Edit and View menu options specific to the Analysis screen are described next.

File Menu


The following File menu options are active:



- **Open Run** - Opens a run file.
- **Add Run** - Opens and views other run files in addition to those that are already open.
- **Close** - Closes the run file currently being viewed.
- **Close All** - Closes all open run files.
- **Save** - Saves changes to the open run file.
- **Save As** - Saves changes to the open run file under a different file name.
- **Export Tables** - Exports the results for all capillaries in the run in .txt format.
- **Export Spectra** - Exports the raw data traces for each capillary in the run in .txt or .cdf format.
- **Exit** - Closes Compass.

Edit Menu


The following Edit menu options are active:



- **Copy** - Lets you copy data shown in the graph, lane, peaks or capillaries panes. See “Copying Data Views and Results Tables” on page 334 for more information.
- **Analysis** - Displays the analysis settings used to analyze the run data and lets you change them as needed. See “Compass Analysis Settings Overview” on page 372 for more information.
- **Preferences** - Lets you set and save custom preferences for data export, plot colors in the graph and Peggy Sue’s Twitter settings. See Chapter 10, “Setting Your Preferences” for more information.

View Menu

The following View menu options are active:



- **Single View** - Displays data in a per capillary (single) view format.
- **Multiple View** - Displays data in a per 12-capillary (multiple) view format.
- **Standards** - Lets you change the data view to show only the fluorescent standards.
- **Registration** - Lets you change the data view to show only the capillary registrations.
- **Samples** - Lets you change the data view to show only immunodetected sample proteins.
- **Grouping** - Lets you analyze replicates by calculating the mean, standard deviation and CV of named proteins.

- **Filter** - Lets you display data only for specific capillaries or named proteins.
- **View Region** - Lets you change the pI (x-axis) range of the data displayed.
- **Show Hidden**- Shows capillaries that are hidden from the data view.


Opening Run Files

You can open one or multiple run files at a time to compare data between runs.

Opening One Run File

NOTE: If you need to open a run file when another run is executing, launch another instance of Compass first, then open the file.

1. Select **File** in the main menu and click **Open Run**.




2. A list of the last 10 runs opened will display. Select one of these runs or click **Browse** to open the Runs folder and select a different file.


Opening Multiple Run Files

NOTE: If you need to open a run file when another run is executing, launch another instance of Compass first, then open the file.


1. To open the first run file, select **File** in the main menu and click **Open Run**.



2. A list of the last 10 runs opened will display. Select one of these runs or click **Browse** to open the Runs folder and select a different file.
3. To open another run file, select **File** in the main menu and click **Add Run**.



4. A list of the last 10 runs opened will display. Select one of these runs or click **Browse** to open the Runs folder and select a different file. When a run is added, its data will append to the open run file, and will display as a second set of up to 96 capillaries in all screen panes. The second run file name will also appear in the Compass title bar:



5. Repeat the last two steps to open additional runs.

How Run Data is Displayed in the Analysis Screen

Data in the run file is compartmentalized for easy review in one of six panes in the Analysis screen.

Experiment Pane: Assay and Capillary Information

The experiment pane displays assay and capillary information for each of the 96 capillaries in the run. A maximized view of the experiment pane is shown below.

Sample	Primary	Cycle	Cap	S	1	2	3
High phospho HeLa in 5-8	ERK1/2	1	1	A1	I1	A...	
Low phospho HeLa in 5-8	ERK1/2	1	2	A2	I2	A...	
High phospho HeLa in 5-8	ERK1/2	1	3	A3	I3	A...	
Low phospho HeLa in 5-8	ERK1/2	1	4	A4	I4	A...	
High phospho HeLa in 5-8	ERK1/2	1	5	A5	I5	A...	
Low phospho HeLa in 5-8	ERK1/2	1	6	A6	I6	A...	
High phospho HeLa in 5-8	ERK1/2	1	7	A7	I7	A...	
Low phospho HeLa in 5-8	ERK1/2	1	8	A8	I8	A...	
High phospho HeLa in 5-8	ERK1/2	1	9	A9	I9	A...	
Low phospho HeLa in 5-8	ERK1/2	1	10	A...	I10	A...	
High phospho HeLa in 5-8	ERK1/2	1	11	A...	I11	A...	
Low phospho HeLa in 5-8	ERK1/2	1	12	A...	I12	A...	

- **To view all columns** - Click the **Experiment** tab, then use the scroll bar or click **Maximize** in the upper right corner.
- **To resize columns** - Click the **Experiment** tab. Roll the mouse over a column border until the sizing arrow appears, then click and drag to resize.

Experiment pane column descriptions for the Peggy Sue Charge default assay are as follows:

- **Sample** - Sample name. If sample names were entered in the assay template (Assay screen), those names will display here. Otherwise, Sample (default name) will display.

NOTE: Data notification icons will display in the sample column if Compass detected a potential analysis issue or data was manually modified by the user. For more information see "Compass Run Data Notifications and Warnings" on page 318.


- **Primary** - Primary antibody name. If primary antibody names were entered in the assay template (Assay screen), those names will display here. Otherwise, Primary (default name) will display.
- **Cycle** - Run cycle number. There are 12 capillaries in one cycle.

NOTE: Peggy Sue can run up to eight cycles (12 capillaries per cycle) in an experiment. Each cycle is designated in the Experiment summary.

- **Cap** - Capillary number.
- **S** - Well on the assay plate used for sample.
- **1** - Well on the assay plate used for primary antibody.
- **2** - Well on the assay plate used for secondary antibody.
- **3** - Well on the assay plate used for tertiary antibody (if used).

Graph Pane: Electropherogram Data

Click the **Graph** tab to view data for immunodetected sample proteins, fluorescent standards or capillary registrations. Data for samples is shown in the following example, and immunodetected proteins are displayed as peaks:



More Graph view options will be described in more detail in "Changing the Electropherogram View" on page 349.


Image Pane: Capillary Separation Image Data

Click the **Image** tab to view final separation images of immunodetected sample proteins, fluorescent standards or capillary registrations for up to 96 capillaries per experimental run. Image data for samples is shown in the following example:



Lane Pane: Virtual Blot-Like Image Data

Click the **Lane** tab to view data for immunodetected sample proteins, fluorescent standards or capillary registrations as bands in individual lanes. This virtual blot-like image is similar to traditional blot results. Data for samples in the lane view is shown in the following example, and immunodetected proteins are displayed as bands.




To view information for a band, roll the mouse over a band until the info box appears.

NOTE: The reported pI for immunodetected sample proteins in Compass may vary slightly from predicted pIs based on sample, buffer and assay conditions.

Lane data displayed in the virtual blot is automatically aligned by Compass. To view raw, unaligned lane data and learn more about virtual blot viewing options, see "Changing the Virtual Blot View" on page 343.

Peaks Pane: Calculated Results

Click the **Peaks** tab to view tabulated results for immunodetected sample proteins, fluorescent standards or capillary registrations. Each row in the table shows the individual results for each peak detected in each capillary. Data for samples in the peaks table is shown in the following example:



Sample	Primary	Cap	Peak	Name	Position	pI	Height	Area	% Area	Width	S/N
High pho...	ERK1/2	C1:1	1	ppErk1	499	5.468	394.8	6337	58.7	0.0464	517.9
High pho...	ERK1/2	C1:1	2	pErk1	583	5.735	153.9	1411	13.1	0.0265	321.1
High pho...	ERK1/2	C1:1	3	ppERK2	607	5.805	711.2	5864	62.4	0.0239	1538.7
High pho...	ERK1/2	C1:1	4	Erk1	639	5.904	401.0	3054	28.3	0.0220	879.5
High pho...	ERK1/2	C1:1	5	pERK2	693	6.129	75.8	701	7.5	0.0469	159.2
High pho...	ERK1/2	C1:1	6	ERK2	783	6.581	219.5	2831	30.1	0.0562	352.9
Low phos...	ERK1/2	C1:2	1	ppErk1	505	5.504	45.4	477	3.5	0.0300	75.9
Low phos...	ERK1/2	C1:2	2	nErk1	582	5.739	74.9	492	3.6	0.0188	137.0

NOTES:

Peaks that Compass names automatically using the user-defined peak name analysis parameters are color-coded.

The reported pI for immunodetected sample proteins in Compass may vary slightly from predicted pIs based on sample, buffer and assay conditions.

- **To view all rows** - Click the **Peaks** tab, then use the scroll bar or click **Maximize** in the upper right corner.
- **To resize columns** - Click the **Peaks** tab. Roll the mouse over a column border until the sizing arrow appears, then click and drag to resize.

Peak table column descriptions are as follows:

- **Sample** - Sample name. If sample names were entered in the assay template (Assay screen), those names will display here. Otherwise, Sample (default name) will display.
- **Primary** - Primary antibody name. If primary antibody names were entered in the assay template (Assay screen), those names will display here. Otherwise, Primary (default name) will display.
- **Cap** - Cycle and capillary number. For example, C1:3 indicates cycle 1, capillary 3.


NOTE: Peggy Sue runs up to eight cycles in an experiment, 12 capillaries at a time. Capillary and cycle numbers are displayed in the Experiment tab.

- **Peak** - Peak number. Peaks are numbered in order of detection.

- **Name** - Peak name. Displays peaks that Compass named automatically using the user-defined peak name analysis parameters. Cells will be blank if Compass was not able to name the peak or if naming parameters were not entered.
- **Position** - Displays the pixel position of a peak in the image.
- **pl** - Displays the calculated pl for the peak.
- **Height** - Displays the calculated peak height.
- **Area** - Displays the calculated peak area.
- **% Area** - Displays the calculated percent area for the named peak compared to all named peaks. This value results from dividing the individual peak area by the sum of all named peak areas for the capillary and multiplying by 100 (shown for named peak sample data only).
- **Width** - Displays the calculated peak width (shown for sample data only).
- **S/N** - Displays the calculated signal to noise ratio for the peak.

Capillaries Pane: User-Specified Peak Names

Click the **Capillaries** tab to view tabulated results for sample proteins associated with specific primary antibodies in the run data. Compass labels these sample peaks automatically using user-defined peak name settings. Each row in the table shows the individual results for the named peaks detected in each capillary. Data for samples in the capillaries table is shown in the following example.



Sample	Primary	Capillary	ppErk1	pErk1	Erk1	ppERK2	pERK2	ERK2
High pho...	ERK1/2	C1:1	6337	1411	3054	5864	701	2831
Low phos...	ERK1/2	C1:2	477	492	12533	440	275	13959
High pho...	ERK1/2	C1:3	7333	1672	3272	6400	833	3457
Low phos...	ERK1/2	C1:4	591	577	13618	589	388	16427
High pho...	ERK1/2	C1:5	8556	1533	3320	6819	921	3667
Low phos...	ERK1/2	C1:6	710	596	14005	561	330	17124

NOTES:

Peaks that Compass names automatically with user-defined peak name settings are color-coded.

Information displayed for fluorescent standards and capillary registration data will be for identified standards or registration peaks.

- **To view all rows** - Click the **Capillaries** tab, then use the scroll bar or click **Maximize** in the upper right corner.

- **To resize columns** - Click the **Capillaries** tab. Roll the mouse over a column border until the sizing arrow appears, then click and drag to resize.

Capillaries table column descriptions are as follows:

- **Sample** - Sample name. If sample names were entered in the assay template (Assay screen), those names will display here. Otherwise, Sample (default name) will display.
- **Primary** - Primary antibody name. If primary antibody names were entered in the assay template (Assay screen), those names will display here. Otherwise, Primary (default name) will display.
- **Capillary** - Cycle and capillary number. For example, C1:3 indicates cycle 1, capillary 3.

NOTE: Peggy Sue runs 12 capillaries at a time in a cycle and is able to run up to eight cycles in an experiment. The Information on cycle and capillary number is displayed in the Experiment tab.

- **Peak Name Columns** - An individual column per peak name will display for every peak identified by name in the run data. Cells for capillaries in these columns will be blank if Compass did not find peaks automatically using the user-defined peak name analysis parameters (or none were entered).
 - **To view % area in the peak name columns** - Select % in the upper right corner of the pane. This displays the calculated percent area for the named peak compared to all named peaks. This value results from dividing the individual peak area by the sum of all named peak areas for the capillary and multiplying by 100.

Sample	Primary	Capillary	ppErk1	pErk1	Erk1	ppERK2	pERK2	ERK2
High pho...	ERK1/2	C1:1	58.7	13.1	28.3	62.4	7.5	30.1
Low phos...	ERK1/2	C1:2	3.5	3.6	92.8	3.0	1.9	95.1
High pho...	ERK1/2	C1:3	59.7	13.6	26.7	59.9	7.8	32.3
Low phos...	ERK1/2	C1:4	4.0	3.9	92.1	3.4	2.2	94.4
High pho...	ERK1/2	C1:5	63.8	11.4	24.8	59.8	8.1	32.1
Low phos...	ERK1/2	C1:6	4.6	3.9	91.5	3.1	1.8	95.0

- **To view peak area in the peak name columns (default)** - Deselect % in the upper right corner of the pane. This displays calculated peak area for the individual peak only.

Viewing Run Data


Each run file contains the following data for up to 96 capillaries:

- **Sample data** - For the immunodetected proteins in the sample.
- **Standards data** - For the fluorescent standards run with each sample.
- **Registration data** - For tracking capillaries as they are moved for various assay steps.

Data can be viewed for one, all, or individually selected capillaries.

Switching Between Sample, Standards and Registration Data Views

You can switch between viewing sample, standards and registration data in a run file using the View bar or View menu:




Data buttons in the View bar:



Show Standards




Show Registrations



Show Samples

- **To view sample data** - Click **Show Samples** in the View bar or select **View** in the main menu and click **Samples**:




- Data in this view is for immunodetected sample proteins only.
- Graph view data displays electropherograms in chemiluminescence units (y-axis) and pI (x-axis).
- Lane view data displays immunodetected sample proteins only.
- Image view data displays immunodetected sample proteins only.
- Results for each immunodetected protein are shown in the peaks and capillaries tables.

NOTE: The reported pI for immunodetected sample proteins in Compass may vary slightly from predicted pIs based on sample, buffer and assay conditions.

For information on checking and identifying sample peaks, see “Step 4 – Checking Samples” on page 325.


- **To view standards data** - Click **Show Standards** in the View bar or select **View** in the main menu and click **Standards**:



- Data in this view is for the fluorescent standards only. These standards are premixed with each sample prior to analysis.
- Graph view data displays electropherograms in fluorescence (y-axis) and position (x-axis).
- Lane view data displays fluorescent standards only.
- Image view data displays fluorescent standards only.
- Standards are highlighted in both the peaks and capillaries tables and marked with an **S**.

For information on checking and identifying standards peaks, see “Step 2 – Checking Fluorescent Sizing Standards” on page 320.

- **To view registration data** - Click **Show Registrations** in the View bar or select **View** in the main menu and click **Registration**:



- Data in this view shows capillary registration only.
- Graph view data displays electropherograms in fluorescence (y-axis) and position (x-axis).
- Lane view data displays capillary registration only.
- Image view data displays capillary registration only.
- Registrations are highlighted in both the peaks and capillaries tables and marked with an R.

Because capillaries must be moved multiple times during the run, their positions are tracked to account for potential changes in imaging position as they are moved to and from the separation block. To do this, one of the standards is used as a capillary registration. An image of the standard used for registration is taken at the end of the separation step prior to chemiluminescent imaging. The position of this peak in the registration image is then compared to the position of the same standard peak in the standards image. The shift value listed in the peaks table is the change in pixel position between the two images. Compass data is then aligned to account for any changes in capillary position.

For information on checking and identifying registration peaks, see “Step 3 – Checking Capillary Registrations” on page 324.

Selecting and Displaying Capillary Data


You can choose to view data from one, multiple, or all capillaries at once.

- **To look at data for one capillary** - Click a row in the experiment pane. Data for just the row selected will display in all data views and tables. The following example shows sample data displayed in the lane and peaks panes when one sample is selected:

The screenshot shows the Compass software interface with the 'Lane' view selected. The 'Peaks' table is displayed below the lane view, showing data for six peaks from the selected sample.

Sample	Primary	Cap	Peak	Name	Position	pI	Height	Area	% Area	Width	S/N	Basel
Low p...	ERK1/2	C1:6	1	ppERK1	510	5.518	68.6	710	4.6	0.0294	135.9	55
Low p...	ERK1/2	C1:6	2	pERK1	583	5.739	85.9	596	3.9	0.0197	187.4	55
Low p...	ERK1/2	C1:6	3	ppERK2	607	5.813	40.7	561	3.1	0.0392	46.1	55
Low p...	ERK1/2	C1:6	4	ERK1	637	5.901	1588.3	14005	91.5	0.0251	3436.7	55
Low p...	ERK1/2	C1:6	5	pERK2	691	6.117	31.7	330	1.8	0.0544	60.3	55
Low p...	ERK1/2	C1:6	6	ERK2	778	6.569	1195.4	17124	95.0	0.0639	1639.0	55

- To look at data for multiple non-sequential capillaries** - Hold the **Ctrl** key and select the rows you want to view in the experiment pane. Data for only the rows selected will display in all data views and results tables. The following example shows sample data displayed in the graph (in multiple view) and capillaries panes when multiple rows are selected:




- **To look at data for multiple sequential capillaries** - Hold the **Shift** key and select the first and last rows you want to view in the experiment pane. Data for only the selected rows will display in all data views and results tables. The following example shows standards data displayed in the image and peaks panes when multiple sequential rows are selected:

The screenshot displays the Compass software interface for a run titled "DemoData_Charge_HeLa_ERK12 - Compass". The interface is divided into several panes:

- Experiment Pane:** A table listing experimental runs. Rows 3 through 6 are selected, corresponding to the data shown in the peaks table below.
- Graph Pane:** A chromatogram showing peaks. Three horizontal yellow bars highlight the peaks corresponding to the selected rows in the experiment pane.
- Peaks Table:** A table showing the results for the selected peaks. The table has columns for Sample, Primary, Cap, Peak, Name, Position, pI, Height, Area, % Area, Width, S/N, and Be.


Sample	Primary	Cap	Peak	Name	Position	pI	Height	Area	% Area	Width	S/N	Be
Low p...	ERK1/2	C1:2	1	ppERK1	505	5.504	45.4	477	3.5	0.0300	75.9	
Low p...	ERK1/2	C1:2	2	pERK1	582	5.739	74.9	492	3.6	0.0188	137.0	
Low p...	ERK1/2	C1:2	3	ppERK2	607	5.815	28.5	440	3.0	0.0442	29.3	
Low p...	ERK1/2	C1:2	4	ERK1	636	5.903	1466.5	12533	92.8	0.0244	2663.7	
Low p...	ERK1/2	C1:2	5	pERK2	690	6.121	24.7	275	1.9	0.0572	40.4	
Low p...	ERK1/2	C1:2	6	ERK2	778	6.575	1061.0	13959	95.1	0.0583	1362.0	

- **To look at data for all capillaries** - Hold the **Shift** key and select the first and last rows in the experiment pane. Data for all rows will display in all data views and results tables. The following example shows sample data displayed in the lane and peaks panes when all rows are selected. The pane scroll bars can be used to view all 8 cycles:



Switching Between Single and Multiple Views of the Capillaries

You can switch between displaying run data in the graph, image and lane panes in a per capillary (single) format or up to 96 capillary format. This view is selected in the View bar or the View menu.



Capillary view buttons in the View bar:




Single View




Multiple View

- **To view data in a per capillary format** - Click **Single View** in the View bar or select **View** in the main menu and click **Single View**:



Data for the row(s) selected in the experiment pane will display as follows:


- Electropherograms in the graph pane display overlaid or stacked (depending on the graph option chosen).
- Results for only the selected row(s) will display in the peaks and capillaries tables.
- The selected row(s) are highlighted in the image pane:



- Lanes for only the selected row(s) are displayed in the lane pane:




- **To view data up to 96-capillary format** - Click **Multiple View** in the View bar or select **View** in the main menu and click **Multiple View**:




Data for the row(s) selected in the experiment pane will display as follows:

- Electropherograms corresponding to the selected row(s) are highlighted and displayed in the graph pane.
- Results for selected capillaries display in the peaks and capillaries tables and are highlighted.
- The selected row(s) are highlighted in the image pane:



- All selected lanes display in the lane pane, and lanes corresponding to the selected row(s) are highlighted.




Hiding Capillary Data

Capillary data can be hidden from view if needed. To do this:


1. Click the **Experiment** tab.

- Select the rows you want to hide, then right click on one of the selected rows and click **Hide**:



Data for selected rows will be hidden in all data views and results tables (except for the image pane).

- To view hidden rows** - Select **View** in the main menu and click **Show Hidden**. Hidden rows will become visible again in all panes, and hidden rows will be marked with an X in the experiment pane:




- To unhide rows** - Select the hidden row(s). Right click on one of the selected rows and click **Unhide**.

Setting Run Data Display Filters

The filter lets you auto-hide specific capillaries in the run file data, or only show data for peaks that Compass identified automatically using the user-defined peak name analysis parameters.

NOTE: When more than one run file is open, filter settings will be applied to all files.

- To filter data to show specific capillaries only** - Select **View** in the main menu and click **Filter**. Uncheck the boxes for the capillaries you do not want shown, then click **OK**.



Only data for the checked capillaries will display in data views and results tables, unchecked capillaries will be hidden.

- **To filter data to show named peaks only** - Select **View** in the main menu and click **Filter**. Select **Show named peaks only** then click **OK**. Only data for peaks that Compass identified automatically using the user-defined peak name analysis parameters will display in data views and results tables, all other peak data will be hidden.

Compass Run Data Notifications and Warnings

If Compass detects a potential data issue, a notification or warning will display next to the row in the experiment pane. A list of warnings, notifications and corrective actions are as follows:

- 
Manual correction of sample data notification - Indicates the user modified sample data manually, such as adding or removing a sample peak. Rolling the mouse over the icon displays the type of data modification made.

Sample	Primary	Cycle	
✓ High ...	ERK1/2	1	1
Low p...	ERK1/2	1	...
High ...	ERK1/2		Baseline Manual

- 
Standards warning - Indicates that one of the standards may not be identified properly. This can be resolved by manually identifying the standards. Please refer to “Step 2 – Checking Fluorescent Sizing Standards” on page 320 for details. Rolling the mouse over the icon displays warning details.

Sample	Primary	5	2
Sample	Primary	5	3

Standards Warning: Low Confidence

- 
Manual correction of standards data notification - Indicates the user modified the standards data manually. Rolling the mouse over the icon displays the type of data modification made.
- 
Registrations warning - Indicates that a capillary registration was not identified properly. This can be resolved by manually identifying the registration peak. Please refer to “Step 3 – Checking Capillary Registrations” on page 324 for details. Rolling the mouse over the icon displays warning details.

ERK Hi...	ERK1/2	1	
ERK Hi...	ERK1/2	1	

Registration Warning: Large Registration Shift

-  • **Manual correction of registrations notification** - Indicates the user modified the capillary registration data manually. Rolling the mouse over the icon displays the type of data modification made.
-  • **Peak fit warning** - Indicates that a peak cannot be fit properly. This is typically caused either by a very small peak in the electropherogram or a peak that is very close to the end of the pI range. This can be resolved by removing the peak. Please refer to “Step 4 – Checking Samples” on page 325 or “Step 4 – Checking Samples” on page 325 for details. Rolling the mouse over the icon displays warning details.

	Kit low-pho...	anti-H...	2	4
	Kit low-pho...	anti-E...	2	5
Peak Fit Warning: Too many iterations				


Checking Your Results

Compass detects immunodetected proteins, fluorescent standards and capillary registrations and reports results automatically with no user-interaction. However, we recommend you review their data as outlined in this section as good general practice to ensure accurate results. If a data warning is presented in the experiment pane, the following steps will also help you identify and correct any issues.

Step 1 – Review the Fluorescent Sizing Standards Movie

All capillaries should have three fluorescent sizing standards. To verify the standards separated in the capillary correctly:

1. When the run has completed, click the **Run Summary** screen tab.
2. Click the **Separation** tab and play the movie (this will be the fluorescent standards separation).



- For each capillary, verify that the standards visibly separated. Each fluorescent point in the capillary represents one of the three sizing standards.


Step 2 – Checking Fluorescent Sizing Standards

To verify the standards are identified correctly:

- Click the **Analysis** screen tab.
- Click **Show Standards** in the View bar. Verification that the standards have been correctly identified can be done in either the graph or lane panes.

Graph Pane:

- Click **Single View** in the View bar.
- Click on the first row in the experiment pane, then click the **Graph** tab. Check that the electropherogram has the appropriate number of fluorescent pl standard peaks for the pl Standard Ladder you are using. They will also be identified with a green **S** in the peaks table. The pl standards at the low and high end of the pl range in the electropherogram are at higher concentrations as they are also used for capillary registration. In the following example, the pl standards shown are those for pl Standard Ladder 3 (P/N 040-646).



If standards are not identified correctly, they can be manually corrected as follows:

- **If an incorrect peak is identified as a standard** - Right click the peak in the electropherogram or peaks table and select **Not a Standard**. Compass should correctly reassign the remaining peaks as standards, and update the peaks table.


- **To set an unidentified peak as a standard** - Right click the peak in the electropherogram or peaks table and select **Force Standard**. Compass will assign the peak as a standard, and correctly reassign the remaining standard peaks. A lock icon indicating the standard was set manually will display next to the peak in the peaks table.

*NOTE: To remove standards peak assignments that were made manually, right click on the peak in the electropherogram or peaks table and click **Clear**.*

- c. Repeat the previous steps for the remaining rows in the experiment pane to make sure all standards are identified correctly.

Lane Pane:

- a. Click **Multiple View** in the View bar.
- b. Click on the first row in the experiment pane, then click the **Lane** tab. Standards will be bands and identified with a green outline. Check that each lane has the appropriate number of fluorescent pI standard bands for the pI Standard Ladder you are using. They will also be identified with a green **S** in the peaks table. The pI standard bands at the low and high end of the pI range in each lane will display greater intensity as they are at higher concentrations and are also used as registration standards. In the following example, the pI standards shown are those for pI Standard Ladder 3 (P/N 040-646). To view band labels, roll the mouse over the individual bands.




If standards are not identified correctly, they can be manually corrected as follows:

- **If an incorrect band is identified as a standard** - Right click the band in the lane or peaks table and select **Not a Standard**. Compass should correctly reassign the remaining bands as standards.
 - **To set an unidentified band as a standard** - Right click the band in the lane or peaks table and select **Force Standard**. Compass will assign the band as a standard, and correctly reassign the remaining standard bands.
- c. Repeat the previous steps for the remaining rows in the experiment pane to make sure all standards are identified correctly.

Step 3 – Checking Capillary Registrations

All capillaries should have a single capillary registration peak. To verify the registrations are identified correctly:

1. Click the **Analysis** screen tab.
2. Click **Show Registrations** and **Single View** in the View bar.
3. Click the **Graph** tab.
4. Click on the first row in the experiment pane. The registration peak will be the first and largest peak in the electropherogram. Check that the two registration peaks are identified and labeled Reg 1 and Reg 2 in the electropherogram. They will also be identified with a purple **R** in the peaks table.



5. If the registration peak is not identified correctly, right click the correct registration peak in the electropherogram or peaks table and select **Registration Force**. Compass will assign the new peak as the capillary registration. A lock icon indicating the registration was set manually will display next to the peak in the peaks table.

*NOTE: To remove registration peak assignments that were made manually, right click on the peak in the electropherogram or peaks table and click **Registration Clear**.*

6. Repeat the previous steps for the remaining rows in the experiment pane to make sure all registrations are identified correctly.

Step 4 – Checking Samples

All immunodetected sample proteins in the graph and lane panes will be labeled automatically with the calculated protein pI.


NOTE: The reported pI for immunodetected sample proteins in Compass may vary slightly from predicted pIs based on sample, buffer and assay conditions.

To verify that sample proteins are identified correctly:

1. Click the **Analysis** screen tab.
2. Click **Show Samples** in the View bar. Verification that sample proteins have been correctly identified can be done in either the graph or lane panes, but manual corrections must be done in the graph pane.

Graph Pane:

- a. Click **Single View** in the View bar.
- b. Click on the row in the experiment pane that contains the sample you wish to check, then click the **Graph** tab.



If sample peaks are not identified correctly, they can be manually corrected as follows:


- **If an incorrect peak is identified as a sample peak** - Right click the peak in the electropherogram or peaks table and select **Remove peak**. Compass will no longer identify it as a sample peak in the electropherogram and the peak data will be removed in the results tables.
- **To set an unidentified peak as a sample peak** - Right click the peak in the electropherogram or peaks table and select **Add Peak**. Compass will calculate and display the results for the peak in the results tables and identify the peak in the electropherogram.

*NOTE: To remove sample peak assignments that were made manually and go back to the original view of the data, right click in the electropherogram and click **Clear All**.*

- Repeat the previous steps for the remaining sample rows in the experiment pane to make sure all sample proteins are identified correctly.

Lane Pane:

- a. Click either **Single View** or **Multiple View** in the View bar.
- b. Click on the row in the experiment pane that contains the sample you wish to check, then click the **Lane** tab. To view sample protein band labels, roll the cursor over the individual bands. If sample bands are not identified correctly, they must be corrected in the graph pane as described in the previous section.



Step 5 – Assigning Peak Names (Optional)

Compass can identify and automatically name sample proteins associated with specific primary antibodies in the run data using user-specified peak name settings. For more information on how to do this see “Peak Names Settings” on page 391.

Group Statistics

The Grouping View is used to analyze replicates by calculating the mean, standard deviation, %CV and SEM for named proteins (see "Peak Names Settings" on page 391 for detailed information on entering named proteins). Statistics for each protein are also plotted for easy comparison between samples, antibodies, and proteins.

Using Groups


1. A group is automatically created for capillaries with the same sample and primary antibody name, so to use this feature, this information must first be entered into the Assay template. To do this:
 - a. Click the **Assay** tab and go to the Template pane.
 - b. Enter sample names and primary antibody names as described in "Step 6 - Add Assay Plate Annotations (Optional)" on page 90. Be sure to enter the same sample and/or primary antibody names for the groups of samples you want to calculate statistics for.

In this example there are two sample types, **High phospho HeLa** and **Low phospho HeLa** which were run with two different antibodies, **ERK 1/2 Primary 1** and **ERK 1/2 Primary 2**.

	1	2	3	4	5	6	7	8	9	10	11	12
A	High pho...	Low phos...	High pho...	Low phos...	High pho...	Low phos...	High pho...	Low phos...	High pho...	Low phos...	High pho...	Low phos...
B	High pho...	Low phos...	High pho...	Low phos...	High pho...	Low phos...	High pho...	Low phos...	High pho...	Low phos...	High pho...	Low phos...
C			ERK 1/2 Primary 1					ERK 1/2 Primary 2				
D					Goat anti-Rabbit HRP							
E					Luminol/Peroxide							

Each of the two samples were run with the two antibodies in every cycle, and the ERK1/2 antibody generates 6 named peaks.

2. Select the **Analysis** tab. Group your results and view the associated statistics by selecting **View** and clicking on **Grouping**. In the box, click **Enable Grouping**.



3. Select a grouping option by clicking the box next to the option. These options allow you to group capillaries in multiple ways:



- **Group across runs** - Groups capillaries from multiple runs. This option creates fewer groups and generates statistics across multiple runs.
- **Group across cycles** - Groups capillaries run in different cycles.
- **No option selected** - When only one run is open, groups will only contain capillaries from the same cycle. When more than one run is open, groups are created for each run.

Viewing Statistics

Peak and Capillary Groups

The Peak Groups pane reports statistics for each named protein in every group. Each group shows the statistics for named proteins which includes average area, standard deviation, %CV and SEM. The number in parenthesis after the sample name indicates the number of capillaries in the group.

Sample	Primary	Cap	Name	Area	Std.Dev	CV	SEM
Hi-phospho HeLa (24)	ERK1/2 Primary 1		ERK2	3204	972	5.3	198
Hi-phospho HeLa (24)	ERK1/2 Primary 1		Erk1	2841	516	8.2	105
Hi-phospho HeLa (24)	ERK1/2 Primary 1		pERK2	1308	256	9.5	52
Hi-phospho HeLa (24)	ERK1/2 Primary 1		pErk1	887	766	6.4	156
Hi-phospho HeLa (24)	ERK1/2 Primary 1		ppERK2	5105	1368	6.8	279
Hi-phospho HeLa (24)	ERK1/2 Primary 1		ppErk1	8716	1438	6.5	293
Hi-phospho HeLa (24)	ERK1/2 Primary 2		ERK2	2625	1157	4.1	236
Hi-phospho HeLa (24)	ERK1/2 Primary 2		Erk1	2597	680	6.2	139
Hi-phospho HeLa (24)	ERK1/2 Primary 2		pERK2	1116	372	3.3	76
Hi-phospho HeLa (24)	ERK1/2 Primary 2		pErk1	678	396	8.5	81
Hi-phospho HeLa (24)	ERK1/2 Primary 2		ppERK2	3973	1839	6.3	375
Hi-phospho HeLa (24)	ERK1/2 Primary 2		ppErk1	8109	2309	8.5	471
Low-phospho HeLa (24)	ERK1/2 Primary 1		ERK2	9745	1538	5.8	314
Low-phospho HeLa (24)	ERK1/2 Primary 1		Erk1	7707	1444	8.7	295
Low-phospho HeLa (24)	ERK1/2 Primary 1		pERK2	1879	306	6.3	62
Low-phospho HeLa (24)	ERK1/2 Primary 1		pErk1	248	71	8.5	14
Low-phospho HeLa (24)	ERK1/2 Primary 1		ppERK2	397	92	3.2	19
Low-phospho HeLa (24)	ERK1/2 Primary 1		ppErk1	710	194	7.3	40
Low-phospho HeLa (24)	ERK1/2 Primary 2		ERK2	9120	2117	3.2	432
Low-phospho HeLa (24)	ERK1/2 Primary 2		Erk1	7774	1601	6.6	327
Low-phospho HeLa (24)	ERK1/2 Primary 2		pERK2	1878	419	2.3	86
Low-phospho HeLa (24)	ERK1/2 Primary 2		pErk1	225	83	6.9	17
Low-phospho HeLa (24)	ERK1/2 Primary 2		ppERK2	326	139	2.7	28
Low-phospho HeLa (24)	ERK1/2 Primary 2		ppErk1	693	202	9.2	41

Clicking the arrow next to a group lists the individual capillaries in the group and reported data for each capillary:


Sample	Primary	Cap	Name	Area	Std.Dev	CV	SEM
Hi-phospho HeLa (24)	ERK1/2 Primary 1		ERK2	3204	972	5.3	198
Hi-phospho HeLa (24)	ERK1/2 Primary 1		Erk1	2841	516	8.2	105
Hi-phospho HeLa	ERK1/2 Primary 1	C1:1	Erk1	2792			
Hi-phospho HeLa	ERK1/2 Primary 1	C1:3	Erk1	4132			
Hi-phospho HeLa	ERK1/2 Primary 1	C1:5	Erk1	3735			
Hi-phospho HeLa	ERK1/2 Primary 1	C2:2	Erk1	3115			
Hi-phospho HeLa	ERK1/2 Primary 1	C2:4	Erk1	3704			
Hi-phospho HeLa	ERK1/2 Primary 1	C2:6	Erk1	3397			
Hi-phospho HeLa	ERK1/2 Primary 1	C3:1	Erk1	2565			
Hi-phospho HeLa	ERK1/2 Primary 1	C3:3	Erk1	3428			
Hi-phospho HeLa	ERK1/2 Primary 1	C3:5	Erk1	3034			
Hi-phospho HeLa	ERK1/2 Primary 1	C4:2	Erk1	2663			
Hi-phospho HeLa	ERK1/2 Primary 1	C4:4	Erk1	3015			
Hi-phospho HeLa	ERK1/2 Primary 1	C4:6	Erk1	2821			
Hi-phospho HeLa	ERK1/2 Primary 1	C5:1	Erk1	2230			
Hi-phospho HeLa	ERK1/2 Primary 1	C5:3	Erk1	2695			
Hi-phospho HeLa	ERK1/2 Primary 1	C5:5	Erk1	2569			
Hi-phospho HeLa	ERK1/2 Primary 1	C6:2	Erk1	2536			
Hi-phospho HeLa	ERK1/2 Primary 1	C6:4	Erk1	2762			
Hi-phospho HeLa	ERK1/2 Primary 1	C6:6	Erk1	2713			
Hi-phospho HeLa	ERK1/2 Primary 1	C7:1	Erk1	2352			
Hi-phospho HeLa	ERK1/2 Primary 1	C7:3	Erk1	2278			
Hi-phospho HeLa	ERK1/2 Primary 1	C7:5	Erk1	2419			
Hi-phospho HeLa	ERK1/2 Primary 1	C8:2	Erk1	2122			
Hi-phospho HeLa	ERK1/2 Primary 1	C8:4	Erk1	2639			
Hi-phospho HeLa	ERK1/2 Primary 1	C8:6	Erk1	2474			
Hi-phospho HeLa (24)	ERK1/2 Primary 1		pERK2	1308	256	9.5	52
Hi-phospho HeLa (24)	ERK1/2 Primary 1		pErk1	887	766	6.4	156
Hi-phospho HeLa (24)	ERK1/2 Primary 1		ppERK2	5105	1368	6.8	279
Hi-phospho HeLa (24)	ERK1/2 Primary 1		ppErk1	8716	1438	6.5	293
Hi-phospho HeLa (24)	ERK1/2 Primary 2		ERK2	2625	1157	4.1	236
Hi-phospho HeLa (24)	ERK1/2 Primary 2		Erk1	2597	680	6.2	139
Hi-phospho HeLa (24)	ERK1/2 Primary 2		pERK2	1116	372	3.3	76
Hi-phospho HeLa (24)	ERK1/2 Primary 2		pErk1	678	396	8.5	81

The **Capillary Groups** pane pivots the **Peak Groups** results to show statistics for named protein peaks in individual columns.

Sample	Primary	Capillary	ppErk1	Std.Dev	% CV	SEM	pErk1	Std.Dev	% CV	SEM	Erk1	Std.Dev	% CV	SEM
Hi-phospho HeLa (24)	ERK1/2 Primary 1		8716	1438	6.5	293	887	766	6.4	156	2841	516	8.2	105
Hi-phospho HeLa (24)	ERK1/2 Primary 2		8109	2309	8.5	471	678	396	8.5	81	2597	680	6.2	139
Low-phospho HeLa (24)	ERK1/2 Primary 1		710	194	7.3	40	248	71	8.5	14	7707	1444	8.7	295
Low-phospho HeLa (24)	ERK1/2 Primary 2		693	202	9.2	41	225	83	6.9	17	7774	1601	6.6	327

Group Plots


The mean values for named peaks in each group are plotted in bar graphs with error bars showing the standard deviation. The plots compare different antibodies for the same sample and different samples for the same antibody to allow a choice of presentation.



Hiding or Removing Capillaries in Group Analysis

Hidden capillaries are not included in groups. However, hiding capillaries provides an easy way to reject individual capillaries from the statistical analysis. See "Hiding Capillary Data" on page 315 for details on how to do this.

Larger groups of capillaries can also be excluded from analysis. To do this, select **View** and click **Filter**.




Uncheck the box next to the cycles or capillaries you wish to remove and click **OK**. This data will now be removed from the grouped view statistics.

Copying Data Views and Results Tables

You can copy and paste data and results tables into other documents, or save a data view as a graphic file.

Copying Data Views

1. Click in the graph or lane pane.
2. Select **Edit** in the main menu and click **Copy**, or right click and select **Copy**.
3. If you selected copy from the graph pane the following window will display, click **Copy**. If you selected copy from the lane pane skip to the next step.




4. Open a document (Microsoft® Word®, Excel®, PowerPoint®, etc.). Right click in the document and select **Paste**. A graphic of the copied data view will be pasted into the document.

Copying Results Tables

1. Click in the peaks or capillaries pane.
2. Select one or multiple rows.
3. Select **Edit** in the main menu and click **Copy**, or right click and select **Copy**.
4. Open a document (Microsoft® Word®, Excel®, PowerPoint®, etc.). Right click in the document and select **Paste**. Data for the rows selected will be pasted into the document.

Saving the Graph View as an Image File

1. Click in the graph pane.
2. Select **Edit** in the main menu and click **Copy**, or right click and select **Copy**.
3. Select an image option (EMF, PNG or PDF) in the pop-up window, then click **Save**.



4. Select a directory to save the file to and enter a file name, then click **OK**.

Exporting Run Files

Results tables and raw plot data can be exported for use in other applications.

Exporting Results Tables

To export the information in the peaks and capillaries tables:


1. Click **File** in the main menu and click **Export Tables**.
2. Select a directory to save the files to and click **OK**. Data will be exported in .txt format.

NOTE: To exclude export of standards data or export results table data in .csv format, see "Setting Data Export Options" on page 413.

Exporting Raw Sample Electropherogram Data

To export raw sample plot data:

1. Click **File** in the main menu and click **Export Spectra**.




- **To export data in .txt format** - Select **Text Format**. Plots will be exported in one file for all capillaries.
 - **To export data in .cdf format** - Select **Andi Format**. Plots will be exported in one file per capillary.
2. Select a directory to save the files to and click **OK**. Data will be exported in the selected format.

Changing Sample Protein Identification

Compass allows you to customize what sample proteins are reported in the results tables by making manual adjustments in the electropherogram or peaks table.

Adding or Removing Sample Data

1. Click **Show Samples** in the View bar.
2. Click **Single View** in the View bar.
3. Click on the row in the experiment pane that contains the sample you wish to correct, then click the **Graph** tab.
 - **To remove a peak from the data** - Right click the peak in the electropherogram or peaks table and select **Remove peak**. Compass will no longer identify it as a sample peak in the electropherogram and the peak data will be removed in the results tables.



- **To add an unidentified peak to the data** - Right click the peak in the electropherogram or peaks table and select **Add Peak**. Compass will calculate and display the results for the peak in the results tables and identify the peak in the electropherogram.

NOTES:


To remove sample peak assignments that were made manually and go back to the original view of the data, right click in the electropherogram and click **Clear All**.

Virtual blot data in the lane pane will also update to reflect changes made in the graph pane.


Hiding Sample Data

You can hide the results for a sample protein in the results tables without completely removing it from the reported results. To do this:

1. Click **Show Samples** in the View bar.
2. Click **Single View** in the View bar.
3. Click on the row in the experiment pane that contains the sample you wish to correct, then click the **Graph** tab.
4. Right click the peak in the electropherogram or peaks table and select **Hide peak**. Compass will hide the peak data in the results tables.



- To view hidden peak data, click **View** in the main menu and click **Show Hidden**. Hidden peak data will display in the results table and be marked with an **X**.




- To unhide a peak, right click on the peak in the electropherogram or peaks table and select **Unhide Peak**.

Changing Peak Names for Sample Data

If Compass did not automatically name a sample protein associated with a specific primary antibody, this can be adjusted manually. To do this:

- Click **Show Samples** in the View bar.
- Click **Single View** in the View bar.
- Click on the row in the experiment pane that contains the sample you wish to correct, then click the **Graph** tab.

- Right click the peak in the electropherogram or peaks table and click **Name**, then click a name in the list. Compass will change the peak name in the electropherogram and results tables, and adjust peak names for other sample proteins accordingly.



NOTES:


For details on how to specify peak name settings, see "Peak Names Settings" on page 391.

Virtual blot data in the lane pane will also update to reflect changes made in the graph pane.


Displaying Sample Data for Named Peaks Only

You can adjust the sample data to display results only for user-specified named peaks. To do this:

1. Click **Show Samples** in the View bar.
2. Click **View** in the main menu and click **Filter**.
3. Check the **Show Named Peaks only** box and click **OK**.



Compass will hide peak data in the electropherogram and results tables for all unnamed peaks, and instead only display data for named peaks in the electropherogram and results tables.



Changing the Virtual Blot View

Options in the lane pane let you change the contrast or invert the virtual blot, remove baseline noise, change lane labels or overlay standards data on sample lanes.

The lane pane toolbar has the following options:

-  Contrast Adjustment
-  Invert
-  Edit Labels
-  Remove Baseline
-  Overlay Standards Data

Adjusting the Contrast


1. Click the **Contrast Adjustment** button. The contrast tool will display:



2. Click the bar and drag it up or down to adjust the contrast.
3. When finished, click **X** to close the tool.

Inverting the Virtual Blot

1. Click the **Invert** button. The virtual blot image will invert:




2. Click the **Invert** button again to return to the default view.

Selecting Lane Labels

The labels shown above the lanes in the virtual blot can be customized. To do this:

1. Click the **Edit Labels** button. The label box will display:




2. Check one or multiple label boxes, and uncheck those you don't want to display. To remove labels completely, uncheck all boxes. The following label options are available:
 - **Sample** - Sample name. If sample names were entered in the assay template (Assay screen), those names will display here. Otherwise, Sample (default name) will display.
 - **Primary** - Primary antibody name. If primary antibody names were entered in the assay template (Assay screen), those names will display here. Otherwise, Primary (default name) will display.
 - **Secondary Ab** - Secondary antibody name. If a secondary antibody name was entered in the assay template (Assay screen), those names will display here. Otherwise, Secondary (default name) will display.
 - **Capillary** - Cycle and capillary number. For example, C1:3 indicates cycle 1, capillary 3

NOTE: Peggy Sue runs up to eight cycles (12 capillaries at a time), so the cycle number displayed will always be C1-C8 depending on the number of cycles programmed.

- **Attributes** - Attribute text. If attribute information was entered for any of the rows or wells in the assay template (Assay screen), up to four can be selected and displayed as lane labels.

Viewing the Uncorrected Sample Baseline

1. Click the **Remove Baseline** button (active for sample data only). This will remove the automatic baseline correction.



2. Click the **Remove Baseline** button again to return to the default, baseline corrected view.

Overlaying Standards Data on Sample Lanes

Data for the standards can be overlaid on the sample data in the virtual blot so you can view the raw, unaligned (uncorrected) lane data. To do this:

1. Click the **Overlay Standards Data** button (active for sample data only). An overlay of the raw sample and standards data will display:




2. Click the **Overlay Standards Data** button again to return to the default, auto-aligned view.

Moving Lanes in the Virtual Blot View

Lanes in the virtual blot view can be moved to make it easier to compare specific lanes of data. To do this:

1. Click a lane in the virtual blot. Hold the mouse button down and drag the lane to a new position.



2. Release the mouse button. The lane will now be repositioned in the virtual blot view.

Changing the Electropherogram View

Options in the graph pane let you zoom and scale electropherograms, overlay or stack plots and change the peak and plot information displayed.

The graph pane toolbar has the following options:



Auto Scale



Graph Options




Stack the Plots




Overlay the Plots

Autoscaling the Electropherogram

Click the **Autoscale** button to scale the y-axis to the largest peak in the electropherogram.




Click the **Autoscale** button again to return to default scaling.



Stacking Multiple Electropherograms

Electropherogram data for multiple capillaries can be stacked vertically in the graph pane for comparison. To do this:

1. Click **Single View**.
2. Select multiple rows in the experiment pane.
3. Click the **Stack the Plots** button. The individual electropherograms for each row selected will stack in the graph pane.




You can customize the colors used for the stacked plot display. For more information see "Selecting Custom Plot Colors for Graph Overlay" on page 414.

Overlaying Multiple Electropherograms

The electropherogram data for multiple capillaries can also be overlaid on top of each other for comparison in the graph pane. To do this:


1. Click **Single View**.
2. Select multiple rows in the experiment pane.
3. Click the **Overlay the Plots** button. The individual electropherograms for each row selected will overlay in the graph pane.




You can customize the colors used for the overlay plot display. For more information see "Selecting Custom Plot Colors for Graph Overlay" on page 414.

Zooming

To zoom in on a specific area of the electropherogram, hold the mouse button down and draw a box around the area with the mouse:




To return to default scaling, right click in the electropherogram and click **Zoom Out**.




Customizing the Data Display

You can customize electropherogram peak labels, plot labels and display options. To do this, select the **Graph Options** button.




Peak Labels

You can customize the labels used to identify peaks in the electropherogram with these options:




Matching Peak Names - Checking this box will draw vertical lines through each named peak. Using this option with Stack the Plots or Overlay the Plots features is useful for visual comparison of named peaks across multiple capillaries.



- **Peak Names** - Checking this box will display peak name labels on all named peaks in the electropherogram.

NOTE: If more than one peak label option is selected, peak name labels will always be used for named peaks.




- **Peak Values** - Checking this box will display the pI labels on all peaks in the electropherogram.

NOTES:

When viewing standards and registration data, this option is called Peak Positions. Labels displayed are peak positions rather than pI.


If more than one peak label option is selected, peak name labels will always be used for named peaks.

The reported pI for immunodetected sample proteins in Compass may vary slightly from predicted pIs based on sample, buffer and assay conditions.




Baseline and Grid Options

You can view the calculated baseline fit, peak integration and show grid lines with these options.




- **Fitted peaks** - Checking this box will display how the peaks were fit by the software.

NOTE: This option is only available for sample data.




- **Baseline Fit** - Checking this box will display the calculated baseline for the peaks. Baseline points will also display for regions of the electropherogram considered to be at baseline.

NOTE: This option is only available for sample data.



- **Grid Lines** - Checking this box will display grid lines in the graph area.



Plot Labels


You can customize the plot labels displayed on the electropherogram with these options.

Plot Label	
<input checked="" type="checkbox"/> Sample	<input type="checkbox"/> Attribute
<input type="checkbox"/> Primary	<input type="checkbox"/> Attribute
<input checked="" type="checkbox"/> Capillary	<input type="checkbox"/> Exposure

Plot labels are shown on the right side of the graph pane.


- **Sample** - Checking this box will display the sample name. If sample names were entered in the assay template (Assay screen), those names will display here. Otherwise, Sample (default name) will display.
- **Primary** - Checking this box will display the primary antibody name. If primary antibody names were entered in the assay template (Assay screen), those names will display here. Otherwise, Primary (default name) will display.
- **Capillary** - Checking this box will display the cycle and capillary number. For example, C1:3 indicates cycle 1, capillary 3.
- **Attributes** - Checking this box will display attribute text. If attribute information was entered for samples or primary antibodies in the assay template (Assay screen), they can be selected as plot labels.

- **Exposure** - Checking this box will display the exposure time(s) used for the data.
The following example shows an electropherogram with all plot labels selected:



Selecting Data Viewing Options


The graph view menu provides multiple options for changing what type of electropherogram data is displayed. To access the view menu, click the down arrow next to the graph pane toolbar:




A check mark next to the menu option indicates it is currently selected. Multiple options can be selected at a time.

NOTE: Unless otherwise noted, graph view menu options are available for sample data only.

- **Sample** - Clicking this option will display raw, uncorrected sample data.




- **Sample Baseline Corrected** - Clicking this option will display sample data with the baseline subtracted (zeroed). This is the default view.




- **Baseline Fit** - Clicking this option will display the calculated baseline for the raw sample data. In the example that follows, both Baseline fit and Sample are selected.

NOTE: This option is selected automatically when Baseline Fit is selected in graph options.




- **Baseline Points** - Clicking this option will display regions of the electropherogram considered to be at baseline. In the example that follows, both Baseline Points and Sample are selected.


NOTE: This option is selected automatically when Baseline Fit is selected in graph options.



- **Fit** - Clicking this option will display the bounding envelope of the fitted peaks as calculated by the software for the raw sample data. In the example that follows, both Fit and Sample are selected.




- **Fit Baseline Corrected** - Clicking this option will display the fitted peaks as calculated by the software for the sample baseline corrected data. In the example that follows, both Fit Baseline Corrected and Sample are selected, the fit plot is on the bottom.



NOTE: When viewing standards or registration data, this option is called Standards Fluorescence or Registration Fluorescence, respectively.

- **Standards** - Checking this box aligns the pI of the raw standards data to the sample data and overlays both electropherograms in the graph pane.




Adding and Removing Baseline Points

Points in the baseline can be added or removed as needed. To do this:

1. Click the **Graph Options** button in the graph pane toolbar and check **Baseline Fit**. This will display baseline points for the raw sample data.
2. Use the mouse to draw a box around the area you want to correct. This will zoom in on the area.


- Right click a baseline point and click **Add Baseline Point** or **Remove Baseline Point**.



*NOTE: To clear the manual addition or removal of baseline points and go back to the original view of the data, right click in the electropherogram and click **Clear All**.*


Selecting the X-Axis pI Range

The pI range used for the x-axis can be changed. To do this, Select **View** in the main menu and click **View Region**. The following pop-up window will display:



- To change the x-axis pI range displayed for the run data** - Enter new values in the Lower or Upper range in pI and click **OK**. Electropherogram and virtual blot data will update to display only the data in the entered range.

- **To see the full x-axis pI range included in the run data** - Check **Full Range**. Electropherogram and virtual blot data will update to display the full range of available data.



NOTE: You can change the default x-axis range that Compass uses. For more information, see "Peak Fit Analysis Settings" on page 384.


Closing Run Files

If more than one run file is open, you can close just one file or all the open files at the same time.


- **To close one of multiple open run files** - In the experiment pane, click on one of the sample rows in the file. Next click **File** from the main menu and click **Close**.
- **To close all open run files** - Select **File** from the main menu and click **Close All**.

Compass Analysis Settings Overview

Compass has a variety of analysis features and settings that you can modify as needed to enhance run data. To access these settings, select **Edit** in the main menu and click **Analysis**. If more than one run file is open, select the run file you want to view analysis settings for from the list:



The following screen will display:




To move between analysis pages in this window, click on an option in the list on the left. The following analysis settings can be user-customized in Compass:

- **Advanced** - Lets you customize analysis settings for samples, standards and image data.
- **Images** - Lets you view the chemiluminescent exposures taken during the run and view data for different exposures in the Analysis screen.
- **Peak Fit** - Lets you customize peak fit settings for sample data.
- **Peak Names** - Lets you enter custom naming settings for sample proteins associated with specific primary antibodies or attributes and have Compass automatically label the peaks in the run data.
- **Standards** - Lets you customize the pI and positions Compass uses to identify fluorescent standards and registration peaks.

Advanced Analysis Settings

The advanced analysis settings page lets you view and change analysis settings for samples, standards and image data. To access these settings, select **Edit** in the main menu and click **Analysis**, then click **Advanced** in the options list.

NOTE: Settings can be modified in an assay prior to starting a run, or in a run file once the run has finished executing. Analysis changes made to an executing run will not be saved to the final run file.



- Click **Import** to import an analysis settings file. This will be explained in more detail in “Importing Analysis Settings” on page 409.
- Click **Export** to export the current analysis settings file. This will be explained in more detail in “Exporting Analysis Settings” on page 410.
- Click **Restore Original** to restore Compass default settings.
- Click **OK** to save changes and exit.
- Click **Cancel** to exit without saving changes.

Standards Settings

- **Peak Width** - The approximate width (at full width half max) used to filter out fluorescence artifacts which improves recognition of standards. The default value is 15.
- **Allowable Drift** - The distance in pixels that standards are expected to move compared to their entered positions in the standards analysis page. This setting assists in recognition of standards. The default value is 100.

Sample Settings

- **Peak Fit Starting Width Ratio** - Focuses the peak fit towards the peak center and aids the overall peak fitting. The default value is 0.5.

Image Settings

- **Median Filter Threshold Ratio** - Pixel ratio used to filter out camera artifacts. The default value is 0.5.
- **Median Filter Threshold Limit** - Pixel threshold value used to filter out camera artifacts. The default value is 100.

Advanced Analysis Settings Groups


Advanced analysis settings are saved as a group, and multiple settings groups can be created. Specific group settings can then be applied to individual capillaries, sample names or attributes in the run data.

NOTES:

We recommend using the Compass default values for advanced analysis settings. These settings are included in the default Advanced group.

Analysis settings are run-file specific. However, settings can be imported or exported for use with other run files. For more information see "Importing and Exporting Analysis Settings" on page 409.

Analysis groups are displayed in the analysis settings box:




The Advanced group shown contains the Compass default analysis settings. You can make changes to this group and create new groups.


To view settings for a group, click on the group name in the analysis settings box.

Creating a New Analysis Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Advanced** in the options list.
2. Click **Add** under the analysis settings box. A new group will be created:



3. Click on the new group and enter a new name.




4. Modify standards, sample or image parameters as needed.
5. To use the new group as the default analysis settings for the run file data, click the arrow in the drop down list next to Default, then click the new group from the list. Analysis settings in the new group will then be applied to the run data.



6. Click **OK** to save changes.

Changing the Default Analysis Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Advanced** in the options list.
2. Click the arrow in the drop down list next to Default, then click a new default group from the list.



3. Click **OK** to save changes. Analysis settings in the group selected will be applied to the run data.

Modifying an Analysis Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Advanced** in the options list.
2. Click on the group in the analysis settings box you want to modify.



3. Modify standards, sample or image parameters as needed.
4. Click **OK** to save changes. The new analysis settings will be applied to the run data.

Deleting an Analysis Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Advanced** in the options list.
2. Click on the group in the analysis settings box you want to delete and click **Remove**.




3. Click **OK** to save changes.

Applying Analysis Groups to Specific Run Data


1. Select **Edit** in the main menu and click **Analysis**, then click **Advanced** in the options list.
2. Click on the group in the analysis settings box you want to apply to specific run data.



3. Application of analysis groups to specific run data is done in the override box. Click **Add** under the override box. A default override data set will be created from sample information found in the run file.




4. Click the cell in the **Apply To** column, then click the down arrow.




5. Select an option from the drop down list. This will apply the settings group selected to specific run data as follows:
- **Sample names** - All sample names entered in the assay template (Assay screen) will display in the drop down list, otherwise the default name of Sample will be shown. Select a sample name to apply group settings to all capillaries that use this sample name in the run file.
 - **Attributes** - All sample attributes entered in the assay template (Assay screen) will display in the drop down list. Select an attribute to apply group settings to all capillaries that use this attribute in the run file.

- **Cycle 1-8** - When this option is selected, group settings will be applied to all capillaries in the cycle.
- **Custom settings** - Lets you choose specific capillaries to apply the group settings to. When this option is selected, the following pop-up box appears to let you enter a specific capillary number or range of capillaries:




6. If you need to change the analysis group used for a data set, click the cell in the **Settings** column and click the down arrow. Select a group from the drop down list.



7. Repeat the previous steps to apply other groups to specific run data.
8. To remove a data set, click on its cell in the **Apply To** column, then click **Remove**.
9. Click **OK** to save changes.

Images Analysis Settings

The Images analysis settings page lets you see what chemiluminescent exposures were taken during the run, and view data for different exposures in the Analysis screen. To access these settings, select **Edit** in the main menu and click **Analysis**, then click **Images** in the options list.



- Click **Import** to import an analysis settings file. This will be explained in more detail in “Importing Analysis Settings” on page 409.
- Click **Export** to export the current analysis settings file. This will be explained in more detail in “Exporting Analysis Settings” on page 410.
- Click **Restore Original** to restore Compass default settings.
- Click **OK** to save changes and exit.
- Click **Cancel** to exit without saving changes.


Exposure Settings

The exposure used for the sample data displayed in the Analysis screen is shown in the All box:

Cycle	Luminescence
All	Exposure 3 240 seconds
1	Exposure 3 240 seconds
2	Exposure 3 240 seconds
3	Exposure 3 240 seconds
4	Exposure 3 240 seconds
5	Exposure 3 240 seconds
6	Exposure 3 240 seconds
7	Exposure 3 240 seconds
8	Exposure 3 240 seconds

NOTE: Peggy Sue runs up to eight cycles (12 capillaries per cycle). The individual exposure setting selected for a specific cycle is applied for all 12 capillaries in that cycle unless selected for All cycles.

- **Exposure #** - Sample data displayed in the Analysis screen is for this specific exposure only. To see the number of exposures and exposure times used for the run data, click the arrow in the drop down list next to **All**.




NOTE: The number of exposures taken and exposure times shown were specified in the assay protocol (Assay screen) prior to executing the run and cannot be modified.

Changing the Sample Data Exposure Displayed

The exposure used for the sample data displayed in the Analysis screen can be changed. To do this:

1. Select **Edit** in the main menu and click **Analysis**, then click **Images** in the options list.
2. Click the arrow in the drop down list next to **All** and select an exposure setting:




3. Click **OK** to save changes and exit. Sample data for the exposure selected will display in the Analysis screen.

Peak Fit Analysis Settings

The peak fit analysis settings page lets you view and change peak fit settings for sample data. To access these settings, select **Edit** in the main menu and click **Analysis**, then click **Peak Fit** in the options list.

NOTE: Settings can be modified in an assay prior to starting a run, or in a run file once the run has finished executing. Analysis changes made to an executing run will not be saved to the final run file.



- Click **Import** to import an analysis settings file. This will be explained in more detail in “Importing Analysis Settings” on page 409.
- Click **Export** to export the current analysis settings file. This will be explained in more detail in “Exporting Analysis Settings” on page 410.
- Click **Restore Original** to restore Compass default settings.
- Click **OK** to save changes and exit.
- Click **Cancel** to exit without saving changes.

Range Settings

- **Minimum** - The pl value below which peaks will not be identified. This value will also be used as the default lower pl range for the data displayed in the electropherogram and virtual blot. The default value is 5.
- **Maximum**: The pl value above which peaks will not be identified. This value will also be used as the default upper pl range for the data displayed in the electropherogram and virtual blot. The default value is 7.

Baseline Settings

- **Threshold** - The variance, or roughness, in a baseline data segment below which a point will be called part of the baseline. The default value is 1.0.
- **Window** - How long baseline data segments are expected to be in pixels. Shorter segments will allow the baseline to follow plateau sections of the signal. The default value is 15.
- **Stiffness** - The amount the baseline is allowed to vary from a straight line. Settings between 0.1 and 1.0 make the baseline fit closer to a straight line. Settings from 1.0 to 10.0 will make the baseline fit follow the data more closely. The default value is 1.0.

Peak Find Settings

- **Threshold** - The minimum signal to noise ratio required for a peak to be identified. A setting of 1.0 will detect many peaks, a setting of 30.0 will detect fewer peaks. The default value is 20.0.
- **Width** - The approximate peak width (at full width half max) in pixels used to detect peaks. The minimum value for this setting is 3.0. Larger widths help to eliminate the detection of shoulder peaks and noise peaks. The default value is 7.0.

Peak Fit Analysis Settings Groups


Peak fit settings are saved as a group, and multiple settings groups can be created. Specific group settings can then be applied to individual capillaries, sample names or attributes in the run data.

NOTES:

We recommend using the Compass default values for peak fit analysis settings. These settings are included in the default Peak Fit group.

Analysis settings are run-file specific. However, settings can be imported or exported for use with other run files. For more information see "Importing and Exporting Analysis Settings" on page 409.

Peak fit groups are displayed in the analysis settings box:




The Peak Fit group shown contains the Compass default analysis settings. You can make changes to this group and create new groups.


To view settings for a group, click on the group name in the analysis settings box.

Creating a New Peak Fit Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Fit** in the options list.
2. Click **Add** under the analysis settings box. A new group will be created:



3. Click on the new group and enter a new name.




4. Modify range, baseline or peak find parameters as needed.
5. To use the new group as the default peak fit settings for the run file data, click the arrow in the drop down list next to Default, then click the new group from the list. Peak fit settings in the new group will then be applied to the run data.



6. Click **OK** to save changes.

Changing the Default Peak Fit Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Fit** in the options list.
2. Click the arrow in the drop down list next to Default, then click a new default group from the list.



3. Click **OK** to save changes. Peak fit settings in the group selected will be applied to the run data.

Modifying a Peak Fit Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Fit** in the options list.
2. Click on the group in the analysis settings box you want to modify.



3. Modify range, baseline or peak find parameters as needed.
4. Click **OK** to save changes. The new peak fit settings will be applied to the run data.

Deleting a Peak Fit Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Fit** in the options list.
2. Click on the group in the analysis settings box you want to delete and click **Remove**.




3. Click **OK** to save changes.

Applying Peak Fit Groups to Specific Run Data


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Fit** in the options list.
2. Click on the group in the analysis settings box you want to apply to specific run data.




3. Application of peak fit groups to specific run data is done in the override box. Click **Add** under the override box. A default override data set will be created from sample information found in the run file.




4. Click the cell in the **Apply To** column, then click the down arrow.



5. Select an option from the drop down list. This will apply the settings group selected to specific run data as follows:
- **Sample names** - All sample names entered in the assay template (Assay screen) will display in the drop down list, otherwise the default name of Sample will be shown. Select a sample name to apply group settings to all capillaries that use this sample name in the run file.
 - **Attributes** - All sample attributes entered in the assay template (Assay screen) will display in the drop down list. Select an attribute to apply group settings to all capillaries that use this attribute in the run file.
 - **Cycle 1-8** - When this option is selected, group settings will be applied to all capillaries in the cycle.
 - **Custom settings** - Lets you choose specific capillaries to apply the group settings to. When this option is selected, the following pop-up box appears to let you enter a specific capillary number or range of capillaries:



- If you need to change the peak fit group used for a data set, click the cell in the **Settings** column and click the down arrow. Select a group from the drop down list.




- Repeat the previous steps to apply other groups to specific run data.
- To remove a data set, click on its cell in the **Apply To** column, then click **Remove**.
- Click **OK** to save changes.

Peak Names Settings

The peak names analysis settings page lets you enter custom naming settings for sample proteins. Compass can automatically label sample peaks in the run data associated with specific blocking reagent names, primary antibody names or attributes. To access these settings, select **Edit** in the main menu and click **Analysis**, then click **Peak Names** in the options list.

NOTE: Settings can be modified in an assay prior to starting a run, or in a run file once the run has finished executing. Analysis changes made to an executing run will not be saved to the final run file.



- Click **Import** to import an analysis settings file. This will be explained in more detail in “Importing Analysis Settings” on page 409.
- Click **Export** to export the current analysis settings file. This will be explained in more detail in “Exporting Analysis Settings” on page 410.
- Click **Restore Original** to restore Compass default settings.
- Click **OK** to save changes and exit.


- Click **Cancel** to exit without saving changes.

Peak Names Analysis Settings Groups

Peak name settings are saved as a group, and multiple settings groups can be created. Specific group settings can then be applied to individual capillaries, blocking reagent names or primary antibody names and attributes in the run data.

NOTE: Analysis settings are run-file specific. However, settings can be imported or exported for use with other run files. For more information see "Importing and Exporting Analysis Settings" on page 409.

Peak name groups are displayed in the analysis settings box:




The Peak Names group shown is a Compass template group. You can make changes to this group and create new groups.

To view settings for a group, click on the group name in the analysis settings box.

Creating a Peak Names Group

1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Names** in the options list.
2. Click on the **Peak Names 1** template group in the analysis settings box.



3. Enter a new name for the group.

Analysis Settings: ERK1

Name	pI	Color	Range	Show
ppErk1	5.5		0.05	<input checked="" type="checkbox"/>

- Click in the first cell in the **Color** column, then click the button.

Analysis Settings: ERK1

Name	pI	Color	Range	Show
ppErk1	5.5	 (0,1) 	0.05	<input checked="" type="checkbox"/>

The color selection box will display:



The color selected will be used to identify the sample protein peak in the peaks and capillaries table in the Analysis screen.

- Click a color or define a custom color and click **OK**. The color selection will update in the table:

Analysis Settings: ERK1

Name	pI	Color	Range	Show
ppErk1	5.5		0.05	<input checked="" type="checkbox"/>

- Click in the first cell in the **Range (%)** column.

Analysis Settings: ERK1

Name	pI	Color	Range	Show
ppErk1	5.5		0.05	<input checked="" type="checkbox"/>

11. Enter a range window for the pI entered. Compass will automatically name peaks found within this percent of the pI. For example, if the pI entered is 5.5 and a 0.1 pI range is used, all peaks between pI 5.4 and 5.6 will be identified with this peak name.

NOTE: The reported pI for immunodetected sample proteins in Compass may vary slightly from predicted pIs based on sample, buffer and assay conditions.

12. Select the checkbox in the first cell of the **Show** column. This will turn peak naming on for the sample protein.

Analysis Settings: ERK1


Name	pI	Color	Range	Show
ppErk1	5.5		0.1	<input checked="" type="checkbox"/>

To turn peak naming off for a particular sample protein, deselect the checkbox in the Show column.

13. To add another sample protein, click **Add** under the analysis settings peak table:

Analysis Settings: ERK1


Name	pI	Color	Range	Show
ppErk1	5.5		0.1	<input checked="" type="checkbox"/>
Peak2	6		0.05	<input checked="" type="checkbox"/>



4. Enter information in the analysis settings peak table as described in “Creating a Peak Names Group” on page 392.
5. Click **OK** to save changes.

Modifying a Peak Names Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Names** in the options list.
2. Click on the group in the analysis settings box you want to modify.



3. Change the information in the analysis settings peak table as described in “Creating a Peak Names Group” on page 392.
4. Click **OK** to save changes.

Deleting a Peak Names Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Names** in the options list.
2. Click on the group in the analysis settings box you want to delete and click **Remove**.




3. Click **OK** to save changes.

Applying Peak Names Groups to Run Data


1. Select **Edit** in the main menu and click **Analysis**, then click **Peak Names** in the options list.
2. Click on the group in the analysis settings box you want to apply to specific run data.




3. Application of peak names groups to specific run data is done in the override box. Click **Add** under the override box. A default override data set will be created.




4. Click the cell in the **Apply To** column, then click the down arrow.




5. Select an option from the drop down list. This will apply the settings group selected to specific run data as follows:
 - **All** - When this option is selected, peak names group settings will be applied to all capillaries in the run file.
 - **Primary antibody names** - All primary antibody names entered in the assay template (Assay screen) will display in the drop down list, otherwise the default name of Primary will be shown. Select a name to apply group settings to all capillaries that use the primary antibody name in the run file.
 - **Attributes** - All primary antibody attributes entered in the assay template (Assay screen) will display in the drop down list. Select an attribute to apply group settings to all capillaries that use this attribute in the run file.
 - **Cycle 1-8** - When this option is selected, group settings will be applied to all capillaries in the cycle.
 - **Custom settings** - Lets you choose specific capillaries to apply the group settings to. When this option is selected, the following pop-up box appears to let you enter a specific capillary number or range of capillaries:



6. If you need to change the peak names group used for a data set, click the cell in the **Settings** column and click the down arrow. Select a group from the drop down list.




- Repeat the previous steps to apply other groups to specific run data.
- To remove a data set, click on its cell in the **Apply To** column, then click **Remove**.
- Click **OK** to save changes. Named peaks will be identified with a peak name label in the electropherogram and virtual blot, and will be color coded in the peaks and capillaries tables:



Standards Settings

The standards analysis settings page lets you view and change the pI and position for fluorescent standards and set the registration peaks. To access these settings, select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list.

NOTE: Settings can be modified in an assay prior to starting a run, or in a run file once the run has finished executing. Analysis changes made to an executing run will not be saved to the final run file.



- Click **Import** to import an analysis settings file. This will be explained in more detail in “Importing Analysis Settings” on page 409.
- Click **Export** to export the current analysis settings file. This will be explained in more detail in “Exporting Analysis Settings” on page 410.
- Click **Restore Original** to restore Compass default settings.
- Click **OK** to save changes and exit.

- Click **Cancel** to exit without saving changes.

Standards Analysis Settings Groups


Standards settings are saved as a group, and multiple settings groups can be created. Specific group settings can then be applied to individual capillaries, sample names or attributes in the run data.

NOTES:

We recommend using the Compass default values for standards analysis settings. These settings are included in the default Standards group.

Analysis settings are run-file specific. However, settings can be imported or exported for use with other run files. For more information see "Importing and Exporting Analysis Settings" on page 409.

Standards groups are displayed in the analysis settings box:




The Std. Ladder groups shown contains the Compass default analysis settings for pI Standard Ladders used with each of the premixes (separation gradients) for charge assays on Peggy Sue. You can select and use one of these default groups, make changes to groups or create new groups.


To view settings for a group, click on the group name in the analysis settings box.

Creating a New Standards Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list.
2. Click **Add** under the analysis settings box. A new group will be created:




3. Click on the new group and enter a new name.




4. Click in the first cell in the **pI** column in the Fluorescent Peaks table.




5. Enter the pI for the fluorescent standard.



6. Click in the first cell in the **Position** column.



7. Enter the position of the fluorescent standard peak.



NOTE: Standards peak positions are relative to each other. Only the difference in their position is used to help identify the standard peaks. When entering standard peak information for the first time, review the standards data in the Analysis screen to find the correct peak position.

8. Repeat the steps above for the remaining standards in the table.
 - **To add another standard** - Click **Add** under the peak table, then modify the information in the new row.
 - **To remove a standard** - Select its row and click **Remove**.
9. Select which standard should be used for capillary registration by clicking the checkbox in the **Registration** column. The first and last standards are typically used for the registration.

New Standards
Fluorescent Peaks

pI	Position	Registration	Fit
4	210	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
4.9	350	<input type="checkbox"/>	<input checked="" type="checkbox"/>
6	520	<input type="checkbox"/>	<input checked="" type="checkbox"/>
6.4	550	<input type="checkbox"/>	<input checked="" type="checkbox"/>
7.3	650	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>


NOTE: In order for Compass to perform data analysis, at least one peak must be selected for registration.

10. Select which standards should be used for pI determination of sample proteins by clicking the checkbox in the **Fit** column. The standards not used for registration are typically also used for fit.

New Standards
Fluorescent Peaks

pI	Position	Registration	Fit
4	210	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
4.9	350	<input type="checkbox"/>	<input checked="" type="checkbox"/>
6	520	<input type="checkbox"/>	<input checked="" type="checkbox"/>
6.4	550	<input type="checkbox"/>	<input checked="" type="checkbox"/>
7.3	650	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>


11. To use the new group as the default analysis settings for the run file data, click the arrow in the drop down list next to Default, then click the new group from the list. Analysis settings in the new group will then be applied to the run data.



12. Click **OK** to save changes.

Changing the Default Standards Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list.
2. Click the arrow in the drop down list next to Default, then click a new default group from the list.



3. Click **OK** to save changes. Analysis settings in the group selected will be applied to the run data.

Modifying a Standards Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list.
2. Click on the group in the analysis settings box you want to modify.



3. Modify fluorescent standards information as described in “Creating a New Standards Group” on page 403.
4. Click **OK** to save changes. The new analysis settings will be applied to the run data.

Deleting an Analysis Group


1. Select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list.
2. Click on the group in the analysis settings box you want to delete and click **Remove**.




3. Click **OK** to save changes.

Applying Analysis Groups to Specific Run Data


1. Select **Edit** in the main menu and click **Analysis**, then click **Standards** in the options list.
2. Click on the group in the analysis settings box you want to apply to specific run data.




- Application of analysis groups to specific run data is done in the override box. Click **Add** under the override box. A default override data set will be created from sample information found in the run file.




- Click the cell in the **Apply To** column, then click the down arrow.



- Select an option from the drop down list. This will apply the settings group selected to specific run data as follows:
 - Sample names** - All sample names entered in the assay template (Assay screen) will display in the drop down list, otherwise the default name of Sample will be shown. Select a name to apply group settings to all capillaries that use this sample name in the run file.
 - Attributes** - All sample attributes entered in the assay template (Assay screen) will display in the drop down list. Select an attribute to apply group settings to all capillaries that use this attribute in the run file.
 - Cycle 1-8** - When this option is selected, group settings will be applied to all capillaries in the cycle.
 - Custom settings** - Lets you choose specific capillaries to apply the group settings to. When this option is selected, the following pop-up box appears to let you enter a specific capillary number or range of capillaries:



- If you need to change the analysis group used for a data set, click the cell in the **Settings** column and click the down arrow. Select a group from the drop down list.



- Repeat the previous steps to apply other groups to specific run data.
- To remove a data set, click on its cell in the **Apply To** column, then click **Remove**.
- Click **OK** to save changes.

Importing and Exporting Analysis Settings

The analysis settings in a run file can be exported as a separate file. This allows the same analysis settings to be imported into other assays or run files at a later time, rather than having to re-enter settings manually.

Importing Analysis Settings

NOTE: Importing an analysis settings file populates the settings in all Analysis pages.


- Open the run file or assay you want to import analysis settings to.

2. Select **Edit** in the main menu and click **Default Analysis** (Assay screen) or **Analysis** (Analysis screen).
3. Click **Import** on any page.
4. Select a settings file (*.settings) and click **OK**. The imported settings will display in all analysis pages.

Exporting Analysis Settings

NOTE: Exporting an analysis settings file exports the settings in all Analysis pages.

1. Open the run file or assay you want to export analysis settings from.
2. Select **Edit** in the main menu and click **Default Analysis** (Assay screen) or **Analysis** (Analysis screen).
3. Click **Export** on any page. The following window displays:



4. The default directory is Compass/Assays. Change the directory if needed.
5. Enter a file name and click **Save**. The settings will be saved as a *.settings file.

Chapter 10:


Setting Your Preferences

Chapter Overview

- Custom Preference Options
- Setting Data Export Options
- Selecting Custom Plot Colors for Graph Overlay
- Setting Up Wes, Sally Sue and Peggy Sue to Send Tweets

Custom Preference Options

You can set and save custom preferences for data export, plot colors in the graph and Twitter communication. To access these settings, select **Edit** in the main menu and click **Preferences**.




To move between preferences pages in this window, click on any option in the list on the left. The following items can be user-customized in Compass:

- **Analysis** - Lets you customize data export options.
- **Graph** - Lets you customize graph color displays.
- **Twitter** - Lets you configure Wes, Sally Sue or Peggy Sue to Tweet run status.

Setting Data Export Options


Select **Edit** in the main menu and click **Preferences**, then click **Analysis** in the options list.



- **Export Standards** - Selecting this option includes data for the standards in each sample when run data is exported. When this option is deselected, only sample data will be exported. This option is selected by default.
- **Export using a comma as the column delimiter** - Selecting this option exports run data in .csv format. When this option is deselected, the data is exported in .txt format.
- Click **Apply** to apply changes to any open run files in Compass.
- Click **Restore Defaults** to restore Compass default settings.
- Click **OK** to save changes and exit.
- Click **Cancel** to exit without saving changes.

Selecting Custom Plot Colors for Graph Overlay

Select **Edit** in the main menu and click **Preferences**, then click **Graph** in the options list.



- **Apply colors to stacked plots** - Selecting this option applies the color scheme to individual electropherograms when the Stack the Plots option is selected in the Analysis screen graph pane.

NOTE: If Apply colors to stack plots is not checked, the colors shown in the preferences screen will be applied only to overlaid electropherograms when the Overlay the Plot option is selected in the graph pane.

- When this option is deselected, plots will use Compass default colors.
- Click **Apply** to apply changes to any open run files in Compass.
- Click **Restore Defaults** to restore Compass default settings.
- Click **OK** to save changes and exit.
- Click **Cancel** to exit without saving changes.

Changing Plot Colors


1. Select **Edit** in the main menu and click **Preferences**, then click **Graph** in the preferences list.
2. Click the color button next to a plot number. The color selection box displays:



3. Select a color or define a custom color and click **OK**. The color button will update to the new color selected.
4. Repeat the steps above for any other plot colors.
5. Check **Apply Colors to Stacked Plots** if you want the new color scheme to also be used for the Stack the Plots option in the graph pane.
6. Click **Apply** to apply changes to plots currently displayed in the graph pane.
7. Click **OK** to save changes and exit. When the Overlay the Plots option is selected in the graph pane, the new color scheme will be used.

Setting Up Wes, Sally Sue and Peggy Sue to Send Tweets

Select **Edit** in the main menu and click **Preferences**, then click **Twitter** in the options list.



- Click **Apply** to apply changes.
- Click **Restore Defaults** to restore Compass default settings.
- Click **OK** to save changes and exit.
- Click **Cancel** to exit without saving changes.

To have Wes, Sally Sue or Peggy Sue tweet a Twitter account:


NOTES:

To set Wes, Sally Sue or Peggy Sue up to tweet, the computer you are using must have an internet connection.


To tweet, Wes, Sally Sue or Peggy Sue must be connected to the internet through a network connection or via the local lab computer.

We recommend setting up a separate Twitter account for Wes, Sally Sue or Peggy Sue. This lets multiple people in the lab follow run progress. It also lets you send tweets directly from Wes, Sally Sue or Peggy Sue to all users, for example to notify others when the instrument is available or when an error has been reset, etc.


1. Click **Set Account**. A set account window will display in Compass and the following browser window will open:




2. Enter a user name or email and password, then click **Authorize app**. A new page will display in the browser with a PIN number.
3. Enter the PIN number in the set account window in Compass and click **OK**:




4. The user name will now appear in the Twitter User Name box. Select one or all of the tweet options in the Tweet When box, then click **Apply**.




5. To confirm the Twitter account is receiving messages, click **Tweet Message**. Enter a test message and click **OK**.



6. If the test Tweet was successful, Compass will display the following message:



7. Click **OK** to save changes and exit. Wes, Sally Sue or Peggy Sue will automatically tweet as the selected options occur, as shown below:



Changing the Twitter Account

To change the Twitter account Wes, Sally Sue or Peggy Sue uses:

1. Select **Edit** in the main menu and click **Preferences**, and click **Twitter** in the preferences list.
2. Click **Clear**.
3. Set up the new account as described in “Setting Up Wes, Sally Sue and Peggy Sue to Send Tweets” on page 416.

Sending Manual Tweets from Wes, Sally Sue and Peggy Sue

You can send tweets directly from Wes, Sally Sue or Peggy Sue. For example, you may want to notify other users that the instrument is available, being serviced or when an error has been cleared. To do this:

1. Select **Edit** in the main menu and click **Preferences** and click **Twitter** in the preferences list.
2. Click **Tweet Message**.
3. Enter a test message and click **OK**. The tweet will be received by any users following the Twitter account Wes, Sally Sue or Peggy Sue uses.

Chapter 11:

Compass Access Control and 21 CFR Part 11 Compliance

Chapter Overview

- Overview
- Enabling Access Control
- Logging In to Compass
- Saving Changes
- Signing Files
- Instrument Command Log
- Run File History
- Troubleshooting Problems and Suggested Solutions
- Authorization Server



Overview


The Compass Access Control feature can be used to help satisfy the 21CFR Part 11 data security requirements when using Simple Western instruments. When Access Control is enabled and the Authorization Server has been installed (see "Authorization Server" on page 432):

- Users are required to log in to Compass when the software is launched
- A history of all actions is maintained
- Data files are signed and encrypted to prevent unauthorized changes (*e.g.*, all files are controlled)
- Each instrument maintains a history of user commands
- Each assay and data file includes a history of signed changes to the file

Compass can be run with or without Access Control enabled. When Access Control is disabled, no user log in is required and files are not encrypted or signed. The instrument history and file history are still maintained but the entries are not signed.

Enabling Access Control

Access Control is enabled in **Preferences**. Select **Edit** in the main menu, click **Preferences**, then select **Access Control**.



To enable Access Control:

1. Check the **Enable** box.
2. Enter the IP address of the Authorization server. Use format X.X.X.X or LocalHost if installing the server on the local machine.


NOTE: Always use the default port setting of 8000, this should not be changed.

3. Close Compass. The next time Compass is launched, a user log in will be required.

*NOTE: Access Control can only be disabled by logging into Compass and deselecting the **Enable** box in the Access Control page of Preferences.*


Logging In to Compass

With Access Control enabled, all users must log in to Compass whenever the software is launched.




Enter your user name and password previously setup by your Compass Administrator.

A successful log in will display the Compass main window with the user information in the lower status bar. The full user name is displayed with the unique user ID in parenthesis:



If there is no activity in Compass for 20 minutes, the user must re-enter their password to perform any controlled actions:




Resolving Log In Issues

Log in failures may occur when:

- The server is temporarily unavailable
- Compass is using the wrong IP address


When this happens, the following message displays:




Click **Disable** to restart Compass with Access Control disabled. Verify or correct the server IP address then close and restart Compass to log in with Access Control enabled.

Saving Changes

When **Save** is selected from the **File** menu, a dialog box will display to allow you to enter a comment before saving the signed file:




The comment is added to the signature entry in the file History:




Signing Files

Select **e-Signature** from the **File** menu to add an electronic signature to a file.



The signed entry will be added to the file History with the meaning of the signature entered in the comment, such as *Approved* or *Verified*.




Date	User Name	Message	Comment
04/15/2013 9:48 AM	b	RunStart2013-04-15_09-49-28_Simon_short	
05/10/2013 4:22 PM	steveg	Saved analysis changes	Changed Baseline fr...
05/10/2013 4:44 PM	steveg	e-Signature applied	Approved


Time: 05/10/2013 4:44 PM User: steveg
Message: e-Signature applied
Comment: Approved

Instrument Command Log

The Instrument Command Log can be viewed at any time by selecting the **Instrument** menu and clicking **Properties**, and then clicking the **Command Log** button:




The Command Log lists all the commands sent to the instrument that were signed by the user who sent the command. If you want to copy the Command Log at any time, right click in the table and select **Copy**, then paste into another document.



Run File History

Select the **Run Summary** screen tab and then the **History** tab to see the file History. To copy the file History, right click in the table and select **Copy**, then paste into another document.



Troubleshooting Problems and Suggested Solutions

If any of the following error messages are encountered, follow the recommended steps below to resolve the issue.

- **Unknown user name or password.**
 - Check if the Caps Lock is on, user name and password are case sensitive.
 - Ask a Compass administrator to confirm your user name. If your password is unknown then the administrator can reset your password (see "Resetting User Passwords" on page 440 for more information).

- **Server not available.**
 - From the **Edit** menu, click **Preferences** and then **Access Control** to confirm the server address is set to the correct Authorization server address. Compass must be able to reach the server on the network.
 - The server must have inbound access to port 8000 enabled.
- **Controlled file cannot be opened without log in.** To open a controlled Run file, enable Access Control by clicking **Edit**, then **Preferences** and **Access Control**. Select **Enable**, close Compass, then re-launch Compass with a valid log in.
- **Uncontrolled file cannot be opened when logged in.** To open an uncontrolled Run file, disable Access Control by clicking **Edit**, then **Preferences** and **Access Control**. Deselect **Enable**, close Compass then re-launch the software.

NOTE: Uncontrolled Assay files can be opened when Compass Access Control is enabled (controlled mode).

- **Command disabled.** Certain commands are only available when a user with the correct permissions is logged in. To change user permissions, use a web browser to log in to the Authorization server web interface at the address shown on the **Access Control** page in **Preferences**, such as: 10.1.3.231:8000.
- **Compass does not prompt for log in.** Compass will only prompt for a log in on launch when Access Control is enabled in Preferences. Enable Access Control by clicking **Edit**, then **Preferences** and **Access Control**. Select **Enable**, close Compass, then re-launch the software. You should now be prompted for a log in.

Authorization Server

The Authorization Server controls the log in access to Compass. In the simplest configuration, the server is run on the same computer as Compass and only that copy of Compass is controlled. A single server can also be used to control access to multiple copies of Compass running on different computers, so long as they have network access to the server. Multiple copies of the server may be run on the same network, and each server will have its own user database.

To enable Compass to use a particular Authorization Server, click **Edit**, then **Preferences** and **Access Control** and enter the server IP address using format X.X.X.X.

NOTES:

Always use the default port setting of 8000, this should not be changed.

If the server is installed on the same computer as Compass (e.g., the local machine), enter LocalHost instead of the IP address. Contact your local IT Administrator to assist with installing the Authorization Server in your preferred format.


Server Administration

The Authorization Server is configured through a web interface at the IP address of the server on port 8000. To access the Server home page, open any browser and type the IP address on port 8000 in a X.X.X.X:8000 or http://X.X.X.X:8000 format. Use LocalHost instead of the IP address if the Server is installed on the local machine.

The default server administrator is:

- User: admin
- Password: admin


After installing the Authorization Server, the administrator user name and password can be changed.




Adding Non-admin Users

Add a user to the server to allow that user to log in to Compass. To do this:


1. Select **Users** from the Site Administration home page:




2. From the Users page, select **Add User**:



3. Fill in the fields to create a new user:




After adding a new user more information can be added:




Permissions

All users can log in to Compass, but the commands available within Compass are controlled by Permission settings. Commands a user does not have permission to use will be disabled. After user permissions have been changed on the server the user must close and re-open Compass to use the new permissions.


Users can belong to groups that have multiple permissions such as Manager or Scientist:




Use the Groups page to change the permissions in a group or create new groups:



To change permissions for a group click **Change**, then select a group:



Move individual group permissions in or out of the Available Permissions and Chosen Permissions boxes by selecting a permission in either box. Click the **left** or **right** arrow button to move the permission into the other box.



Adding Admin Users

To create a user with administrator permissions:

1. Follow the steps described in "Adding Non-admin Users" on page 434 to create the admin user.
2. Under permissions, select **Staff status** and **Superuser status**:

Permissions	
<input checked="" type="checkbox"/> Active	Designates whether this user should be treated as active. Unselect this instead of deleting accounts.
<input checked="" type="checkbox"/> Staff status	Designates whether the user can log into this admin site.
<input checked="" type="checkbox"/> Superuser status	Designates that this user has all permissions without explicitly assigning them.


3. Assign the admin user to a group.

NOTE: Selecting Superuser status enables server permissions only. Admin users must be also be assigned to a group to in order to have Compass permissions.

Resetting User Passwords

To reset a user password:

1. Select **Users** from the Site Administration home page, then select the user to change. The following screen displays:



2. Raw passwords are not stored, they must be changed manually. Click the text link to access the password change form:



3. Enter the new password, then click **Change password**.

Encryption Details

Compass uses the SHA1 hash algorithm to generate a 160 bit hash code that is unique for all files. All files saved by Compass are encrypted with a digital key. This key along with the hash codes guarantees the file history is correct and no other edits were made. All changes saved to a file have the electronic signature of the user who saved the file. The **e-Signature** command allows a user to sign off on a state such as approved or verified.

There is no individual ownership of files, all users who log into Compass can open any file.

