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名称	供应商	货号
10%/mL NALM Cells (included in the Milo Training Kit)	ProteinSimple	035-119
scWest Kit	ProteinSimple	K600
Milk-Free Antibody Diluent	ProteinSimple	043-524
Glycine-HCl	Sigma-Aldrich	G2879-500G
Rabbit anti- β -Actin Antibody	Cell Signaling Technology	8457
Mouse anti- β -Actin Antibody	Cell Signaling Technology	3700
Rabbit anti- β -Actin Antibody	Cell Signaling Technology	4970
Mouse anti- β -Tubulin Antibody	Genscript	A01717
Rabbit anti- β -Tubulin Antibody	Novus Biologicals	NB600-936
Goat anti-GAPDH Antibody	Novus Biologicals	NB300-320
Rabbit anti-AKT1 Antibody	Cell Signaling Technology	29385
Donkey anti-Goat IgG (H+L) Cross-Adsorbed Secondary Antibody, Alexa Fluor 488	Thermo Fisher Scientific	A-11055
Donkey anti-Rabbit IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor 555	Thermo Fisher Scientific	A-31572
Donkey anti-Mouse IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor 647	Thermo Fisher Scientific	A-31571

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8. 紫外捕获结束后，按照标准程序进行芯片清洗和免疫杂交。
- The screenshot shows the Milo software interface with the 'Assay' tab selected. On the left, there's a 'protein simple' logo and a large 'RUN' button. Below it is a graph showing 'Current (mA)' vs 'Time (h)' with a y-axis from 2.4 to 2.54 and an x-axis from 0 to 2. A status bar at the bottom left says 'Program Started Connecting to DAQ-6001...'. In the center, there are four custom scripts (1, 2, 3, 4) and a table for their parameters:

SCRIPT 1	SCRIPT 2	SCRIPT 3	SCRIPT 4
10 s	0 s	0 s	0 s
60 s	0 s	0 s	0 s
240 V	240 V	0 V	0 V
0 s	240 s	0 s	0 s

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1. Protein Expression Heterogeneity with Milo, the First Single-Cell Western System, Application Note, ProteinSimple, a Bio-Techne Brand
2. Simple Western and Single-Cell Western are Proven Technologies, Publication Spotlight, ProteinSimple, a Bio-Techne Brand

WHERE SCIENCE INTERSECTS INNOVATION™

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For research use or manufacturing purposes only. Trademarks and registered trademarks are the property of their respective owners.

AN_Signal-Enhancement-on-Milo_STRY0086829