

Wes/Jess 试剂耗材价格表 2020

Item	Part No.	Description	RMB Price
Protein Normalization Assay Module for Jess 总蛋白归一化试剂盒	AM-PN01	Contains 12-230 kDA separation module (8x 25-capillary) and Protein Normalization reagents.	¥12,534
Total Protein Detection Module for Jess /Wes , Peggy Sue or Sally Sue 总蛋白检测试剂盒	DM-TP01	Use with a Wes or Sally Sue/Peggy Sue Separation Module. Includes : Luminol-S, Peroxide, Antibody Diluent 2, Total Protein Reconstitution Agent 2, Total Protein Streptavidin-HRP.	¥2,767
Protein Normalization Module for Jess Jess 总蛋白归一化试剂盒	DM-TP02	Use with a Wes/Jess Separation Module. Includes: Protein Normalization Reagent and Reconstitution Agent	¥2,390
Anti-Rabbit Detection Module for Jess /Wes, Peggy Sue or Sally Sue 兔抗检测试剂盒	DM-001	Use with a Wes or Sally Sue/Peggy Sue Separation Module. Includes : Luminol-S, Peroxide, Antibody Diluent 2, Anti-Rabbit Secondary Antibody, Streptavidin-HRP.	¥2,767
Anti-Mouse Detection Module for Jess /Wes, Peggy Sue or Sally Sue 鼠抗检测试剂盒	DM-002	Use with a Wes or Sally Sue/Peggy Sue a Separation Module. Includes : Luminol-S, Peroxide, Antibody Diluent 2, Anti-Mouse Secondary Antibody, Streptavidin-HRP.	¥2,767
No secondary Detection Module for Jess /Wes, Peggy Sue or Sally Sue 无二抗检测试剂盒	DM-003	Use with a Wes or Sally Sue/Peggy Sue a Separation Module. Includes : Luminol-S, Peroxide, Antibody Diluent 2, Streptavidin-HRP.	¥2,219
Biotin Detection Module for Jess /Wes, Peggy Sue or Sally Sue 生物素检测试剂盒	DM-004	Use with a Wes or Sally Sue/Peggy Sue a Separation Module. Includes: Secondary Streptavidin HRP, Luminol-S, Peroxide, Antibody Diluent 2, and Streptavidin-HRP (for biotinylated ladder).	¥2,767
Anti-Human IgG Detection Module for Jess /Wes, Peggy Sue or Sally Sue 抗人 IgG 检测试剂盒	DM-005	Use with a Wes or Sally Sue/Peggy Sue a Separation Module. Includes: Anti-Human IgG Secondary HRP Antibody, Luminol-S, Peroxide, Antibody Diluent 2, Streptavidin-HRP	¥2,767
Anti-Goat Detection Module for Jess /Wes, Peggy Sue and Sally Sue 抗羊检测试剂盒	DM-006	Use with a Wes or Sally Sue/Peggy Sue Separation Module. Includes: Anti-Goat Secondary HRP Antibody, Luminol-S, Peroxide, Milk-free Antibody Diluent (use to dilute goat primary antibody and as blocking agent), and Streptavidin-HRP. Upon receipt, store at 2-8°C.	¥2,767
Anti-Rabbit NIR Detection Module for Jess 兔抗 NIR 检测试剂盒	DM-007	Use with a Wes/Jess Separation Module. Includes: Anti-Rabbit NIR Secondary Antibody, Antibody Diluent 2 and Streptavidin-NIR	¥1,400

Anti-Rabbit IR Detection Module for Jess 兔抗 IR 检测试剂盒	DM-008	Use with a Wes/Jess Separation Module. Includes: Anti-Rabbit IR Secondary Antibody, Antibody Diluent 2 and Streptavidin-NIR	¥1,400
Anti-Mouse NIR Detection Module for Jess 鼠抗 NIR 检测试剂盒	DM-009	Use with a Wes/Jess Separation Module. Includes: Anti-Mouse NIR Secondary Antibody, Antibody Diluent 2 and Streptavidin-NIR	¥1,400
Anti-Mouse IR Detection Module for Jess 鼠抗 IR 检测试剂盒	DM-010	Use with a Wes/Jess Separation Module. Includes: Anti-Mouse IR Secondary Antibody, Antibody Diluent 2 and Streptavidin-NIR	¥1,400
RePlex Module Jess 洗脱试剂盒	RP-001	This RePlex Module provides all the reagents needed to remove primary and secondary antibodies in a RePlex assay with two immunoassays or immunoassay with Total Protein assay in a single run. Please note this module can only be used with Jess. Includes RePlex Reagent 1 (store at 18–28 °C) and RePlex Reagent 2 (store at 2–8 °C).	¥1,067
12-230 kDa Jess /Wes Separation Module, 8 x 13 capillary cartridges 12- 230 kDa Wes 分离试剂盒, 8×13 份	SM-W002	"Use with a Detection Module. Reagents for 104 data points includes two parts: Separation box with 8 Capillary Cartridges (13 capillaries per cartridge), 8 Wes 12-230 kDa Pre-filled Plates (pre-filled with Separation Matrix 2, Stacking Matrix 2, Split Running Buffer 2 and Matrix Removal Buffer), Sample Buffer, Wash Buffer (Store at RT), plus Standard Pack 1 clamshell with 8 ea. DTT, FL Standard master mix with 29 kDa System Control protein (Store at 4C)"	¥6,992
12-230 kDa Jess /Wes Separation Module, 8 x 25 capillary cartridges 12- 230 kDa Wes 分离试剂盒, 8×25 份	SM-W004	"Use with a Detection Module. Reagents for 200 data points includes two parts: Separation box with 8 Capillary Cartridges (25 capillaries per cartridge), 8 Wes 12-230 kDa Pre-filled Plates (pre-filled with Separation Matrix 2, Stacking Matrix 2, Split Running Buffer 2 and Matrix Removal Buffer), Sample Buffer, Wash Buffer (Store at RT), plus Standard Pack 1 clamshell with 8 ea. DTT, FL Standard master mix with 29 kDa System Control protein. (Store at 4C)"	¥10,780
66-440 kDa Jess /Wes Separation Module, 8 x 13 capillary cartridges 66-440 kDa Wes 分离试剂盒, 8×13 份	SM-W006	"Use with a Detection Module. Reagents for 104 data points includes two parts: Separation box with 8 Capillary Cartridges (13 capillaries per cartridge), 8 Wes 66-440 kDa Pre-filled Plates (pre-filled with Separation Matrix 3, Stacking Matrix 2, Split Running Buffer 3 and Matrix Removal Buffer), Sample Buffer, Wash Buffer (Store at RT), plus Standard Pack 3 clamshell with 8 ea. DTT, FL Standard master mix with 90 kDa System Control protein (Store at 4C)"	¥6,992

66-440 kDa Jess/Wes Separation Module, 8 x 25 capillary cartridges 66-440 kDa Wes 分离试剂盒, 8×25 份	SM-W008	"Use with a Detection Module. Reagents for 200 data points includes two parts: Separation box with 8 Capillary Cartridges (25 capillaries per cartridge), 8 Wes 66-440 kDa Pre-filled Plates (pre-filled with Separation Matrix 3, Stacking Matrix 2, Split Running Buffer 3 and Matrix Removal Buffer), Sample Buffer, Wash Buffer (Store at RT), plus Standard Pack 3 clamshell with 8 ea. DTT, FL Standard master mix with 90 kDa System Control protein (Store at 4C)"	¥10,780
2-40 kDa Jess /Wes Separation Module, 8 x 13 capillary cartridges 2-40 kDa Wes 分离试剂盒, 8×13 份	SM-W010	Requires Compass v3.0 or higher. Use with a Detection Module. Reagents for 104 data points includes two parts: Separation box with 8 Capillary Cartridges (13 capillaries per cartridge), 8 Wes 2-40 kDa Pre-filled Plates (pre-filled with Separation Matrix 2, Stacking Matrix 2, Split Running Buffer 4 and Matrix Removal Buffer), Sample Buffer, Wash Buffer (Store at RT), plus Standard Pack 5 clamshell with 8 ea. DTT, FL Standard master mix with 26 kDa System Control protein (Store at 4C)	¥6,992
2-40 kDa Jess /Wes Separation Module, 8 x 25 capillary cartridges 2-40 kDa Wes 分离试剂盒, 8×25 份	SM-W012	"Requires Compass v3.0 or higher. Use with a Detection Module. Reagents for 200 data points includes two parts: Separation box with 8 Capillary Cartridges (25 capillaries per cartridge), 8 Wes 2-40 kDa Pre-filled Plates (pre-filled with Separation Matrix 2, Stacking Matrix 2, Split Running Buffer 4 and Matrix Removal Buffer), Sample Buffer, Wash Buffer (Store at RT), plus Standard Pack 5 clamshell with 8 ea. DTT, FL Standard master mix with 26 kDa System Control protein (Store at 4C)"	¥10,780
单一组分:			
8 pk, EZ Standard Pack 1 标准分子量标记	PS-ST01EZ-8	Includes 8 packets with one each of the Single Use Biotinylated Ladder (MW ladder spans 12-230 kDa) and Single Use FL Standard (29 kDa system control protein) vials. Reagents for 200 Simple Western data points. For use with Jess/Wes, Sally Sue, or Peggy Sue.	¥2,920
8 pk, EZ Standard Pack 2 标准分子量标记	PS-ST02EZ-8	Includes 8 packets with one each of the Single Use Biotinylated Ladder (MW ladder spans 12-230 kDa) and Single Use FL Standard (180 kDa system control protein) vials. Reagents for 200 Simple Western data points. For use with Jess/Wes Sally Sue, or Peggy Sue.	¥2,920
8 pk, EZ Standard Pack 3 大分子量标记	PS-ST03EZ-8	Includes 8 packets with one each of the Single Use Biotinylated Ladder (MW ladder spans 66-440 kDa) and Single Use FL Standard (90 kDa system control protein) vials. Reagents for 200 Simple Western data points. For use with Jess/Wes Sally Sue, or Peggy Sue.	¥2,920
8 pk, EZ Standard Pack 4 大分子量标记	PS-ST04EZ-8	Includes 8 packets with one each of the Single Use Biotinylated Ladder (MW ladder spans 66-440 kDa) and Single Use FL Standard (200 kDa system control protein) vials. Reagents for 200 Simple Western data points. For use with Jess/Wes Sally Sue, or Peggy Sue.	¥2,920
8 pk, EZ Standard Pack 5 小分子量标记	PS-ST05EZ-8	Includes 8 packets with one each of the Single Use Biotinylated Ladder (MW ladder spans 2-40 kDa) and Single Use FL Standard (26 kDa system control protein) vials. Reagents for 200 Simple Western data points. For use with Jess/Wes Sally Sue, or Peggy Sue.	¥2,920

Sample Buffer 样品缓冲液	042-195	10X Sample Buffer (440 uL)	¥590
Luminol-S 鲁米诺发光底物	043-311	Luminol-S (2 mL)	¥1027
Peroxide 过氧化物	043-379	Peroxide (2 mL)	¥1027
Anti-Mouse Secondary Antibody 二抗（抗小鼠）	042-205	Goat Anti-Mouse Secondary HRP Conjugate (For Wes, Jess, Sally Sue, and Peggy Sue only)	¥1740
Anti-Rabbit Secondary Antibody 二抗（抗兔）	042-206	Goat Anti-Rabbit Secondary HRP Conjugate (for Wes, Jess, Sally Sue, and Peggy Sue only).	¥1740
Milk-free Antibody Diluent 抗体稀释液	043-524	Milk-free Antibody Diluent for Wes, Peggy Sue and Sally Sue (20 mL). We recommend its use for dilution of goat primary, anti-goat secondary and as blocking agent as it has been shown to be superior in reducing background compared to other diluents tested. Ships at room temperature. Upon receipt, store at 2-8°C.	¥786
Anti-Goat Secondary HRP Antibody 抗羊二抗	043-522-2	Donkey Anti-Goat Secondary HRP Conjugate (ready-to-use reagent). For Wes, Sally Sue, and Peggy Sue (2 x 1 mL). Use our Milk-free Antibody Diluent (043-524) to dilute goat primary and as a blocking agent. Ships at room temperature. Upon receipt, store at 2-8°C.	¥974
Wes Antibody Diluent 2 抗体稀释液	042-203	Antibody Diluent 2 for Wes (20 mL). Ships overnight at room temperature. Antibody Diluent is recommended for dilution of primary and secondary antibodies. Upon receipt, store at 4C.	¥695