

Wes/Abby/Jess耗材试剂2022促销价格

Item	Part No.	Description	价格	促销价
Total Protein Detection Module for Jess /Wes , Peggy Sue or Sally Sue 总蛋白检测试剂盒	DM-TP01	Use with a Wes or Sally Sue/Peggy Sue Separation Module. Includes : Luminol-S, Peroxide, Antibody Diluent 2, Total Protein Reconstitution Agent 2, Total Protein Streptavidin-HRP.	¥2,767	¥1,998
Anti-Rabbit Detection Module for Jess /Wes, Peggy Sue or Sally Sue 抗兔检测试剂盒	DM-001	Use with a Wes or Sally Sue/Peggy Sue Separation Module. Includes : Luminol-S, Peroxide, Antibody Diluent 2, Anti-Rabbit Secondary Antibody, Streptavidin-HRP.	¥2,767	¥1,998
Anti-Mouse Detection Module for Jess /Wes, Peggy Sue or Sally Sue 抗鼠检测试剂盒	DM-002	Use with a Wes or Sally Sue/Peggy Sue a Separation Module. Includes : Luminol-S, Peroxide, Antibody Diluent 2, Anti-Mouse Secondary Antibody, Streptavidin-HRP.	¥2,767	¥1,998
No secondary Detection Module for Jess /Wes, Peggy Sue or Sally Sue 无二抗检测试剂盒	DM-003	Use with a Wes or Sally Sue/Peggy Sue a Separation Module. Includes : Luminol-S, Peroxide, Antibody Diluent 2, Streptavidin-HRP.	¥2,219	¥1,609
Biotin Detection Module for Jess /Wes, Peggy Sue or Sally Sue 生物素检测试剂盒	DM-004	Use with a Wes or Sally Sue/Peggy Sue a Separation Module. Includes: Secondary Streptavidin HRP, Luminol-S, Peroxide, Antibody Diluent 2, and Streptavidin-HRP (for biotinylated ladder).	¥2,767	¥1,998
Anti-Human IgG Detection Module for Jess /Wes, Peggy Sue or Sally Sue 抗人IgG检测试剂盒	DM-005	Use with a Wes or Sally Sue/Peggy Sue a Separation Module. Includes: Anti-Human IgG Secondary HRP Antibody, Luminol-S, Peroxide, Antibody Diluent 2, Streptavidin-HRP	¥2,988	¥2,189
Anti-Goat Detection Module for Jess /Wes, Peggy Sue and Sally Sue 抗羊检测试剂盒	DM-006	Use with a Wes or Sally Sue/Peggy Sue Separation Module. Includes: Anti-Goat Secondary HRP Antibody, Luminol-S, Peroxide, Milk-free Antibody Diluent (use to dilute goat primary antibody and as blocking agent), and Streptavidin-HRP. Upon receipt, store at 2-8°C.	¥2,767	¥1,998
Anti-Rabbit NIR Detection Module for Jess 抗兔NIR检测试剂盒	DM-007	Use with a Wes/Jess Separation Module. Includes: Anti-Rabbit NIR Secondary Antibody, Antibody Diluent 2 and Streptavidin-NIR	¥1,400	¥999
Anti-Rabbit IR Detection Module for Jess 抗兔IR检测试剂盒	DM-008	Use with a Wes/Jess Separation Module. Includes: Anti-Rabbit IR Secondary Antibody, Antibody Diluent 2 and Streptavidin-NIR	¥1,400	¥999
Anti-Mouse NIR Detection Module for Jess 抗鼠NIR检测试剂盒	DM-009	Use with a Wes/Jess Separation Module. Includes: Anti-Mouse NIR Secondary Antibody, Antibody Diluent 2 and Streptavidin-NIR	¥1,400	¥999
Anti-Mouse IR Detection Module for Jess 抗鼠IR检测试剂盒	DM-010	Use with a Wes/Jess Separation Module. Includes: Anti-Mouse IR Secondary Antibody, Antibody Diluent 2 and Streptavidin-NIR	¥1,400	¥999
RePlex Module Jess洗脱试剂盒	RP-001	This RePlex Module provides all the reagents needed to remove primary and secondary antibodies in a RePlex assay with two immunoassays or immunoassay with Total Protein assay in a single run. Please note this module can only be used with Jess. Includes RePlex Reagent 1 (store at 18–28 °C) and RePlex Reagent 2 (store at 2–8 °C).	¥1,067	¥588
Chemiluminescent Substrate 化学发光底物	PS-CS01	The chemiluminescent Substrate includes Luminol-S (P/N: 042-233) and Peroxide (P/N: 042-234). Each order contains two bottles of 2.0mL of substrate/bottle.	¥1,870	¥1,018
12-230 kDa Jess /Wes Separation Module, 8 x 13 capillary cartridges (12- 230 kDa Wes分离试剂盒, 8×13份)	SM-W002	Use with a Detection Module. Reagents for 104 data points includes two parts: Separation box with 8 Capillary Cartridges (13 capillaries per cartridge), 8 Wes 12-230 kDa Pre-filled Plates (pre-filled with Separation Matrix 2, Stacking Matrix 2, Split Running Buffer 2 and Matrix Removal Buffer), Sample Buffer, Wash Buffer (Store at RT), plus Standard Pack 1 clamshell with 8 ea. DTT, FL Standard master mix with 29 kDa System Control protein (Store at 4C)"	¥6,992	¥4,298
12-230 kDa Jess /Wes Separation Module, 8 x 25 capillary cartridges (12- 230 kDa Wes分离试剂盒, 8×25份)	SM-W004	Use with a Detection Module. Reagents for 200 data points includes two parts: Separation box with 8 Capillary Cartridges (25 capillaries per cartridge), 8 Wes 12-230 kDa Pre-filled Plates (pre-filled with Separation Matrix 2, Stacking Matrix 2, Split Running Buffer 2 and Matrix Removal Buffer), Sample Buffer, Wash Buffer (Store at RT), plus Standard Pack 1 clamshell with 8 ea. DTT, FL Standard master mix with 29 kDa System Control protein. (Store at 4C)"	¥10,780	¥6,628
66-440 kDa Jess /Wes Separation Module, 8 x 13 capillary cartridges (66-440 kDa Wes分离试剂盒, 8×13份)	SM-W006	Use with a Detection Module. Reagents for 104 data points includes two parts: Separation box with 8 Capillary Cartridges (13 capillaries per cartridge), 8 Wes 66-440 kDa Pre-filled Plates (pre-filled with Separation Matrix 3, Stacking Matrix 2, Split Running Buffer 3 and Matrix Removal Buffer), Sample Buffer, Wash Buffer (Store at RT), plus Standard Pack 3 clamshell with 8 ea. DTT, FL Standard master mix with 90 kDa System Control protein (Store at 4C)"	¥6,992	¥4,298
66-440 kDa Jess/Wes Separation Module, 8 x 25 capillary cartridges (66-440 kDa Wes分离试剂盒, 8×25份)	SM-W008	Use with a Detection Module. Reagents for 200 data points includes two parts: Separation box with 8 Capillary Cartridges (25 capillaries per cartridge), 8 Wes 66-440 kDa Pre-filled Plates (pre-filled with Separation Matrix 3, Stacking Matrix 2, Split Running Buffer 3 and Matrix Removal Buffer), Sample Buffer, Wash Buffer (Store at RT), plus Standard Pack 3 clamshell with 8 ea. DTT, FL Standard master mix with 90 kDa System Control protein (Store at 4C)"	¥10,780	¥6,628
2-40 kDa Jess /Wes Separation Module, 8 x 13 capillary cartridges (2-40 kDa Wes分离试剂盒, 8×13份)	SM-W010	Requires Compass v3.0 or higher. Use with a Detection Module. Reagents for 104 data points includes two parts: Separation box with 8 Capillary Cartridges (13 capillaries per cartridge), 8 Wes 2-40 kDa Pre-filled Plates (pre-filled with Separation Matrix 2, Stacking Matrix 2, Split Running Buffer 4 and Matrix Removal Buffer), Sample Buffer, Wash Buffer (Store at RT), plus Standard Pack 5 clamshell with 8 ea. DTT, FL Standard master mix with 26 kDa System Control protein (Store at 4C)	¥6,992	¥4,298
2-40 kDa Jess /Wes Separation Module, 8 x 25 capillary cartridges (2-40 kDa Wes分离试剂盒, 8×25份)	SM-W012	Requires Compass v3.0 or higher. Use with a Detection Module. Reagents for 200 data points includes two parts: Separation box with 8 Capillary Cartridges (25 capillaries per cartridge), 8 Wes 2-40 kDa Pre-filled Plates (pre-filled with Separation Matrix 2, Stacking Matrix 2, Split Running Buffer 4 and Matrix Removal Buffer), Sample Buffer, Wash Buffer (Store at RT), plus Standard Pack 5 clamshell with 8 ea. DTT, FL Standard master mix with 26 kDa System Control protein (Store at 4C)"	¥10,780	¥6,628
Protein Normalization Module	DM-PN02	Compatible with Jess™ systems only. Includes 8x Protein Normalization Reagent and Protein Normalization Reconstitution Agent. Compatible with Separation Module - sold separately.	¥2,730	¥1,583